



PROFESSIONAL MEDICAL PRODUCTS

OXYGEN CONCENTRATOR 5 L

User manual

ATTENTION: *The operators must carefully read and completely understand the present manual before using the product.*

GIMA 34582



CONTEC MEDICAL SYSTEMS CO., LTD
No.112 Qinhuang West Street,
Economic & Technical Development Zone,
Qinhuangdao, Hebei Province, PEOPLE'S
REPUBLIC OF CHINA
cms@contecmed.com.cn
Made in China



CONTEC21



Prolinx GmbH Brehmstr. 56, 40239,
Duesseldorf, Germany



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

USER NOTICE

Dear users, thanks for purchasing the Oxygen Concentrator (hereinafter referred to as machine).

This Manual is written and compiled in accordance with the council directive IEC 60601-1:2020 and ISO 80601-2-69:2020. The information contained in this document is subject to change without notice.

The Manual describes, in accordance with the machine's features and requirements, main structure, performance, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. As well as the safety procedures to protect both the user and machine. Refer to the respective chapters for details.

Please read the Manual carefully before using this machine. These instructions describe the operating procedures to be followed strictly. This Manual will tell you the operating procedures, possible abnormal operation, possible damage to this machine and danger to personal injury which must be noticed during using this machine. Our company is not responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Owing to the forthcoming renovation, the machine you received may not be totally in accordance with the description of this Manual. We would sincerely regret for that.

This machine is a medical device, which can be used repeatedly.

Warning:

- This Manual only provides technical instruction, please follow doctor's advice for oxygen healthcare.
- The machine doesn't suit for surgery or the salvage of no spontaneous respiratory.
- Keep the machine away from strong magnetic field or electromagnetic interference source.
- Open flames during oxygen therapy are dangerous and is likely to result in fire or death. Do not allow open flames within 2 m of the machine or any oxygen carrying accessories.
- Avoid placing the machine in environment with contaminants or fumes.
- Please refer to the correlative medical literature about the clinical restrictions and contraindications.
- The machine is only used for supplying oxygen, don't use it for first-aid treatment or sustaining life.
- There is a risk of fire associated with oxygen enrichment during oxygen therapy. Do not use the oxygen concentrator or accessories near sparks or open flames.
- To ensure receiving the therapeutic amount of oxygen delivery according to your medical condition: the machine must be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels. And the machine must be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.

Our company reserves the final elucidative right.

Contents

Important information.....	1
Chapter 1 Overview.....	4
1.1 Features.....	4
1.2 Applied range.....	4
1.3 Environment	4
Chapter 2 Principle	5
2.1 Basic principles	5
2.2 Oxygen swing absorbers(PSA).....	5
2.3 Uncertainty of parameter measurement	5
Chapter 3 Technical characteristic	6
3.1 Main performance.....	6
3.2 Main parameters	6
3.3 Safety categories.....	7
Chapter 4 Introduction for parts and functions.....	7
4.1 Parts name.....	7
4.2 Alarm	9
4.3 Functions for accumulating time	9
4.4 Filter	9
4.5 Function of timing shutdown.....	10
4.6 Accessories	10
4.7 Software information.....	10
4.8 Packing.....	10
Chapter 5 Operation	10
Chapter 6 Maintenance, Transport and Storage	11
6.1 Maintenance.....	12
6.2 Transport and storage.....	13
Chapter 7 Troubleshooting	14
Chapter 8 Symbol Meanings	15
Appendix 1 Alarm information.....	17
Appendix 2 EMC guidance and manufacturer's declaration	17

Important information

Please read the Manual carefully before using the machine. Do not use this product or any available optional equipment without first completely reading and understanding these instructions. If you are unable to understand the warnings, cautions or instructions, contact a healthcare personnel before attempting to use this equipment—otherwise, injury or damage may occur.

The safety prompt symbols (such as warning, attention, etc.) described in the Manual apply to all dangerous operations which may result in the loss of personal property, as defined below:

Symbol	Meanings
	Warning: high hazard, improper operation may cause injury and death to person or loss to property.
	Attention: potential hazard, improper operation may cause injury to person or loss to property.

Warning:

-  No modification of this equipment is allowed.
-  Do not modify this equipment without authorization of the manufacturer.
-  If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.
-  Additional monitoring or attention may be required for patients using this device who are unable to hear or see alarms or communicate discomfort.
-  The machine is only applicable for use following the applied range described in the Manual, it can't be used for first-aid treatment or sustaining life.
-  Oxygen therapy in certain circumstances can be hazardous, seeking medical advice before using the machine is advisable.
-  The concentrator should always be kept in the upright position to prevent cabinet damage while being transported.
-  Avoid severe vibration and lying upside down during transportation.
-  To avoid electroshock hazard, maintenance to the machine only can be performed by the personnel appointed or authorized by manufacturer. Users are not permitted to maintain it by themselves.
-  Such as power supply voltage instability, beyond AC220 V±22 V range, please install the regulator before use.
-  Do not lubricate fittings, connections, tubing, or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
-  The oil, grease or such substances contacted with oxygen at a certain pressure may produce strong self-ignition. So they should be placed far away from the machine, pipeline, connector and other equipment related to the oxygen. No lubricant can be used in the machine without prior written consent of the manufacturer.
-  Oxygen is a kind of combustion-supporting material, so don't use the machine in environment with inflammable and explosive goods. And smoking or naked flames (including the electrostatic spark caused by friction) is prohibited during using.
-  Don't place the machine in wet environment, and avoid dripping water or other liquid into the machine.
- ◆ Don't immerse the machine into any liquid. If the machine is splashed or coagulated by water, please stop operating.
- ◆ Avoid using it during bathing. If the patient needs to be provided oxygen continuously, please use

it in the other room which is far away from the bathroom at least 2.5 m.

- ◆ If the machine drops into water or other liquid carelessly, please don't touch, shut off the power immediately and contact the authorized dealer or manufacturer.

⚠ The machine should be used in a well-ventilated room, and it should be placed far away from the wall, furniture or similar goods for more than 10 cm.

⚠ NEVER block the air opening of the product or place it on a soft surface, such as a bed or couch, where the air opening may be blocked. Keep the opening free from lint, hair and the like.

⚠ Keep the machine clean, and don't drip or insert any substance into the outlet.

⚠ To reduce the accidental risk, please obey the following operations:

- ◆ Don't move the machine under the state of working.
- ◆ After connecting to the power, the machine needs to be looked after always.
- ◆ Keep the power cord away from the objects generated heat or heating objects.
- ◆ When the machine is not used, please cut off power. Don't use the device in a place where it is difficult to disconnect the power supply.
- ◆ In order to cut off from the network power supply, the plug must be unplugged when the device is not used.
- ◆ Don't drag the power cord to avoid electric shock hazard.

⚠ It may influence the machine performance when using it near the portable communication equipment.

⚠ Do not use the machine near the device with the frequency of 430 MHz~470 MHz, such as wireless intercom equipment, or it may cause an unexpected interference to the machine and once occurred, the machine need a restart.

⚠ There is a risk of inaccurate results or unexpected interference to use the machine during specific investigations or treatments.

⚠ Please operate the machine following doctor's advice or operation steps described in the Manual, if lack of oxygen or insufficient oxygen concentration, please contact the doctor or dealer immediately, don't adjust it by yourselves.

⚠ Check the machine periodically to make sure that there is no visible damage that may affect patient's safety or monitoring performance. It is recommended that the machine should be inspected weekly at least. When there is obvious damage, or any one of the following situations appears: (1) the power cord or the plug is broken, (2) the machine can't work normally, (3) the machine is broken, please contact the engineer for maintaining.

⚠ Please don't connect the machine with other concentrator or oxygen therapeutic equipment in series or parallel.

⚠ Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.

⚠ It is recommended by the manufacturer that keep the machine working for thirty minutes at least after switching on it, and avoid switching on/off the machine frequently, otherwise it will shorten the life of the machine.

⚠ The oxygen delivery settings of the oxygen concentrator should periodically reassessed for the effectiveness of the therapy.

⚠ Medical disposable nasal cannula is sterilized by ethylene oxide gas. Please don't use it if the package is damaged.

⚠ DO NOT use the machine while examining by MRI and CT, as the induced current may cause burn.

-
- ⚠ The device cannot be used in the MRI environment.
 - ⚠ The disposal of scrap machine and its accessories, packings, wastes and residue should comply with the corresponding national laws and regulations, to avoid polluting to the local environment. And the packaging materials must be placed in the region where the children are out of reaching.
 - ⚠ Waste molecular sieves can be sent back to your distributor or manufacturer for disposal, or disposed of according to relevant laws and regulations.
 - ⚠ Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located. If you intend to smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned the oxygen concentrator off before smoking.
 - ⚠ Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.
 - ⚠ Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions, if the oxygen concentrator is turned on, but not in use; the oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.
 - ⚠ Geriatric, paediatric or any other patient unable to communicate discomfort can require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to the responsible care giver to avoid harm.
 - ⚠ If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
 - ⚠ The oxygen delivery setting has to be determined for each patient individually with the configuration of the equipment to be used, including accessories.
 - ⚠ The proper placement and positioning of the PATIENT interface is critical to the effectiveness of the therapy.
 - ⚠ It is dangerous for children to play with the accessories, ensure the accessories are placed where the children are out of reaching.

Attention:

- 🔔 Keep the machine away from dust, vibration, corrosive or flammable substances, and higher or lower temperature and humidity.
- 🔔 The person who is allergic to silicone, PVC, TPU, TPE or ABS can not use this machine.
- 🔔 Accessories must be routed and secured properly. Don't position the tubes or cords around the neck. Ensure the patient can move freely while wearing the cannula. Ensure the tiny parts be placed away from children to avoid accident swallow.
- 🔔 Keep children and pets away from nasal cannula and tubing to avoid choking or strangulation or in case they may cause an unexpected change on the controller.
- 🔔 Please check the packing before use to make sure the machine and accessories are totally in accordance with the packing list, or else the machine may have the possibility of working abnormally.
- 🔔 Please use the accessories properly, for example, use of a paediatric cannula on an adult patient may cause adverse effect on the therapy.
- 🔔 Ensure the oxygen concentrator, its parts and accessories are specified for use at rated flowrate.
- 🔔 Incompatible parts or accessories can result in degraded performance.

-
-  The responsible organization can be accountable for ensuring the compatibility of the oxygen concentrator and all of the parts or accessories used to connect to the patient before use.
 -  When the machine is carried from cold environment to warm or humid environment, please do not use it immediately, wait four hours at least is recommended.
 -  Never try to sterilize the machine by high temperature or high pressure or steam sterilizing process, please refer to relative chapters in the Manual for clean.
 -  The machine doesn't suit all users, if you can't get satisfactory result, please stop using it.
 -  Date of manufacture: see the label.
 -  During the test normal operation of the OXYGEN CONCENTRATOR will deplete the ambient oxygen inside the environmental chamber if the gas output leaves the environmental chamber. An external air source is required to compensate and monitoring of the oxygen concentration inside the chamber is recommended.
 -  Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
 -  If necessary, our company can provide some information(such as circuit diagrams, component lists, illustrations, calibration methods, etc.), so that the qualified technical personnel of the user can repair the machine components designated by our company.

Chapter 1 Overview

The machine is consisted of main unit, humidifier bottle, flowrate controller, etc. It takes the molecular sieve as the adsorbent, adopts pressure swing absorbers(PSA) directly concentrate the medical oxygen from the air.

The machine takes the advantages of small in volume, light in weight, convenient in moving(as it has the turning truckle), stable in performance, high in safety, easy in operating, low in noise, safe in working, which complies with requirements of medical device.

1.1 Features

- a. Easy and convenient to operate.
- b. Small in volume, light in weight and low in power consumption.
- c. Concentrate the oxygen with high concentration from the air directly, easy to operate.

1.2 Applied range

1.2.1 Intended purpose

The device can be used in medical institutions for supplying oxygen.

1.2.2 Patient Population

Adult and child.

1.2.3 Intended users

Professional medical staff.

1.2.4 Medical indications

The device can be used for adjuvant therapy of diseases (including cardiovascular disease, respiratory system disease, hypoxic disease).

1.2.5 Contraindications

Patients with oxygen poisoning or oxygen allergy.

1.3 Environment

- a. Temperature: 10 °C~40 °C

- b. Relative humidity: $\leq 80\%$
- c. Atmospheric pressure: 860 hPa~1060 hPa
- d. Power supply: AC220 V, 50 Hz
- e. Input power: 450 VA

Attention:

- ① If stored or used the machine outside the temperature and humidity specified by the manufacturer, the system may not achieve the performance standards declared.
- ② The temperature and atmospheric pressure range of oxygen concentration status indicator(OCSI) is consistent with the machine.
- ③ Use of this device at an altitude above 1900m is expected to adversely affect the flowrate and the percentage of oxygen and consequently the quality of the therapy.

Chapter 2 Principle

2.1 Basic principles

The device takes air as the material, adopts pressure swing adsorption (PSA) to generate the oxygen (concentration 90 % ~ 96 %, oxygen 93 percent for short).

2.2 Oxygen swing absorbers(PSA)

Pressure swing adsorption, atmospheric pressure desorption. The compressed air passed the air filter respectively enters the separated electromagnetic valve, then the nitrogen, carbon dioxide, vapour in air is selectively adsorbed by the molecular sieve, and oxygen enrichment passed the separation unit of nitrogen and oxygen to form product gas. When the molecular sieve in the separation unit of nitrogen and oxygen adsorbs near the status of saturation, the compressed air enters the other molecular sieve regenerated to continue adsorbing oxygen. The saturated tower makes the molecular sieve desorb and regenerate by decompressing to atmospheric pressure and introducing some oxygen for cleaning the beds of molecular sieve, to prepare for next adsorption. The separation unit of nitrogen and oxygen consisted in series or parallel achieves the purpose of getting continuous production of oxygen through the PLC system controlling.

Sieve is a porous filtering material and is considered a wear item. Some factors that could affect sieve material life include humidity, temperature, particulates, air contaminates, air intake, vibration and other environmental conditions. The frequency or intensity of use may also affect the effective service life.




Figure 1

2.3 Uncertainty of parameter measurement

Serial number	Test parameter	Value	Device error	Measurement uncertainty
1	Maximum recommended flowrate	5 L/min	$\pm 10\%$	$\pm 2\%$
2	concentration	93%	$\pm 3\%$	$\pm 0.4\%$

Note: The tolerances declared in the documents include the uncertainty of the measurement.

Chapter 3 Technical characteristic

3.1 Main performance

- 1) With turning truckle, easy to move.
- 2) Remove the impurity by built-in filter.
- 3) OCSI function.
- 4) Function of accumulating time.
- 5) Function of timing shutdown.
- 6) Alarm for power failure.
- 7) Alarm for low power supply.
- 8) Alarm for high/low pressure protection.
- 9) Alarm for low flowrate.
- 10) Alarm for high temperature.
- 11) With working indicator.

3.2 Main parameters

- 1) Maximum recommended flowrate: 5 L/min
- 2) Flowrate range: 0.5 L/min~5 L/min
- 3) Oxygen concentration(reach the stated oxygen concentration after turning on the machine for about 30 minutes): $93\% \pm 3\%$ (percentage of volume) when the flowrate at the range of 0.5 L/min~5 L/min.
- 4) Flowrate range when the outlet nominal pressure is 0 and 7 kPa: 0.5 L/min ~ 5 L/min
- 5) Flowrate change in maximum recommended flowrate when back pressure of 7 kPa is applied: <0.5 L/min.
- 6) Outlet pressure: 20 kPa~50 kPa
- 7) Input power: 450 VA
- 8) Working voltage: AC220 V $\pm 10\%$, 50 Hz ± 1 Hz
- 9) Operating noise: ≤ 55 dB(A), at the flowrate of 5 L/min
- 10) When the outlet nominal pressure is 0, the function diagram for oxygen concentration and flowrate is shown as Figure 2:




Figure 2

- 11) The relationship between oxygen concentration and altitude In the plateau environment, with the increase of altitude, the atmospheric pressure gradually decreases, the oxygen intake rate also decreases. At the same flowrate, the output oxygen concentration in the plateau environment is lower than that in the plain environment.




Figure 3

12) Dimension: 508 mm×260 mm×530 mm

13) Weight: 21 kg

3.3 Safety categories

- Device class: class II equipment
- The degree of protection against electroshock: type BF applied part (medical disposable nasal cannula)
- The degree of protection against ingress of water: IP21
- Operation mode: continuous operation
- Applied part for protecting defibrillator discharge effect: no
- Signal input/output parts: no

Chapter 4 Introduction for parts and functions

4.1 Parts name




Figure 4 Front panel view




Figure 5 Rear panel view

- Flowrate controller: anti-clockwise rotate the knob to increase the oxygen flowrate(the maximum recommended is 5 L/min); clockwise rotate it to decrease the oxygen flowrate(the minimum recommended is 0.5 L/min). And the LCD displays the flowrate value.
Note:There will be a delay in the flow of electronic display. Please adjust the flow after the display is stable.
- Power switch: switch on the power to open the machine, then it can work normally, switch off the power, then it stops working. The concentrator is isolated from the SUPPLY MAIN by the Power switch.
- Overload protector: when the operating current exceeds the current value limited by the overload protector, the overload protector disconnects and the machine turns off, then turn off the power switch and cut off the power immediately, after excluding the over-current reasons, press it, the machine will return to normal.
- Power indicator: it is green when the power is turned on; in power failure state, switch on the power switch, the indicator is not light.
- Low concentration indicator: the yellow indicator is light when the oxygen concentration is below 82 %(+3 %), and the red indicator is light when the oxygen concentration is below 60 %(± 5 %).
- Setting button for timing shutdown: press "+", "-" to adjust the time for timing shutdown, range: 0~120 minutes, step: 5 minutes.
- Oxygen outlet: when the machine starts working, it exports high-concentration oxygen with constant speed from this outlet.
- Over-pressure relief valve: when the pressure in the humidifier bottle is too high which is arisen from occlusion or bending of nasal oxygen cannula, the over-pressure relief valve will reduce the pressure in the bottle by releasing the oxygen in the humidifier bottle, to ensure that the machine can work normally.
- Humidifier bottle: connect the top screw cap with the oxygen outlet of the machine. The oxygen entered into the humidifier bottle is humidified, then output the oxygen from the oxygen outlet.
- Air intake filter: be used to purify the air inhaled by the machine, and it is placed in air intake filter box.
- Power socket: the power cord connects with the machine by this socket, to provide stable AC for the machine. Plug power cord to an electrical outlet otherwise you can not use the concentrator.
- Exhaust port: be used for ventilation, heat radiation and waste gas output. It should be kept smooth.



Use the accessories not specified by the manufacturer may affect the machine performance.

4.2 Alarm

A. Alarm for power interruption

When the machine works normally, cut off the mains power, the machine gives an auditory alarm, which will stop after continuing for a moment.

When the alarm occurs, turn off the machine at first. Then make sure the power plug is connected well and there is no power failure. Turn on the machine again, if alarm still exists, please turn off and contact the dealer.

This alarm can be triggered by unplugging the power cord during normal operation.

B. Alarm for low power supply

When the power supply falls below the value necessary to maintain normal operation of the concentrator, it will give a visual (the red indicator flashing (frequency:2 Hz) for alarm) and auditory alarm.

When alarm occurs, better to use a voltage regulator as power supply to the concentrator.

C. Alarm for high/low pressure protection

When the abnormalities (the pressure in the oxygen tank is too low or high) appear, the machine will give a visual (the red indicator flashing (frequency:2 Hz) for alarm) and auditory alarm.

When the alarm occurs, make sure the flowrate is adjusted within normal range. Rolling the knob clockwise to decrease it or anti-clockwise to increase it. Once corrected, turn it off for 60 seconds and then turn power back on.

D. Alarm for low flowrate

When the machine works normally, if the flowrate is lower than 0.5 L/min, the machine will give a visual (the red indicator flashing (frequency:2 Hz) for alarm) and auditory alarm.

When the alarm occurs, rolling the knob anti-clockwise to increase it, and wait with a few minutes.

This alarm can be triggered by setting the flowrate below 0.5 L/min.

E. Alarm for high temperature

When the machine works abnormally, and the air temperature is higher than default setting, the machine will give a visual (the red indicator flashing (frequency:2 Hz) for alarm) and auditory alarm.

When the alarm occurs, check if the cabinet filters requires replacement. Make sure the oxygen concentrator is at least three inches away from walls, draperies or furniture. And don't use it within a short time to make sure the machine cools down adequately.

F. Alarm for low oxygen concentration

When the machine works normally, if the oxygen concentration is less than 82 % (+3 %), the yellow indicator is light, then please contact the suppliers, users can continue to use it, but make sure that there is spare oxygen nearby. When the oxygen concentration is less than 60 % (± 5 %), the red indicator flashes, please turn it off immediately, use the spare oxygen and contact the suppliers in time.

This alarm can be triggered by setting the flowrate to maximum, the yellow indicator light illuminates first, and then the red.

4.3 Functions for accumulating time

When the machine works normally, the LCD displays the accumulated time, unit is hour, when it arrives 99999, the machine stops accumulating.

When using the function of timing shutdown, the LCD can't display the accumulated time. Turn it on again, the LCD will accumulate the time automatically.

4.4 Filter

A filter is placed between the oxygen tank and oxygen outlet, it can filter such particles whose diameter is more than 1 μm , which ensures the oxygen quality.

4.5 Function of timing shutdown

The machine has the function of timing shutdown, user can set the working time (range: 0~120 minutes) according to requirements.

4.6 Accessories

No.	Object	Model
1	Humidifier bottle	M14*1.5
2	Nasal Oxygen Cannula	MB01

4.7 Software information

Software name: CONTEC21

Software model: no

Version: V1

Naming rule for version: <Major enhance software upgrade>.<Minor enhance software upgrade>

The software version can be found in the machine.

4.8 Packing

- 1) An oxygen concentrator
- 2) A User Manual
- 3) A power cord
- 4) A humidifier bottle
- 5) A medical disposable nasal oxygen tube(gift for test machine)
- 6) A filter sponge

Chapter 5 Operation

A. Check the filter sponge

Before turning on the machine, please check the intake filter sponge to make sure that it is clear and dry. Refer to section 6.1 for its maintenance.

B. Connect with humidifier bottle




Figure 6 Humidifier bottle

- 1) Clockwise screw the knob on the top cover of the humidifier bottle, to take out the humidifier bottle.
- 2) Remove the top cover of the humidifier bottle. Inject distilled water or cold water to the humidifier bottle to the position between "MINIMUM" and "MAXIMUM", then tighten the top cover of humidifier bottle.
- 3) Connect the humidifier bottle to the oxygen outlet of the machine in counter clockwise direction.
- 4) Connect the oxygen tube to the oxygen outlet of the machine.

 Please inject water to the humidifier bottle following the Manual, the water volume should not exceed the maximum scale marked on the humidifier bottle, otherwise it may spilled over from the cannula which can cause the user choked.

C. Connect the power cord

Turn off the power switch, connect the machine to the power socket by the power cord.

D. Startup for the machine

Turn on the power switch, the power indicator is light, then the machine can work normally. After the the start-up interface showed as figure 7, it enters the main interface, as figure 8. It displays the status “Normal” and the accumulation time(as “00004h”) normally.




Figure 7 start-up interface



Figure 8 main interface

Adjust the flow controller clockwise to decrease it or anti-clockwise to increase it, to make it in the location recommended by the medical personnel, then the machine will export continuous and stable oxygen.

When it is abnormal, the status displayed in the screen would change, like “EL” “EH” “ELL” “ET”“EP”. And other visual and auditory alarm could occur in the meantime.

E. Start absorbing oxygen after connecting the oxygen tube

Place the cannula over your ears and position the prongs in your nose as instructed.




Figure 9

Gas flow at the outlet of the cannula can be checked while the concentrator is warming up. Wave your hand in front of the nasal prongs. You should be able to hear and feel the flow of gas. If you do not feel the gas flow, check the cannula connection for leaks. Or place the end of the nasal cannula under the surface of half-full cup of water and look for the bubbles.

F. Setting for timing shutdown

After turning on the machine, the accumulated time displays on the LCD, press “+” or “-” to enter timing interface, the original value is 30 minutes, it increases five minutes each pressing “+” once, the maximum is 120 minutes. The timing time decreases five minutes each pressing “-” once, 0 minutes at least. Set a value as medical personnel recommend, when the machine runs to the set time, it will turn off automatically.



Figure 10 original timing time



Figure 11 after pressing button“-”

G. Finish absorbing oxygen

After stopping absorbing oxygen, turn off the power switch, pull out the power plug to cut off the power.

Chapter 6 Maintenance, Transport and Storage

 Don't maintain the machine when operating.

6.1 Maintenance

 Please cut off the power before clean.

A. Clean the humidifier bottle

Rinse the outer surface of the humidifier bottle with clean water, and brush the inside thoroughly until no visible impurities, then install the bottle cap and let it air dry for later use. The clean should be performed by healthcare professionals. If replacement is required, please use the model provided by our company.

Note: Clean before and after each use or clean once a day during continuous use, or clean when dust falls on the outer surface of the humidifier bottle.

B. Clean/replacement for the intake filter sponge



Figure 12



Figure 13

Remove the filter sponge, soak and clean it with neutral cleanser, such as Ruhof band cleanser mixed with clean water. After that, rinse it with clean water until there are no bubbles and impurities left. Let it air dry for later use. The clean should be performed by healthcare professionals. If replacement is required, please use the model provided by our company.

Note: it is recommended that the filter sponge should be cleaned once per 100 hours, reinstall it to the machine after fully drying. Please replace another new filter sponge if the one cleaned is not completely dry if the device need to be used immediately.

 It may cause harm to the machine if the intake filter sponge is not installed to the machine or turning on the machine when the intake filter sponge is not fully dried.

C. Replacement for the intake filter



Figure 14



Figure 15

Check the internal filter material by the plastic shell of the intake filter, it is recommended to replace it when the black area reaches 80% or more. Press the filter cover by hand, open the cover after unlocking, then unplug the intake filter to replace it.

It is recommended that the intake filter should be replaced once per 1000 hours at most. Or the lifetime of the machine could be influenced to a certain degree, even the oxygen concentration could be reduced when the filter is clogged severely.

 The intake filter is a dissipative part which is not applicable for long-term use.

D. Replacement of nasal oxygen tube

The nasal oxygen tube is a sterile and disposable product, do not use it repeatedly or crosswise, or it will be insanitary and may be harmful to your health. User can purchase the nasal oxygen tube certificated by self.

 To avoid cross-infection, it is recommended that each person uses the oxygen tube solely. The oxygen tube is disposable, repeated use have a risk of infection.

The nasal oxygen tube is the only parts or accessories of the oxygen concentrator intend for single use.

Cannula requirements: length of 4-25 ft including all oxygen tubing, with crush-proof and single lumen tubing. For example, adult standard flow (rated for up to 6L/min continuous flow) for lengths up to 7 ft.

Note: The nasal oxygen tube may become contaminated with body fluids or expired gases during both NORMAL CONDITION and SINGLE FAULT CONDITION.

Note: Patients should use disposable nasal oxygen tubes of appropriate specifications. The following action is prohibited: Use of a paediatric cannula on an adult patient.

E. Other Maintenance opints

The machine has 3-year service life, and sieve beds for a period of one year. It is recommended to replace the molecular sieve when it reaches up to the expected hours or the low oxygen concentration alarm appears. Keep the machine working for 30 minutes at least after switching it on, and avoid switching on/off the machine frequently, otherwise it will shorten the life of the machine.

Note: No lubricants other than those recommended by the manufacturer are to be used.

6.2 Transport and storage

- 1) The packed machine can be transported by ordinary conveyance or according to transport contract. During transportation, avoid strong shock, vibration and splashing with rain or snow, and to avoid damage, please keep upright, no turning over.
- 2) The packed machine should be stored in room with no corrosive gases and good ventilation. Temperature: -20 °C ~ +55 °C; Relative Humidity: ≤95 %.
- 3) Unpackaged machine should be stored in a dry area, humidity of 40 % or below is recommended. High humidity may influence the life of the sieve beds.

Chapter 7 Troubleshooting

Trouble	Possible Reason	Solution
Turn on the power switch, the machine doesn't work or the indicator isn't light.	1) The power plug is not inserted well. 2) No live with the power socket. 3) The control board of the machine is broken.	1) Check the power plug. 2) Check the power. 3) Please contact the dealer.
The machine stops after working for a period of time, or the oxygen concentration descends.	1) Occlusion for the air intake or exhaust port 2) The intake filter sponge is dirty. 3) The air intake filter is dirty. 4) Too high temperature. 5) Too low voltage of the machine. 6) The fan doesn't work. 7) The molecular sieve failure.	1) Place the machine to the clear location in room, or check the air intake or exhaust port. 2) Install the intake filter sponge via clean and drying. 3) Replace the filter. 4) Place the machine to the ventilated location. 5) Collocate the manostat to ensure the voltage is in the range of 220 V±10 %. 6) Replace the fan. 7) Contact the suppliers to replace the molecular sieve.
The machine can't export oxygen, and there is no bubble in the humidifier bottle.	1) The knob of flow controller is close. 2) The oxygen tube is kinked or damaged. 3) The top cover of humidifier bottle is not screwed. 4) The machine failure.	1) Open the flow controller and check. 2) Unwrap the knot or replace the oxygen tube. 3) Screw down the top cover of the humidifier tube. 4) Please contact the dealer.
Big noise with abnormal sound.	The machine failure.	Please contact the dealer.
LCD display "ELL"	low flow alarm: kinked or blocked tubing, cannula or humidifier.	1) Inspect for kinks or blockages. Correct, clean or replace item. Inspect the flowrate controller and LCD, if it is lower than 0.5 L/min, rolling the knob anti-clockwise to increase it. 2) Please contact the dealer.
LCD display "EL"	System low pressure	1) Inspect the flowrate controller and LCD, if it is higher than 5 L/min, rolling the knob clockwise to decrease it. 2) Please contact the dealer.
LCD display "EH"	System high pressure	1) Inspect for kinks or blockages. Inspect the flowrate controller, if it is lower than 0.5 L/min, rolling the knob anti-clockwise to increase it. 2) Please contact the dealer.
LCD display "ET"	Overheat alarm: unit overheating due to blocked air intake.	1) Remove and clean cabinet filters. 2) Ensure oxygen concentrator at least three inches away from walls, draperies or furniture. 3) Wait at least an hour with the machine off, then turn on again. 4) Please contact the dealer.
LCD display "EP"	Power supply falls below the value necessary to maintain normal operation.	A voltage regulator is suggested when the power supply is not stable.

Chapter 8 Symbol Meanings

Symbol	Meanings
	Refer to instruction manual/booklet
	O ₂ between 60 % to 82 % (Label on the machine panel)
	Alarm
	No smoking
	No open flame: Fire, open ignition source and smoking prohibited
I/O	Power indicator
	"ON"(power)
	"OFF"(power)
	Decrease the time for timing shutdown
	Increase the time for timing shutdown
0 ⁵ --120min	Timing range: 0 ~ 120 minutes, step: 5 minutes
+	Increase the flowrate
-	Decrease the flowrate
	Overload protection (250 V,3 A)
	Type BF applied part
	Class II applied
IP21	Covering Protection rate
	Manufacturer
SN	Serial number
STERILE EO	Sterilized using ethylene oxide
	Do not re-use
LOT	Lot number

	This way up
	Fragile, handle with care
	Keep in a cool, dry place
	Stacking limit is 2 layers
	Atmospheric pressure limitation
	Temperature limitation
	Humidity limitation
	WEEE disposal
	Use-by date Note: this symbol shall be accompanied by a date to indicate that the device should not be used after the end of the year, month or day shown.
	Date of manufacture
	REGULATION(EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017
	Medical device
	Authorized representative in the European Community
	Oxygen output
	Alternating current
~220V 50Hz	AC power supply, 220V,50Hz
	Imported by
	Unique Device Identifier

Appendix 1 Alarm information

	State	Screen display	Grouping of alarm conditions	Delay time	Priority
1	Power failure	— —	Technical alarm	No delay	Low
2	Low power supply	EP	Technical alarm	2 s	High
3	High temperature	ET	Technical alarm	1 s	High
4	High pressure	EH	Technical alarm	8 s	High
5	Low oxygen concentration	Normal	Technical alarm	10 s	High
		Normal	Technical alarm	10 s	Low
6	Low flowrate	ELL	Technical alarm	10 s	High
7	Low pressure	EL	Technical alarm	8 s	High

Audio characteristic for alarm:

High-level alarm: pulse burst composed of 10 pulses, interval: 3 s, fundamental frequency: 750 Hz. Alarm volume: centered around the device, the alarm sound measured at a radius of 1 m shall not be less than 55 dB (A).

Low-level alarm: single pulse, non-repetitive, fundamental frequency: 750 Hz. Alarm volume: centered around the device, the alarm sound measured at a radius of 1 m shall not be less than 55 dB (A).

Response frequency for buzzer in power interruption: about 0.32 Hz.

Verification of the functionality of the Alarm system:

Operator can verify the functionality of the Alarm system until oxygen concentrator stabilizes(2 minutes or more after machine starts).

For example, to trigger the alarm for low flowrate, operator can adjust the flowrate of the oxygen concentrator below 0.5 L/min, and this will cause a high-level alarm with “ELL” display, as indicated above or in chapter 4.2. Another example, alarm for low oxygen concentration can be triggered by setting the flowrate to maximum, the yellow indicator light will illuminate after several seconds(low-level alarm), and it may become a high-level one if the oxygen concentrate is low enough, that's depend on the actual situation of the machine.

Appendix 2 EMC guidance and manufacturer's declaration

Warnings

- The machine is subject to special EMC precautions and it must be installed and used in accordance with these guidelines.
- The electromagnetic field can affect the machine performance, so other equipment used near the machine must meet the corresponding EMC requirements. Mobile phones, X-rays or MRI devices are possible interference source, as they can emit high-intensity electromagnetic radiation.
- The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

The following cable types must be used to ensure that they comply with interference radiation and immunity standards:

Name	Cable length(m)
Power cord	1.8 m

Cable introduction

- The machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, it should be observed to verify normal operation in the configuration in which it will be used.
- Active medical devices are subject to special EMC precautions and they must be installed and used in accordance with these guidelines.
- Portable and mobile RF equipment may affect the use of medical electrical equipment.
- Operation of the machine or system below the minimum amplitude or minimum value may result in inaccurate results.
- Devices or systems may still be interfered by other equipment, even if other equipment meets the requirements of the corresponding national standard.

Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
The CONTEC21 is intended for use in the electromagnetic environment specified below. The customer or the user of the CONTEC21 should assure that it is used in such an environment.	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

Table 2

Guidance and manufacturer's declaration - electromagnetic Immunity		
The CONTEC21 is intended for use in the electromagnetic environment specified below. The customer or the user of the CONTEC21 should assure that it is used in such an environment		
Immunity Test	IEC 60601 Test level	Compliance level
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz
Radiated RF IEC61000-4-3	3 V/m 80 MHz – 2,7 GHz	3 V/m 80 MHz – 2,7 GHz
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CONTEC21 is used exceeds the applicable RF compliance level above, the CONTEC21 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CONTEC21.		
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity		
The CONTEC21 is intended for use in the electromagnetic environment specified below. The customer or the user of the CONTEC21 should assure that it is used in such an environment		
Immunity Test	IEC 60601 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV signal input/output	±2 kV for power supply lines Not Applicable
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz	30 A/m 50Hz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 4

Guidance and manufacturer's declaration - electromagnetic Immunity							
The CONTEC21 is intended for use in the electromagnetic environment specified below. The customer or the user of the CONTEC21 should assure that it is used in such an environment							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380 – 390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400 –	Bluetooth, WLAN,	Pulse modulation b)	2	0,3	28

		2570	802.11 b/g/n, RFID 2450, LTE Band 7	217 Hz			
	5240	5100	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5500	—					
	5785	5800					
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.							
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.							
The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = \frac{6}{d} \sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.							

Table 5

Guidance and manufacturer's declaration - electromagnetic Immunity			
The CONTEC21 is intended for use in the electromagnetic environment specified below. The customer or the user of the CONTEC21 should assure that it is used in such an environment			
Radiated fields in close proximity IEC61000-4-39 (Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields)	Test Frequency	Modulation	IMMUNITY Test Level (A/m)
	134.2 kHz	Pulse modulation 2.1 kHz	65
	13.56 MHz	Pulse modulation 50 kHz	7.5

 Warning

- Don't near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CONTEC21, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note:

- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in

another environment to ensure its accuracy



Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment.*

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

OXY-50 BLUETOOTH PULSE OXIMETER with software

Use and maintenance book

ATTENTION: Operators must read and understand this manual completely before using the product.

GIMA 35103

CONTEC MEDICAL SYSTEMS CO., LTD
No. 112 Qinhuang West Street,
Economic & Technical Development Zone,
Qinhuangdao, Hebei Province,
PEOPLE'S REPUBLIC OF CHINA
cma@contemed.com.cn
Made in China

REF CMS60D

EU REP Prolinx GmbH
Brehmstr. 56, 40239, Duesseldorf, Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessato (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

CMS2.782.089(BT,USB)(CE)ESS/1.4 1.4.01.51.322 2025.05

User Notice

Dear users, thank you very much for purchasing the Pulse Oximeter (hereinafter referred to as device).

This Manual is written and compiled in accordance with the council directive MDD93/42/EEC for medical devices and harmonized standards. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

It is a medical device, which can be used repeatedly.

The Manual describes, in accordance with the device's features and requirements, main structure, functions, specifications, correct methods for transportation, assembly, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and device. Refer to the respective chapters for details.

Please read the User Manual carefully before using this device. The User Manual which describes the operating procedures should be followed strictly. Failure to follow the User Manual may cause measuring abnormality, device damage and human injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, human injury and device damage due to users' negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Owing to the forthcoming renovation, the specific products you received may not be totally in accordance with the description of this User Manual. We would sincerely regret for that.

Our company has the final interpretation to this manual. The content of this manual is subject to change without prior notice.

Warnings

Remind that it may cause serious consequences to tester, user or environment.

- Explosive hazard—DO NOT use the device in environment with inflammable gas such as anesthetic.
- DO NOT use the device while examining by MRI or CT, as the induced current may cause burn.
- Do not take the information displayed on the device as the sole basis for clinical diagnosis. The device is only used as an auxiliary means in diagnosis. And it must be used in conjunction with doctor's advice, clinical manifestations and symptoms.
- The maintenance to the device can only be performed by qualified service personnel specified by manufacturer, Users are not permitted to maintain or refit the device by themselves. Unauthorized modification of the device would result unacceptable risk.
- Uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation disturbance users. It is not recommended that the sensor is used on the same finger for more than 2 hours.
- For some special users who need a more careful inspection on the test site, please don't place the device on the edema or tender tissue.
- Please do not stare at the red and infrared light emitter (the infrared light is invisible) after turning on the device, including the maintenance staff, as it may be harmful to the eyes.
- Each part of the device is firmly fixed, if accidental falling leads to the small parts such as a button to fall off, avoid swallowing of these parts, it may cause suffocation.
- The device contains silicone, PVC, TPU, TPE and ABS materials, whose biocompatibility has been tested in accordance with the requirements in ISO 10993-1, and it has passed the recommended biocompatibility test. The person who is allergic to silicone, PVC, TPU, TPE or ABS can not use this device.
- Do not wrap the SpO₂ probe or USB cable around neck to avoid an accident.
- The disposal of scrap device, its accessories and packaging should follow the local laws and regulations, to avoid polluting to the local environment. And the packaging materials must be placed in the region where the children are out of reaching.
- The device can not be used with the equipment not specified in the Manual. Only the accessories appointed or recommended by the manufacturer can be used, otherwise it may cause injury to the tester and operator or damage to the device.
- The SpO₂ probe accompanied is only suitable for using with the device. The device can only use the SpO₂ probe described in the Manual, so the operator has the responsibility to check the compatibility between the device and the SpO₂ probe before using, incompatible accessories may cause device performance degradation, device damage or patient injury.
- Do not reprocess the accompanying SpO₂ probe.
- Check the device before use to make sure that there is no visible damage that may affect user's safety and device performance. When there is obvious damage, please replace the damaged parts before use.
- When the message "Sensor Off" or "Sensor Fault" appears on the screen, it indicates that the SpO₂ probe is disconnected or line fault occurs. Check the connection of the SpO₂ probe and whether there is damage for the probe, if necessary, please replace the probe to avoid risks. The probe fault will not result in a safety hazard.
- Functional testers can not be used to assess the accuracy of the SpO₂ probe and Pulse Oximeter.
- Some functional testers or patient simulators can be used to verify whether the device works normally, for example, INDEX-2LFE Simulator (software version: 3.00), please refer to the Manual for the detailed operation steps.
- Some functional testers or patient simulators can measure the accuracy of the device copied calibration curve, but they can not be used to evaluate the device accuracy.
- When using the device, please keep it away from the equipment which can generate strong electric field or strong magnetic field. Using the device in an inappropriate environment may cause interference to the surrounding radio equipment or affect its working.
- When storing the device, keep it away from children, pets and insects to avoid affecting its performance.
- Do not place the device in places exposed to direct sunlight, high temperature, humidity, dust, cotton wool or easy to splash water, to avoid affecting its performance.
- The measured accuracy will be affected by the interference of electrosurgical equipment.
- When several products are used on the same patient simultaneously, danger may occur which is arisen from the overlap of leakage current.

- CO poisoning will appear excessive estimation, so it is not recommended to use the device.
- This device is not intended for treatment.
- The intended operator of the device may be a patient.
- Avoid maintaining the device during using.
- Users should read the product manual carefully before use and operate according to the requirements.

1 Overview

The oxygen saturation is the percentage of HbO₂ in the total Hb in the blood, so-called the O₂ concentration in the blood, it is an important physiological parameter for the respiratory and circulatory system. A number of diseases related to respiratory system may cause the decrease of SpO₂ in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body, and the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. Serious symptoms might bring danger to human's life. Therefore, prompt information of patients' SpO₂ is of great help for the doctor to discover the potential danger, and is of great importance in the clinical medical field.

Insert the finger when measuring, the device will directly display the SpO₂ value measured, it has a higher accuracy and repeatability.

1.1 Features

- A. Easy to use.
- B. Small in volume, light in weight, convenient to carry.
- C. Low power consumption.

1.2 Intended purpose

The Pulse Oximeter is a non-invasive device intended for the spot-check or continuously monitor of oxygen saturation of arterial hemoglobin (SpO₂) and the pulse rate of adult, pediatric and neonate patients through finger in home and hospital environments (including clinical use in internist/surgery, anesthesia, intensive care etc.). The device is not intended for single use and out-of-hospital transport use.

1.3 Environment requirements

Storage Environment

- a) Temperature: -40 °C ~ + 60 °C
 - b) Relative humidity: ≤ 95%
 - c) Atmospheric pressure: 500 hPa ~ 1060 hPa
- Operating Environment
- a) Temperature: +10 °C ~ + 40 °C
 - b) Relative Humidity: ≤ 75%
 - c) Atmospheric pressure: 700 hPa ~ 1060 hPa

1.4 Precautions

1.4.1 Attention

Point out conditions or practices that may cause damage to the device or other properties.

- Before using the device, make sure that it locates in normal working state and operating environment.
- In order to get a more accurate measurement, it should be used in a quiet and comfortable environment.
- When the device is carried from cold or hot environment to warm or humid environment, please do not use it immediately, wait four hours at least is recommended.
- If the device is splashed or coagulated by water, please stop operating.
- DO NOT operate the device with sharp things.
- High temperature, high pressure, gas sterilizing or immersion disinfection for the device is not permitted. Refer to User Manual in the relative chapter (6.1) for cleaning and disinfection. Please take out the internal battery before cleaning and disinfection.
- The device is suitable for adult.
- The device may not be suitable for all users, if you can't get a satisfactory result, please stop using it.
- Data averaging and signal processing have a delay in the upgrade of SpO₂ data values. When the data update period is less than 30 seconds, the time for obtaining dynamic average values will increase, which is arisen from signal degradation, low perfusion or other interference, it depends on the PR value.
- The device has 3-year service life, date of manufacture see the label.
- The expected service life of the attached parts or accessories of the equipment is two year.
- If the shelf life is less than the expected service life, the shelf life of the attached parts or accessories of the equipment is two year.
- The shelf life shall prevail if the accessory or part has a shelf life that is shorter than its intended service life.
- The device does not provide over-limit alarm function for SpO₂ and PR, so it is inapplicable for using in the place where need such function.

- This device has the function of prompting, users can check on this function according to chapter 5.5.1 as a reference.
- The device has the function of limits prompting, when the measured data is beyond the highest or lowest limit, the device would start prompting automatically on the premise of the prompting function is on.
- The device has the function of prompting, this function can either be paused, or closed for good. This function could be turned on through menu operation if you need. Please check the chapter 5.5.1as a reference.
- The device hasn't low-voltage alarm function, it only shows the low-voltage, please change the battery when the battery voltage is used up.
- The maximum temperature at the SpO₂ probe -tissue interface should be less than 41°C which is measured by the temperature tester.
- During measuring, when abnormal conditions appear on the screen, please pull out your finger and reinsert it to measure again.
- If some unknown error appears during measuring, remove the battery to terminate operating.
- Do not contort or drag the wire of the device.
- The plethysmographic waveform is not normalized, as a signal inadequacy indicator, when it is not smooth and stable, the accuracy of the measured value may degrade. When it tends to be smooth and stable, the measured value read is the optimal and the waveform at this time is also the most standard.
- If necessary, please visit our official website to get the information about SpO₂ probe that can be used with this device.
- If the device or component is intended for single-use, then the repeated use of these parts will pose risks on the parameters and technical parameters of the equipment known to the manufacturer.

- If necessary, our company can provide some information (such as circuit diagrams, component lists, illustrations, etc.), so that the qualified technical personnel of the user can repair the device components designated by our company.
- Operators can contact our company to obtain the modified Bland and Altman plot.
- The measured results will be influenced by the external colouring agent (such as nail polish, colouring agent or color skin care products, etc.), so don't use them on the test site.
- Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO₂ and pulse rate.
- As to the fingers which are too cold or too thin or whose fingernail is too long, it may affect the measured results, so please insert the thicker finger such as thumb or middle finger deeply enough into the probe when measuring.
- The finger should be placed correctly (see Attached figure 6), as improper assembly or improper contact position for sensor will influence the measurement.
- The light between the photoelectric receiving tube and the light-emitting tube of the device must pass through the subject's arteriole. Make sure the optical path is free from any optical obstacles like rubberized fabric, to avoid inaccurate results.
- Excessive ambient light may affect the measured results, such as surgical light (especially xenon light sources), bilirubin lamp, fluorescent lamp, infrared heater and direct sunlight, etc. In order to prevent interference from ambient light, make sure to place the sensor properly and cover the sensor with opaque material.

- Frequent movement (active or passive) of the subject or severe activity can affect the measured accuracy.
- The SpO₂ probe should not be placed on a limb with the blood pressure cuff, arterial ductus or intraluminal tube.
- The measurement value may be inaccurate during defibrillation and in a short period after defibrillation, as it has not defibrillation function.
- The device has been calibrated before leaving factory.
- The device is calibrated to display functional oxygen saturation.
- The equipment connected with the Oximeter interface should comply with the requirements of IEC 60601-1.

1.4.2 Clinical restriction

- A. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, Raynaud's syndrome, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- B. The measurement will be influenced by intravascular staining agents (such as indocyanine green or methylene blue), skin pigmentation.
- C. The measured value may be normal seemingly for the tester who has anemia or dysfunctional hemoglobin(such as carboxyhaemoglobin (COHb), methaemoglobin (MetHb) and sulphaemoglobin (SuHb)), but the tester may appear hypoxia, it is recommended to perform further assessment according the clinical situations and symptoms.
- D. Pulse oxygen only has a reference meaning for anemia and toxic hypoxia, as some severe anemia patients still show better pulse oxygen measured valued.
- E. Contraindication:
 - a. The person who is allergic to silicone, PVC, TPU TPE or ABS.
 - b. The damaged skin tissue.
 - c. During cardiopulmonary resuscitation.
 - d. When the patient is hypovolemic.
 - e. For assessing the adequacy of ventilatory support.
 - f. For detecting worsening lung function in patients on a high concentration of oxygen.

1.5 Clinical indications

The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger.

2 Principle




Figure 1. Operating principle

An experience formula of data processing is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO₂) in red light & near-infrared light zones. On the basis of the principle of Photoelectric Oxyhemoglobin Inspection Technology and Photoplethysmography technology, it uses two light beams of different wavelengths to irradiate the human fingertip to obtain the measurement information from the photosensitive element, after processed by the electronic circuits and microprocessor, displays the measured results on the screen.

3 Functions

- A. SpO₂ value display
- B. PR value and bar graph display
- C. Pulse waveform display
- D. Low-battery indication: low-battery indication appears when the battery voltage is too low to work
- E. Adjustable screen brightness
- F. PR sound indication
- G. Voice prompt for over-limit, sensor off /finger-out and low battery
- H. With SpO₂ value and pulse rate value record function, the stored data can be uploaded to computer.
- I. It can be connected with an external oximeter probeReview function
- J. Clock function

4 Product Introduction

4.1 Appearance




Figure 2. Front View

4.2 Battery assembly




Figure 3. Batteries Assembly

- A. Refer to Figure 3. Use a screwdriver to unscrew the two screws from the battery compartment on the back of the product and open the back cover of the battery compartment.
- B. Insert the two AA size batteries properly in the right direction.
- C. Close the battery back cover, screw on the screw.

Please take care when you insert the batteries, for the improper insertion may damage the device.

Please replace two new batteries of the same kind at the same time.

4.3 Probe assembly

Inserting the SpO₂ probe of the pulse oximeter in the upper jack(see Figure 4). (The probe is limited to be produced by our company; never replace it with the similar one by other manufacturers).




Figure 4. Probe Installation




Figure 5. USB Port

when inserting zhe probe, make the protruding part of the probe plug correspond to the groove of the probe socket. Pull out the probe directly and don't rotate the probe.

4.4 USB port

It is used to connect a personal computer to export the trend data(see Figure 5).

4.5. Structure, accessories and software description

- A. Structure: main unit, SpO₂ probe, USB cable.
 - B. Accessories: one SpO₂ probe, two AA size batteries (optional), one USB cable, one User Manual.
- Please check the device and accessories according to the list to avoid that the device can not work normally.
- C. Software description
- Release version: V2

5 Operating

5.1 Application method

- A. Insert the finger into the probe as shown in Figure 6.




Figure 6. Sketch map for finger placement

(The appearance of actual probe may be different with the one shown as Figure 6, please refer to the actual probe.)

when inserting the finger, the light emitting from the sensor must be directly irradiated to the side of the fingernail.

In the process of using the tested finger had better not shake, the human body also had better not be in motion state.

- Long press the "power on/off " button, until the device turns on.
- Do not shake the finger and keep the user in a stable state during the measurement process.
- Wait a few seconds, the data can be read directly from the screen in the measure interface.

5.2 Pause sound prompt

- A. Sound prompt, including: over-limit, low-battery, finger out, sensor off and sensor fault.
- B. Under the measurement interface, turn on the sound prompt, when the sound prompt occurs, short press the button to pause the sound prompt, and it will resume automatically after about 60s.
- C. If you want to turn off the sound prompt permanently, please set it in menu.

5.3 Review Interface

- A. In the measure interface, press "up button" to enter the Review Interface 1 directly, as shown in Figure 7:




Figure 7-1. Review Interface 1

Figure 7-2. Review Interface 2

- B. In review interface, press "menu button" to switch between Review Interface 1 and Review Interface 2.
- C. In Review Interface 1, the user can observe the trend waveform composed by storage data. Each screen can show storage data for 105 seconds. The yellow line shows the SpO₂ trend waveform, and the red line shows the PR trend waveform. The time underside shows the starting time of displaying the date in the screen, press the "left button" or "right button" to view the information on the previous or next page of the stored data trend chart.
- D. The Review Interface 2 shown based in Review Interface 1, the stored SpO₂ value and PR value in each second can be observed here, the underside date from left to right marks time, SpO₂ value, PR value. Press "left button" or "right button" to display the blood oxygen and pulse of the previous or next second; Long press the "left button" or "right button", and the pulse and blood oxygen will be display with a data interval of 10 seconds.
- E. Press "up button" to exit the review Interface, return to the measure interface.

5.4 Clock interface

In the measure interface, press the "right button" can enter the clock interface of Figure 8. Press the "right button" again can return to the measure interface.




Figure8. Clock interface




Figure 9. Main Menu




Figure 10. Setting for sound prompt

5.5 Menu operations

In the measure interface, press the "menu button "can enter the menu of Figure 9. Users can adjust the setting through the main menu, such as the sound prompt, record, clock., system, etc. can be set, methods are as follows:

5.5.1 Sound prompt setting

Under main menu, press the "up button" or "down button" to select "Prompt", then press the "left button" or "right button" to enter its setting interface shown in Figure 10.

Press the "up button" or "down button" to select the option to be adjusted, then press the "left button" or "right button" to change the value.

“SpO₂ HI(%)”: upper limit prompt for SpO₂ over-limit
“SpO₂ LO(%)”: lower limit prompt for SpO₂ over-limit
“PR HI(bpm)”: upper limit prompt for PR over-limit
“PR LO(bpm)”: lower limit prompt for PR over-limit
“Prompt Sound”: prompt for over-limit, low-battery, finger out, sensor off and sensor fault, “off”: close, “on”: open.
“Pulse Sound”: PR sound, “off”: close, “on”: open.
Lower limit can not exceed the upper limit, and the upper limit can not be lower than the lower limit when adjusting the values. SpO₂ range: 0 % ~ 100 %, PR range: 0 ~ 254 bpm
The values displayed in Figure 10 are the initial values of over-limit prompt.

After setting, press the "menu button" to exit the Prompt Settings Menu interface, and return to “Main Menu” interface.
5.5.2 Data storage
Under the main menu, press the "up button" or "down button" to select “Record”, then press the "left button" or "right button" to enter the “Record Menu” interface as shown in Figure 11.
Press the "up button" or "down button" to select the option to be adjusted, then press the "left button" or "right button" to change the value.
It indicates that the device is storing when the red dot “REC●” in measurement interface flickers
“Mode”: record mode selection, including: “Auto” and “Manual” mode. Under “Manual” mode, select to turn on / off memory by “Record”.

Auto record: start recording after stable data appear, pull out the finger to finish recording a group of data (99 group of data at most), the total duration does not exceed 72 hours.
Manual record: after manual storage is started, the storage state needs to be terminated manually to complete a group of store, store up to 24-hour data.

When the memory is full, it will display “Memory is full!”, then it will enter the standby mode after several seconds. When exiting the standby mode next time, it will display “Memory is full!” to prompt user that the memory has been full, then it will enter the measure interface.
△Under manual mode, when “Record” is “ON”, the device will prompt to clear the data stored last time.
It will display “Recording...” when there is no operation under record state for 15s, then it will enter energy saving mode after several seconds, pressing the "power on/off button", the device would return to the former interface; pressing any button(power on/off excluded), it will display “Recording...”.
△ Under data recording state, after the display screen turns off automatically, in order to save power, pulse sound indication will turn off automatically.
“Seg”: data segment.
After setting, press the "menu button" to exit storage menu, return to main menu.
“Delete All”: delete all records (auto record mode is shown as Figure 11).
△ Please upload data in time after recording, otherwise the data may be covered when the storage space is full.
△ The historical data will be deleted once switching the mode. Under record state, the record mode can not be switched; under manual mode, the record mode can be switched only when turning off recording firstly.

Record Menu	
Mode	Auto
Seg	12
Delete All	

Figure 11. Record menu

Clock Menu	
Set Time	no
Set Year	2019
Set Month	01
Set Day	01
Set Hour	03
Set Minute	00
MENU→EXIT	

Figure 12. Clock menu

System Menu	
Hard.Ver.	2.0.0
Soft.Ver.	2.0.2
ID	user
Demo	off
Sound Volume	3
Brightness	1
MENU→EXIT	

Figure 13. System menu

5.5.3 Clock setting
a. Connect the master device to synchronize device time
Under the PC software interface, after search for the device (refer to relative chapter (5.6) for the connection method), then can synchronize the device time.
b. Set device time manually
Under main menu, press the "up button" or "down button" to select “Clock”, then press the "left button" or "right button" to enter its setting interface shown in Figure 12.
Press the "up button" or "down button" to select the option to be adjusted, then press the "left button" or "right button" to change the value.
“Set Time”: set the time, “yes”: allow, “no”: prohibit
“Set Year”: set the year
“Set Month”: set the month
“Set Day”: set the day
“Set Hour”: set the hour
“Set Minute”: set the minute

Adjustable range for year: 2015 ~ 2045, month: 1 ~ 12, day: 1 ~ 30 (when there are 31 days in a month, it is 1 ~ 31), hour: 1 ~ 23, minute: 1 ~ 59.
After setting, press the "menu button" to exit clock menu, return to main menu.
5.5.4 System setting and other options introduction
Under main menu, press the "up button" or "down button" to select “System”, then press the "left button" or "right button" to enter the interface as shown in Figure 13.
Press the "up button" or "down button" to select the option to be adjusted, then press the "left button" or "right button" to change the value.
“Hard.Ver”: hardware version
“Soft.Ver”: software version
“ID”: user name
“Demo”: set the Demo mode, “on”: turn on the Demo mode, “off”: turn off the Demo mode.
“Sound Volume”: set the sound volume, adjustable range: 1 ~ 3
“Brightness”: set the screen brightness, adjustable range: 1 ~ 4
After setting, press the "menu button" to exit system setting menu, return to main menu.

5.5.5 Exit main menu
Under main menu, press the "menu button" to exit the main menu and return to the measurement interface.
5.6 Data upload
Connect the device to the computer by the USB cable, upload the data after connecting the PC software properly, refer to “Software operating instruction” for details.
△The PC software can be downloaded from our official website.
5.7 Power off
Long press the "power on/off" button, until the device turns off.
△when the device is in storing , it can't be turned off.

6 Maintain, Transport and Storage
6.1 Cleaning and disinfection
The device must be turned off before cleaning, and it should not be immersed into liquid.
Please take out the internal battery before cleaning, do not immerse it into liquid.
Use 75% alcohol to wipe the device enclosure, nature dry or clean it with clean and soft cloth. Do not spray any liquid on the device directly, and avoid liquid penetrating into the device.

6.2 Maintenance
A. Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient’s safety and monitoring performance. It is recommended that the device should be inspected weekly at least. When there is obvious damage, stop using it.
B. Please clean and disinfect the device before/after using it according to the User Manual (6.1).
C. Please replace the batteries in time when low-battery appears.
D. Please take out the batteries if the device is not used for a long time.
E. The device need not to be calibrated during maintenance.
6.3 Transport and Storage
A. The packed device can be transported by ordinary conveyance or according to transport contract. During transportation, avoid strong shock, vibration and splashing with rain or snow, and it can not be transported mixed with toxic, harmful, corrosive material.
B. The packed device should be stored in room with no corrosive gases and good ventilation.

Trouble	Possible Reason	Solution
The values can not be displayed normally or stably.	1) The finger is not properly inserted. 2) The finger is shaking or the patient is moving. 3) The device is not used in environment required by the manual. 4) The device works abnormally.	1) Please insert the finger properly and measure again. 2) Let the patient keep calm. 3) Please use the device in normal environment. 4) Please contact the after-sales.
The device can not be turned on	1) The battery is drained away or almost drained away. 2) The battery is assembled incorrectly. 3) The device’s malfunction.	1) Please change batteries. 2) Please assemble the battery again. 3) Please contact the local service center.
The display disappears suddenly.	1) The device enters into the energy saving mode. 2) Low battery. 3) The device works abnormally.	1) Normal. 2) Please change batteries. 3) Please contact the after-sales.
The data can not be stored.	1) The device is not operated according to the manual. 2) The device works abnormally.	1) Please operate the device according to the manual. 2) Please contact the after-sales.

Symbols	Meaning	Symbols	Meaning
	Refer to instruction manual/booklet for information related to safety.	PRbpm	Pulse rate (bpm)
	Type BF applied part	%SpO ₂	Pulse oxygen saturation (%)
	Manufacturer		The battery power is full
	Two grid of the battery		One grid of the battery
	The lack of battery power.(Please change batteries in time for exact measuring)		Serial number
	Use-by date		Recycling garbage WEEE (2012/19/EU)
	USB		Battery anode
	Battery cathode		It means this pulse oximeter is protection against ingress of solid foreign objects of 12.5mm and greater, protection against ingress of water dripping when device tilted up to 15°.
	Temperature limitation		Humidity limitation
	Atmospheric pressure limitation		This way up
	Fragile, handle with care		Keep away from rain
	Close the sound prompt		Pause the sound prompt
	Open the sound prompt		Close the pulse sound indication
	Open the pulse sound indication		Power on/off button
	Menu button		left button/prompt pause button
	Right button/clock button		Up button/replay button
	down button		Recyclable
	Manufacture Date	Sensor Off	The probe is disconnected.
Finger Out	The finger is not inserted.	Sensor Fault	Probe failure
REC●	Recording	P/N	Material code
	Batch No.		1. The finger clip falls off (no finger inserted) 2. Probe error 3.Signal inadequacy indicator
	Alarm inhibit		Medical device
	Authorized representative in the European Community		This item is compliant with Directive 93/42/EEC of 14 june 1993 concerning medical devices; Including, at 21 march 2010, the amendments by Council Directive 2007/47/EC.
	unique device identifier		Stacking limit by number
	Caution: read instructions (warnings) carefully		Imported by

Note: Your device may not contain all the following symbols.

9 Specification	
SpO ₂ [see note 1]	
Display range	0% ~ 100%
Measured range	0% ~ 100%
Accuracy [see note 2]	70%~100%: ±2%; 0%~69%: unspecified.
Resolution	1%
PR	
Display range	30 bpm ~ 250 bpm
Measured range	30 bpm ~ 250 bpm
Accuracy [see note 3]	±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm
Resolution	1 bpm
Accuracy under low perfusion [see note 4]	Low perfusion 0.4%: SpO ₂ : ±4%; PR: ±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm.
Light interference	Under normal and ambient light conditions, the SpO ₂ deviation ≤ 1%
Pulse intensity	Continuous bar graph display, the higher display indicates the stronger pulse.
Upper and lower limit of measured values	
SpO ₂	0% ~ 100%
PR	0 bpm ~ 254 bpm
Optical sensor [see note 5]	
Red light	Wavelength: about 660 nm, optical output power: < 6.65 mW
Infrared light	Wavelength: about 905 nm, optical output power: < 6.75 mW
Memory	Up to 99 group of data under auto mode, total duration does not exceed 72 hours. Up to 24-hour data under manual mode.
Safety class	Internally powered equipment, type BF applied part
International Protection	IP22
Working voltage	DC 2.6 V - 3.6V
Working current	≤ 100 mA
Power supply	Dry battery (2AA)
Operation time	The device can continuously work for 20 hours when it was powered by two new batteries within the warranty period.
Dimension and Weight	
Dimension	110(L) × 60(W) × 24(H) mm
Weight	About 120g (with Dry battery(2AA))

Note 1: the claims of SpO₂ accuracy shall be supported by clinical study measurements taken over the full range. By artificial inducing, get the stable oxygen level to the range of 70 % to 100 % SpO₂, compare the SpO₂ values collected by the secondary standard pulse oximeter equipment and the tested equipment at the same time, to form paired data, which are used for the accuracy analysis.
There are 12 healthy volunteers (male: 6. female: 6; age: 18~50; skin color: dark black: 3, medium dark: 1, light: 7, white: 1) data in the clinical report.
Note 2: because pulse oximeter equipment measurements are statistically distributed, only about two-thirds of pulse oximeter equipment measurements can be expected to fall within ±Arms of the value measured by a CO-OXIMETER.
Note 3: Patient simulator has been used to verify the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator.
Note 4: percentage modulation of infrared signal as the indication of pulsating signal strength, patient simulator has been used to verify its accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions, compare them with the known SpO₂ and PR values of input signal.
Note 5: optical sensors as the light-emitting components, will affect other medical devices applied the wavelength range. The information may be useful for the clinicians who carry out the optical treatment.For example, photodynamic therapy operated by clinician.

Appendix		
State	Prompt condition delay	Prompt signal generation delay
Low voltage prompt	1s	20ms
SpO ₂ prompt	330ms	20ms
Pulse rate prompt	330ms	20ms
Probe error prompt	16ms	20ms

EMC	
This equipmen is suitable for professional healthcare facility environments and home healthcare environments.	

Warning:

- Don’t near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the this equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Note:

- this equipment needs special precautions regarding EMC and needs to put into service according to the EMC information provided below.
- The basic performance: SpO₂ measured range: 70% ~ 100%, absolute error: ±2%; PR measured range: 30 bpm ~ 250 bpm, accuracy:±2 bpm during the pulse rate range of 30 bpm ~ 99 bpm and ±2% during the pulse rate range of 100 bpm ~ 250 bpm.
- When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in another environment to ensure its accuracy.
- Other devices may affect this device even though they meet the requirements of CISPR.

Product configuration:		
Serial number	Name	Cable length
1	SpO ₂ probe	1.5m
2	USB cable	1m

Table 1	
Guidance and Declaration - Electromagnetic Emissions	
Emissions test	Compliance
Radiated RF EMISSIONS CISPR 11	Group 1
Radiated RF EMISSIONS CISPR 11	Class B
Harmonic distortion IEC 61000-3-2	Not applicable
Voltage fluctuations and flicker IEC 61000-3-3	Not applicable

Table 2		
Guidance and Declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	Not applicable
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable
Voltage dips and Voltage interruptions IEC 61000-4-11	0 % U _r ; 0,5 cycle. At0°,45°,90°,135°, 180°,225°,270°and315°. 0 % U _r ; 1 cycle and 70 % U _r ; 25/30 cycles; Single phase: at 0°. 0 % U _r ; 250/300 cycle	Not applicable
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15MHz - 80 MHz 6 V in ISM and amateur radio bands between 0,15MHz to 80 MHz 80%AM at 1kHz	Not applicable
Radiated RF IEC61000-4-3	10V/m 80 MHz-2,7GHz 80%AM at 1kHz	10V/m 80 MHz-2,7GHz 80%AM at 1kHz
NOTE U _r is the a.c.mains voltage prior to application of the test level		

Table 3						
Guidance and manufacturer’s declaration - electromagnetic immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC60601-1-2 Test level (V/m)	Compliance level (V/m)
	385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	27	27
	450	430 - 470	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	28	28
	710	704 - 787	LTE Band 13,17	Pulse modulation b) 217 Hz	9	9
	745					
	780					
	810	800–960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	28	28
	870					
	930					
	1720	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	28	28
	1845					
	1970					
	2450	2400–2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	28	28
	5240					
	5500	5100–5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	9	9
	5785					

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment

GIMA WARRANTY TERMS
The Gima 12-month standard B2B warranty applies

