EVOLVE LINE

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Dear customers, thanks for choosing our Gima cardiology stethoscope. Please read this manual thoroughly and carefully before using this product and follow the instructions.



Modern stethoscope

1. High acoustics sensitivity
2. Floating diaphragm for monitoring low and high frequencies
3. Colorful combination chestpiece with non chill ring
4. Wonderful and durable design
5. Soft sealing eartips

Changing frequencies using the floating diaphragm

Your Gima stethoscope is equipped with a floating diaphragm that enables you to listen to both low and high frequency sounds without turning over the chestpiece.

To listen to low and high frequency sounds without removing and turning over the chestpiece, simply alternate between light and firm pressure on the floating diaphragm.

Headset adjustment

1. Gima stethoscope is equipped with comfortable sealing eartips. They are designed to conform to the ear canal so as to block out ambient noise. The eartips are fixed onto the eartubes for a secure connection.
2. To reduce spring tension, grasp each eartube at its bend and gently pull outwards as shown in Fig 1.
3. To increase spring tension, squeeze the eartubes together as shown in Fig 2.
4. When the eartips are in your ears, the eartubes should be pointed forward as shown in Fig 3. Improperly wearing the stethoscope headset can result in poor acoustical seal, and in some cases, complete sound blockage.


Fig. 1


Fig. 2

 Make sure that the eartubes are pointed slightly forward into the ear canal.

Fig. 3


Correct

Incorrect

Removing the diaphragm and cleaning the chestpiece

With the diaphragm side up, grasp the rim with thumbs and index fingers of both hands and roll the rim off the edge of the chestpiece. Remove the diaphragm from the rim and clean the parts in soapy water or wipe with alcohol. Chestpiece surfaces can be wiped with alcohol or soapy water. Dry all parts and surfaces thoroughly before reassembly.


Assemble and replacement of the floating diaphragm

1. Using your fingertip, apply a small amount of talcum powder to the inside surface of the flexible edge of the diaphragm. This will facilitate assembly and maintain smooth and quiet low/high frequency alternation. Pat the diaphragm with fingers to remove excess talc.
2. Insert the flexible edge of the diaphragm into the groove of the rim. This is best accomplished by

starting with the rim positioned above the legible side of diaphragm.

3. Visually examine the ring to ensure that the flexible edge is smoothly engaged inside the rim. If necessary, flex the assembly by pinching the outer rim edge between the thumb and fingers. Repeat this flexing procedure after rotating the assembly one quarter.
4. To attach the diaphragm assembly to the chestpiece, engage the groove of the rim [with the diaphragm attached] around the chestpiece at one point and hold it in place with your thumbs.
5. Slowly roll the rim around and over the chestpiece edge using both thumbs, moving in opposite directions around the chestpiece.
6. Visually inspect the edge where the diaphragm engages the rim for uniform containment. Minor adjustments may be made by slightly pulling and rolling the rim away from the diaphragm, allowing the diaphragm to slip into position.

1. 2. 3. 4. 5. 6.

Cleaning

In order to avoid dirt accumulation, always store the stethoscope in its packaging. The outside of the instrument can be cleaned with a damp, soft and fluff-free fabric. For disinfection, the stethoscope can be cleaned with a fabric moistened with alcohol.

Do not use thermal sterilization or soaking.

Warning

Never bring the stethoscope into liquid and make sure that no liquid penetrates the housing! Since the stethoscope is not designed for operations, simple cleaning and/or disinfection with a surface disinfectant with an alcoholic base is sufficient. Eartips can be removed from the eartubes for thorough cleaning.

Avoid extreme heat, cold, solvents and oils.

Keep out of the reach of children

Use on intact skin

Dispose of in accordance with current regulations

Packaging and storage

A robust vibration-proof packaging can prevent the product from compacting outside.

No special storage requirements.

REF	Product code		Consult instructions for use		Caution: read instructions (warnings) carefully
LOT	Lot number		Keep in a cool, dry place		Medical Device compliant with Regulation (EU) 2017/745
MD	Medical Device		Keep away from sunlight		Manufacturer
UDI	Unique identifier		Authorized representative in the United Kingdom		Manufacturing date
CH REP	Authorized Representative in Switzerland				

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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