
ELETTRODI ELECTRODES

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Self-adhesive electrode
pads for electrical stimulation

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GB - Self-adhesive electrode is intended for use as a reusable, conductive adhesive interface between the patient's skin and the marketed electrical stimulators, i.e. TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscular Stimulation) IF or IFC (Interferential Current) or PGF (Pulsed Galvanic Stimulation) for transmitting electrical current, for OTC (Over-The-Counter) or prescription use.
Patient population: the self-adhesive electrode pads are for adults only

Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 80, 20537 Hamburg, Germany

Shenzhen Roundwhale Technology Co., Ltd.
202,2/F, Building 27, Dafa Industrial park,
Longxi Community, Longgang District, Shenzhen,
China - Made in China

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Application

1. Clean skin thoroughly prior to each application of electrodes. Electrodes do not stick well if any lotions, oil, make-up, dirt, etc, is left on the skin.
2. Turn the stimulator off before you apply the electrodes to your skin.
3. Insert the pins of the lead wire from the device into the electrode wire connectors.
4. Remove the electrodes from the protective liner and apply firmly to the skin.

Removal

1. Turn the stimulator off before you remove the electrodes from the skin.
2. Lifting a corner of electrode, remove the electrode from protective liner, peel it away. Do not peel the electrode by wire.
3. Place the electrodes on the liner and disconnect the lead wire by twisting and pulling at the electrode wire connectors.

Care and storage

1. Between uses, store the electrodes in the resealable bag in a cool, dry place.
2. The life of the electrodes varies depending on skin conditions, storage, amount of use, type of stimulation and stimulation site. Electrode life may be extended by carefully following the application instructions above and below.
3. Avoid electrodes from being exposed to direct sunlight or extreme humidity. Recommended storage conditions: 5°C-27°C. Relative humidity: 40%-70%.
4. DO NOT leave the electrodes uncovered in air as this will decrease the viscosity and performance of the product.

Caution

1. Read instructions carefully prior to use.
2. The electrodes are intended for single patient use only.
3. DO NOT apply to broken skin. Should skin rash occur, discontinue use and contact your physician, DO NOT use while driving a motorized vehicle or while operating machinery.
4. DO NOT disassemble or alter electrodes.
5. For use with electrical stimulation type devices only.
6. DO NOT use on children or infants. These devices have not been evaluated for pediatric use.

UDI

LATEX
FREE

LOT