LiceTreat

• **PETTINE MANUALE
PER PIDOCCHI**

• **MANUAL LICE COMB**

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• **PEINE PARA
PIOJOS MANUAL**

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• **XEIPOKINHTO XTENA
ΨΕΙΡΩΝ**

• **PENTE MANUAL
PARA PIOLHOS**

• **PIPTENE MANUAL
PENTRU PADUCHI**

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UTILIZATION: The long, delicately closed teeth of the lice comb are especially designed for detaching the lice and their nits.

INSTRUCTION: Disentangle the hair before the usage for facilitating the passage of the special teeth. Your hair should be dry. • Little by little slowly comb your hair from the root to the end, pay attention to every spot of host for the lice or nit. • The built in magnifier and the cleaning brush on the cover cap help mousing out the smallest lice and verify a great carding. • Repeat the operation several times if necessary, notably during the period of infection or for the prevention.

CLEANING AND DISINFECTION: Use the brush in order to carefully clean the teeth of the lice comb after the usage. It's recommended to disinfect the metallic part with 70% alcohol after each use.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Gima S.p.A - Via Marconi, 1