# Tegaderm™ High Performance Foam

## DESCRIZIONE

**Novità  
2014!**

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## FORMATI DISPONIBILI

**Formato sacrale**

**90616 large oval**  
 Size: 7 1/2 in x 8 3/4 in / 19 cm x 22,2 cm  
 Foam: 5 1/2 in x 6 3/4 in / 14 cm x 17,1 cm

**90613 medium oval**  
 Size: 5 1/2 in x 6 1/2 in / 14,3 cm x 15,6 cm  
 Foam: 4 in x 4 1/2 in / 10 cm x 11 cm

**90614 mini oval**  
 Size: 2 3/4 in x 3 in / 6,9 cm x 7,6 cm  
 Foam: 1 1/4 in x 1 1/2 in / 3,1 cm x 3,8 cm

**90611 small oval**  
 Size: 4 in x 4 1/2 in / 10 cm x 11 cm  
 Foam: 2 1/2 in x 3 in / 6 cm x 7,6 cm

**90612 square**  
 Size: 5 1/2 in x 5 1/2 in / 14,3 cm x 14,3 cm  
 Foam: 4 in x 4 in / 10 cm x 10 cm

**90610 square**  
 Size: 3 1/2 in x 3 1/2 in / 8,8 cm x 8,8 cm  
 Foam: 2 in x 2 in / 5 cm x 5 cm

**90619 heel/elbow design**  
 Size: 5 1/2 in x 5 1/2 in / 13,97 cm x 13,97 cm  
 Foam: 3 in x 3 in / 7,62 cm x 7,62 cm

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




# Tegaderm™ High Performance Foam

## BIBLIOGRAFIA

**3M™ Tegaderm™ Foam Adhesive Dressing  
in the treatment of leg skin ulcers**

---

M. BUCALOSSI <sup>1,2</sup>, S. SERANTONI <sup>1,3</sup>, L. GAZZABIN <sup>1,2</sup>, G. MOSTI <sup>1,4</sup>, F. MARIANI <sup>1,2</sup>

---

**Aim.** 3M™ Tegaderm™ Foam Adhesive Dressing is indicated in the treatment of leg skin ulcers. The aim of this study was to evaluate the effectiveness and tolerability of this dressing in a group of patients with leg ulcers of various aetiology, often recalcitrant to traditional topical treatment.

**Methods.** This prospective multicentric trial evaluated 73 ulcers of differing aetiology and degree of exudate (61 patients). The first dressings were performed weekly and involved surgical toilet and the laying of 3M™ Tegaderm™ Foam Adhesive Dressing on the lesion.

**Results.** Sixty-three ulcers healed within nine months from the beginning of therapy, of whom one third within one month, while 10 were still in treatment. Each patient showed good treatment compliance and there were no side effects.

**Conclusions.** Our study demonstrated that 3M™ Tegaderm™ Foam Adhesive Dressing is an optimal dressing for the treatment of all types of leg ulcers, assuring rapid healing of all skin lesions treated and an excellent quality of life for the patient.

**Key words:** Foam - Ulcer - Varicose ulcer.

**L**eg skin ulcers are the result of a variety of pathological processes and have become more and more relevant by virtue of their frequency and the treatment problems involved.

*Conflict of interest.*—None.  
Received on November 15, 2010.  
Accepted for publication on December 14, 2010.

Corresponding author: M. Bucalossi, Via Camilla 52, 53100 Siena, Italia. E-mail: matteo.bucalossi@aiuc.it

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The condition affects about 1% of the adult population in the developed nations (with more than 3% in the over-70s) and its nature makes it a veritable social scourge as it does not just involve the patient directly concerned but the whole family entourage.<sup>1</sup> Reported figures, however, are certainly underestimated. It is not infrequent that the first therapeutic approach to be attempted by the patient himself, taking advantage of the experience of friends or neighbours. Where there is no result or the clinical picture worsens, the family doctor is called in and only later is the patient referred to a skin ulcer specialist. It has been recognised that certain factors influence the management of ulcer treatment: the type of dressing, the frequency with which the dressing is changed, its tolerance, the time before the ulcer clears up and the prevention of recurrences.<sup>2</sup> The treatments proposed today are very often complex but their common characteristic is that they produce alternating or indeed unsatisfactory results.

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## BIBLIOGRAFIA

### CHRONIC WOUNDS

<p>Tegaderm™ Adhesive Foam Dressings — VS — Coloplast Biatain® Foam Adhesive Dressing — VS — Smith &amp; Nephew ALLEVYN™ Adhesive Dressing — VS — Medline Optifoam® Adhesive Dressing — VS — Molnlycke Mepilex® Border Dressing — VS — ConvaTec Versiva® XC® Adhesive Dressing — VS — Smith &amp; Nephew ALLEVYN™ Gentle Border Dressing</p>	<p><i>in vivo</i> and <i>in vitro</i> Studies  Surgical Wounds</p>	<p><i>In vitro</i> and Healthy Human Studies to Assess Foam Adhesive Dressing Breathability and Fluid Handling Properties</p>	<p><b>Tegaderm™ Adhesive Foam Dressing demonstrated superior fluid handling, out-performing other leading adhesive foam dressings and provided significantly longer wear times.</b></p> 	<p>Holm, D. Zehrer, C. Walters, S. A.</p>	<p>Presented at 2011 SAWC Spring</p>
<p>Tegaderm™ Adhesive Foam Dressing</p>	<p>Case Series</p>	<p>Use of a Foam Adhesive Dressing on Chronic Wounds</p>	<p><b>Tegaderm™ Adhesive Foam Dressing demonstrated successful clinical outcomes on three distinct wounds.</b> The dressing was easy to use, conformable and comfortable for the patient, and maintained a moist healing environment even on hard to dress and hard to heal wounds.</p>	<p>Brown-Etris, M. et al.</p>	<p>Presented at 2011 SAWC Spring and WOCN</p>
<p>Tegaderm™ Adhesive Foam Dressing</p>	<p>Case Study</p>	<p>Healing an Intractable Wound using a High Performance Foam Dressing</p>	<p><b>Tegaderm™ Adhesive Foam Dressing was used successfully to manage exudate and the wound environment on a chronic, hard to heal wound, allowing the wound to progress to healing and complete closure.</b> The dressing managed wound exudate effectively and reduced dressing change frequency.</p>	<p>Verrall, D. Kerr, A.</p>	<p>Poster presented at 2011 EWMA (in Brussels)</p>

### COST-EFFECTIVE

<p>Tegaderm™ Adhesive Foam Dressing — VS — Coloplast Biatain® Foam Adhesive Dressing — VS — Smith &amp; Nephew ALLEVYN™ Adhesive Dressing — VS — Medline Optifoam® Adhesive Dressing — VS — Molnlycke Mepilex® Border Dressing — VS — ConvaTec Versiva® XC® Adhesive Dressing</p>	<p>Health Economic Study</p>	<p>Impact on Cost: Adhesive Foam Dressing Comparisons in Home Care, Long-Term Care, and Hospital Settings</p>	<p>Tegaderm™ Adhesive Foam Dressing was compared to multiple foam dressings for wear time in order to develop an economic model based on material and labor costs in hospital, long-term care and home care settings. <b>In each care setting, Tegaderm™ Adhesive Foam Dressing was found to be cost effective.</b></p> 	<p>Carlson, A. Loegering, C. Zehrer, C. Walters, S. A.</p>	<p>Poster presented at 2011 SAWC Spring</p>
---	------------------------------	---	---	--	---

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