



GIMA

PROFESSIONAL MEDICAL PRODUCTS

ALTIMETRO DA PARETE

WALL MOUNTED HEIGHT METER

TOISE MURALE

ALTÍMETRO DE PARED

ALTÍMETRO DE PAREDE

WANDBEFESTIGTES ROLLMESSBAND

NAŚCIENNY WYSOKOŚCIOMIERZ

MĚSURÁTOR ÎNĂLȚIME MONTAT PE PERETE

FALRA SZERELHETŐ MAGASSÁGMÉRŐ

SEINÄASENTEINEN PITUUSMITTA

ΕΠΙΤΟΙΧΙΟ ΑΝΑΣΤΗΜΟΜΕΤΡΟ

مقياس طول حائطي التثبيت

РЪСТОМЕР, МОНТИРАН НА СТЕНА

SEINALE PAIGALDATAV PIKKUSMÕÕDIK

MJERITELJ VISINE MONTIRAN NA ZID

ANT SIENOS MONTUOJAMAS AUKŠČIO MATUOKLIS

WANDGEMONTEERDE LENGTEMETER

VÝŠKOMER MONTOVANÝ NA STENU

STENSKI MERILNIK VIŠINE

VÄGGMONTERAD LÄNGDMÄTARE

Incertezza di 9 mm a 2 metri

Uncertainty of 9 mm at 2 m

Au-delà de 2 mètres la valeur peut osciller de 9 mm

Oscilación de 9 mm a una distancia de 2 metros

Incerteza de medição (desvio-padrão) de 9 mm a 2 metros

Schwankung von 9 mm bei einer Distanz von 2 m

Niedokładność 9 mm na 2 m

Diferență de 9 mm la 2 m

2 méteren 9 mm-es pontatlanság előfordulhat

Epätarkkuus 9 mm 2 metrissä

Αβεβαιότητα 9 mm στα 2 μέτρα

نسبة الشك في القياس 9 ملم لكل 2 متر

Неточность 9 мм на 2 м

Mõõtemääramatus 2 m-l 9 mm

Odstupanje od 9 mm na 2m

9 mm paklaida 2m

Tolerantie 9 mm bij 2 m

Neistota 9 mm v 2 m

Odstopanje 9 mm na 2 m

Osäkerhet på 9 mm vid 2m

REF 27335



Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Made in China

CE 0476



FEATURES

GIMA wall mounted height meter (cod. 27335) is an instrument to measure a person's height. It can be secured to the wall and measure up to 2 m. It is space-saving and allows a quick reading.

INSTRUCTIONS



Do not pull the slider abruptly, as the graduated tape could break.

The product should be used by qualified personnel only.



Do not use the equipment if damaged. Refer back to your retailer. Avoid poor repairs. Repairs shall be carried out with original spare parts only, which should be installed according to the intended use.

Since the product is made of corrosion-proof materials suitable for environmental conditions appropriate to its normal use, it does not require special care.

However it must be stored in airtight conditions, making sure that it is protected from dust and dirt to maintain hygiene. Moreover, it is recommended that the product is stored in a place which can be easily reached by personnel in case of necessity.

UNPACKING



Always remember that packaging elements (paper, cellophane, staples, adhesive tape, etc.) can cut and/or injure if they are not handled carefully. They shall be removed with appropriate means by responsible personnel; the same applies to the tools used to remove packaging (scissors, knives, etc.).

After opening the packaging, first check all pieces and parts composing the product. Check that they are all present and in perfect

condition.

INSTALLATION

In order to erect the wall mounted height meter, the cursor should be positioned on the wall at a height of 2 metres exactly from the floor surface and, keeping the product compact, fix the small plate using the pins provided. Now the wall mounted height meter is ready for use.

OPERATION

In order to take a height measurement simply take the cursor with one hand and pull it downwards until it touches the head of the person to be measured, then read the measurement in line with the red bar in the window.

MAINTENANCE

The product has been designed to be durable. We recommend that the instrument is cleaned with a light detergent.

 Periodically clean the wall mounted height meter tape with a cloth dampened with disinfectant. Pay special attention to the slider.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

	<p>IT - Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB - Caution: read instructions (warnings) carefully FR - Attention: lisez attentivement les instructions (avertissements) ES - Precaución: lea las instrucciones (advertencias) cuidadosamente PT - Cuidado: leia as instruções (avisos) cuidadosamente DE - Achtung: Anweisungen (Warnungen) sorgfältig lesen PL - Ostrzeżenie - Zobacz instrukcję obsługi RO - Atenție: Citiți și respectați cu atenție instrucțiunile (avertismentele) de utilizare HU - Figyelem: Figyelmesen olvassa el és kövesse a használati utasításokat (figyelmeztetéseket) FI - Huomio: Lue käyttöohjeet (varoituskeset) ja noudata niitä huolellisesti GR - Προσοχή: Διαβάστε προσεκτικά τις οδηγίες (ενστάσεις) AR - الحذر - قراءة التعليمات (التحذيرات) بعناية</p>	<p>IT - Fabbricante GB - Manufacturer FR - Fabricant ES - Fabricante PT - Fabricante DE - Hersteller PL - Producent RO - Producător HU - Gyártó FI - Valmistaja GR - Παραγωγός AR - الشركة المصنعة</p>
	<p>IT - Leggere le istruzioni per l'uso GB - Consult instructions for use FR - Consulter les instructions d'utilisation ES - Consultar las instrucciones de uso PT - Consulte as instruções de uso DE - Gebrauchsanweisung beachten PL - Przeczytaj instrukcję użytkowania RO - Citiți instrucțiunile de utilizare HU - Olvassa el a használati utasításokat FI - Lue käyttöohjeet GR - Διαβάστε προσεκτικά τις οδηγίες χρήσης AR - اقرأ بدقة وحرص تعليمات الاستخدام</p>	<p>IT - Conservare in luogo fresco ed asciutto GB - Keep in a cool, dry place FR - À conserver dans un endroit frais et sec ES - Conservar en un lugar fresco y seco PT - Armazenar em local fresco e seco DE - An einem kühlen und trockenen Ort lagern PL - Przechowywać w suchym miejscu RO - A se păstra într-un loc răcoros și uscat HU - Száraz, hűvös helyen tárolandó FI - Säilytä kuivassa ja viileässä GR - Διατηρείται σε δροσερό και στεγνό περιβάλλον AR - يحفظ في مكان بارد وجاف</p> <p>IT - Conservare al riparo dalla luce solare GB - Keep away from sunlight FR - À conserver à l'abri de la lumière du soleil ES - Conservar al amparo de la luz solar PT - Guardar ao abrigo da luz solar DE - Vor Sonneneinstrahlung geschützt lagern PL - Przechowywać z dala od światła słonecznego RO - A se păstra ferit de razele soarelui HU - Napfénytől védve tárolandó FI - Säilytä auringonvalolta suojassa GR - Κρατήστε το μακριά από ηλιακή ακτινοβολία AR - يحفظ بعيداً عن ضوء الشمس</p>
	<p>IT - Dispositivo medico conforme alla Direttiva 93/42/CEE GB - Medical Device complies with Directive 93/42/EEC FR - Dispositif médical conforme à la directive 93/42 / CEE ES - Dispositivo médico según a la Directiva 93/42 / CEE PT - Dispositivo médico em conformidade com a Diretiva 93/42/CEE DE - Medizinprodukt gemäß Richtlinie 93/42/CEE PL - Wyrób medyczny zgodny z dyrektywą 93/42 / CEE RO - Dispozitiv medical realizat în conformitate cu prevederile Directivei 93/42/CEE HU - Dispozitiv medical realizat în conformitate cu prevederile Directivei 93/42/CEE FI - Direktiivin 93/42/ETY mukainen lääkekinnällinen laite GR - Ιατρική συσκευή σύμφωνα με την οδηγία 93/42 / CEE AR - جهاز طبي يتوافق مع التوجيه CEE/93/42</p>	<p>IT - Rappresentante autorizzato in Svizzera GB - Authorized Representative in Swiss FR - Représentant autorisé en Suisse ES - Representante autorizado en Suiza PT - Representante autorizado no Suíço DE - Autorisierter Vertreter im Schweizer PL - Autoryzowany Przedstawiciel w Szwajcarski FI - Valtuutettu edustaja Sveitsiläinen RO - Reprezentant autorizat în Elveția HU - Meghatalmazott képviselő az Svájc GR - Εξουσιοδοτημένος Αντιπρόσωπος στο ΕΛΒΕΤΟΣ AR - الممثل المعتمد في سويسرا</p>
	<p>IT - Codice prodotto GB - Product code FR - Code produit ES - Código producto PT - Código produto DE - Erzeugniscode PL - Numer katalogowy RO - Cod produs HU - Termékkód FI - Tuotekoodi GR - Κωδικός προϊόντος AR - كود المنتج</p>	<p>CH REP</p>
	<p>IT - Numero di lotto GB - Lot number FR - Numéro de lot ES - Número de lote PT - Número de lote DE - Chargennummer PL - Kod partii RO - Număr de lot HU - Tételszám FI - Eränumero GR - Αριθμός παρτίδας AR - رقم الدفعة</p>	