

ENGLISH

SIDE GUARDS

WARNINGS

Before performing any operation, visually check the integrity of the device. If the product shows any anomalies, immediately contact the point of sale where the purchase was made. Erregi sas assumes no responsibility for the improper use of the device. In the event of a serious incident related to the use of the device, please report the incident to the manufacturer and the competent authority of the member state where the user and/or patient is located.

DESCRIPTION

Side guards available in multiple versions and with various fastening systems, designed to ensure maximum adaptability to specific requirements. Manufactured from rubberized, water-repellent and fire-retardant fabrics, with polyurethane foam padding.

INTENDED USE

Devices used to prevent accidental impacts against the rigid side rails of the bed and to ensure the safety of adult patients with motor and control impairments.

INDICATIONS AND PRECAUTIONS FOR USE

Always use the device over a garment to protect the patient's skin.

CONTRAINDICATIONS

No contraindications have been detected, except in cases of proven sensitivity to the component materials.

COMPOSITION

OUTER FABRIC: Rubberized polyester/polyurethane, water-repellent and flame-resistant class 1

FILLING: Expanded polyurethane

WASHING INSTRUCTIONS

OUTER FABRIC:

			
WASH AT MAX 95°C	DO NOT BLEACH	DO NOT IRON	TUMBLE DRY AT LOW TEMPERATURES
	USE ANY SOLVENT EXCEPT TETRACHLORETHYLENE		

FILLING:

CLEAN WITH A SOFT DAMP CLOTH.

The manufacturer is not responsible for any changes in the technical specifications of the device in case of non-compliance with the washing instructions.

DISINFECTION

Use products for medical-surgical use.

STORAGE

Store in a place protected from light, heat, and moisture.

DISPOSAL

Dispose of the device in accordance with the applicable environmental protection regulations and separate waste collection.

	Product code		Medical Device
	Lot number		Manufacturer
	Medical Device compliant with Regulation (EU) 2017/745		Consult instructions for use
	Keep away from sunlight		Keep in a cool, dry place
	Date of manufacture		Not made with natural rubber latex
	Authorized Representative in Switzerland		Distributed by
	Unique Device Identifier		

DO NOT remove the label from the Medical Device as it ensures traceability of the supplied product (Manufacturer, UDI, batch number, product code).

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Gima 27638



Erregi S.a.s.

Via Adua 40
21040 Vedano Olona (Va)
commerciale@erregiausilimedici.it
www.erregiausilimedici.it
MADE IN ITALY



18/24



Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in Italy



Élite Med Group Sagi

Via Luigi Lavizzari 4
6850 Mendrisio
Switzerland

