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EN LATEX PROBE COVERS FOR VAGINAL AND ANAL ULTRASOUND SCANNING
INDICATIONS

For the prevention of the spread of infectious microorganisms to successive patients due to cross contamination during a vaginal or anal ultrasound.

USERS
Medical personnel with specific training in ultrasound scanning.

TARGET PATIENTS
Patients of all ages, sex, and condition, who undergo an ultrasound scan to examine, monitor or diagnose a disease in the anorectal area. Patients with female sexual organs who undergo ultrasound scan to examine, monitor, or diagnose a disease in the uterus, ovaries, fallopian tubes, cervix, and/or pelvic area. Patients with female sexual organs who undergo an ultrasound scan to monitor a pregnancy.

INSTRUCTIONS FOR USE
Dispense the coupling gel for ultrasound scanner inside the probe cover. Unroll the probe cover down with care over the probe or transducer of ultrasound scanner.

WARNINGS
The use of the latex cover is complementary and does not exempt from the rest of the measures for hygiene and

infection prevention during the ultrasound scan. Avoid handling the latex probe cover with fingernails, jewellery, sharp objects, etc. because they can damage the probe cover. Handle carefully when you remove it from its package.

If the individual wrapping is obviously damaged, discard the probe. If the latex probe cover is brittle, sticky or shows any visible signs of damage, do not use it. Check the expiry date printed on the package.

The latex probe cover is for one-use only. If the probe cover is used more than once, it increases the risk of tearing and being damaged and the infection of patient. Never wash it. Use only water-based ultrasound gels. Do not use oil-based gels, they may damage the latex cover.

Keep the probe cover in a cool dry place away from sunlight and heat. Do not take the probe cover out of its package until immediately prior to use. Do not use the probe cover if it has not been stored adequately. Discard the used probe cover in a container for refuse of medical devices.

Any serious incident related to the product must be reported to the manufacturer and the competent authority in your country.

CONTRA-INDICATIONS
This product is made from natural rubber latex. Latex may cause allergic reactions, including anaphylactic shock.

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