ENGLISH

Instructions/Instructions for Use for Reusable NON-SURGICAL INSTRUMENTS & Accessories

Introduction: These Instructions apply to Reusable NON-SURGICAL instruments and accessories supplied by S. Jee Enterprises

WARNINGS

Visually inspect for damage or deterioration before each use - Wear Gloves while using the instruments. Be careful of the sharp point and of the moving parts of the instruments, if any. The instruments must be stored in dry place and must be handled with care to prevent damage. Use only for its intended purpose to ensure maximum effectiveness. Steam (moist heat) is the only recommended method. Personal Protective Equipment (PPE) should be worn while using the instruments. Saline and cleaning/disinfection agents containing aldehyde, chloride, active chlorine, bromine, bromide, iodine or iodide are corrosive and should not be used. Do not allow biologic soil to dry on contaminated devices. Metal brushes and scouring pads must not be used during manual cleaning, because they will damage the surface and finish of the instruments. Use only soft bristle nylon brushes with different shapes, lengths and sizes to aid with manual cleaning. Do not process instruments with polymer components at temperatures equal to or greater than 140°C/285°F. Oils or silicone lubricants should not be used on instruments. Dry instruments with a clean, absorbent non-shedding lint free cloth. Clean, filtered compressed air may be used to remove moisture from lumens, holes, cannulas and difficult to access areas. After cleaning, all devices should be thoroughly inspected for residue biologic soil or detergent.

If contamination is still present repeat the cleaning process. Only qualified person should use the healthcare devices.

Although there is no chance of any serious incident due to these medical devices, anyway if you find any serious incident that has occurred in relation to these devices inform us (manufacturer) and relevant competent/regulatory authority.

Cleaning Steps

Place instruments in a suitable validated washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the instruments for maximum cleaning exposure; e.g. place concave instruments on their side or upside down, place heavy devices on the bottom of trays and baskets. If the washer/disinfector is equipped with special racks use them according to the manufacturer's instructions.

Recommended parameters are:

Sterilization

Always consult and follow the sterilizer manufacturer instructions for load configuration and equipment operation. Sterilizing equipment should have demonstrated efficacy (e.g. FDA clearance, EN 13060 or EN 285 compliance). Additionally the manufacturer's recommendations for installation, validation, and maintenance should be followed. Validated exposure times and temperatures to achieve a 10 sterility assurance level (SAL) are listed in the following table.

Cycle	Description
1	Pre-wash -- Cold Softened Tap Water -- 2 Minutes
2	Enzyme Spray & Soak -- Hot Softened Tap Water -- 1 Minute
3	Rinse -- Cold Softened Tap Water
4	Detergent Wash -- Hot Tap Water (64-66°C/146-150°F) -- 2 Minutes
5	Rinse -- Hot Purified Water (64-66°C/146-150°F) -- 1 Minute
6	Hot Air Dry (116°C / 240°C) -- 7 to 30 Minutes

Notes:

- The washer/disinfector manufacturer's instructions should be followed.
- A washer/disinfector with demonstrated efficacy (e.g. FDA approval, validated to ISO 15883) should be used.
- Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfector.

	Cycle Type	Minimum Temperature	Minimum Exposure Time	Pressure
For Europe recommended parameters	Pre-vacuum / Vacuum Pulse	134°C/273°°F	3 Minutes	28 p.s.i.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Simboli - Symbols - Symboles - Symbole - Símbolos - Símbolos - Σύμβολα - Символи - Symboly - Symboler - Sümbolid - Symbolit - Simboli - Szimbólumok - Simboliai - Simboli - Symbolen - Symbolika - Simboluri - Simboli - Symboler

الرموز

REF	**IT** - Codice prodotto **GB** - Product code **FR** - Code produit **DE** - Code produit **ES** - Código de producto **PT** - Código do produto **GR** - Κωδικός προϊόντος **BG** - Код на продукта **CZ** - Kód produktu **DA** - Produktkode **FI** - Tuotekoodi **HR** - Šifra proizvoda **HU** - Termék kódja **LV** - Produkta kods **NL** - Productcode **PL** - Kod produktu **RO** - Cod produs **SL** - Šifra izdelka **SV** - Produktkod **SA** - رمز المنتج		**IT** - Conservare in luogo fresco ed asciutto **GB** - Store in a cool, dry place **FR** - Conserver dans un endroit frais et sec **DE** - An einem kühlen, trockenen Ort aufbewahren **ES** - Conservar en un lugar fresco y seco. **PT** - Armazene em local fresco e seco **GR** - Φυλάσσεται σε δροσερό, ξηρό μέρος **BG** - Съхранявайте на хладно и сухо място **CZ** - Skladujte na chladném a suchém místě **DA** - Opbevares på et køligt, tørt sted **HR** - Čuvati na hladnom i suhom mjestu **HU** - Tárolja hűvös, száraz helyen **LV** - Uzglabāt vēsā, sausā vietā **NL** - Op een koele, droge plaats bewaren **PL** - Przechowywać w chłodnym, suchym miejscu **RO** - A se pastra intr-un loc racoros si uscat **SL** - Hraniti na hladnem in suhem mestu **SV** - Förvara på en sval, torr plats **SA** - تخزينها في مكان بارد وجاف
UDI	**IT** - Identificatore unico dispositivo **GB** - Unique device identifier **FR** - Identifiant unique de l'appareil **DE** - Eindeutige Gerätekennung **ES** - Identificador de dispositivo único **PT** - Identificador exclusivo do dispositivo **GR** - Μοναδικό αναγνωριστικό συσκευής **BG** - Уникален идентификатор на устройството **CZ** - Jedinečný identifikátor zařízení **DA** - Unik enhedsidentifikator **FI** - Ainutlaatuinen laitetunniste **HR** - Jedinstveni identifikator uređaja **HU** - Egyedi eszközazonosító **LV** - Unikāls ierīces identifikators **NL** - Unieke apparaat-ID **PL** - Unikalny identyfikator urządzenia **RO** - Identificator unic de dispozitiv **SL** - Enolični identifikator naprave **SV** - Unik enhetsidentifierare **SA** - معرف الجهاز الفريد		**IT** - Conservare al riparo dalla luce solare **GB** - Store away from sunlight **FR** - Conserver à l'abri du soleil **DE** - Vor Sonnenlicht schützen **ES** - Almacenar lejos de la luz solar **PT** - Armazene longe da luz solar **GR** - Φυλάσσεται μακριά από το φως του ήλιου **BG** - Да се съхранява далеч от слънчева светлина **CZ** - Skladujte mimo dosah slunečního záření **DA** - Opbevares væk fra sollys **HR** - Čuvati dalje od sunčeve svjetlosti **HU** - Tárolja napfénytől távol **LV** - Glabāt prom no saules gaismas **NL** - Uit de buurt van zonlicht bewaren **PL** - Przechowywać z dala od światła słonecznego **RO** - A se păstra departe de lumina soarelui **SL** - Hraniti ločeno od sončne svetlobe **SV** - Förvaras åtskilt från solljus **SA** - يخزن بعيدا عن أشعة الشمس
LOT	**IT** - Numero di lotto **GB** - Batch number **FR** - Numéro de lot **DE** - Chargennummer **ES** - Número de lote **PT** - Número do lote **GR** - Αριθμός παρτίδας **BG** - Партиден номер **CZ** - Číslo šarže **DA** - Batchnummer **FI** - Erän numero **HR** - Broj serije **HU** - Tételszám **LV** - Partijas numurs **NL** - Batchnummer **PL** - Numer partii **RO** - Numărul lotului **SL** - Številka serije **SV** - Batchnummer **SA** - رقم الدفعة	CE	**IT** - Dispositivo medico conforme al regolamento (UE) 2017/745 **GB** - Medical device compliant with Regulation (EU) 2017/745 **FR** - Dispositif médical conforme au Règlement (UE) 2017/745 **DE** - Medizinprodukt gemäß Verordnung (EU) 2017/745 **ES** - Dispositivo médico que cumple con el Reglamento (UE) 2017/745 **PT** - Dispositivo médico em conformidade com o Regulamento (UE) 2017/745 **GR** - Ιατροτεχνολογικό προϊόν συμμορφούμενο με τον Κανονισμό (ΕΕ) 2017/745 **BG** - Медицинско изделие, отговарящо на Регламент (ЕС) 2017/745 **CZ** - Zdravotnický prostředek v souladu s nařízením (EU) 2017/745 **DA** - Medicinsk udstyr i overensstemmelse med forordning (EU) 2017/745 **HR** - Medicinski uređaj usklađen s Uredbom (EU) 2017/745 **HU** - Az (EU) 2017/745 rendeletnek megfelelő orvostechnikai eszköz **LV** - Medicīnas ierīce, kas atbilst Regulai (ES) 2017/745 **NL** - Medisch hulpmiddel dat voldoet aan Verordening (EU) 2017/745 **PL** - Wyrób medyczny zgodny z Rozporządzeniem (UE) 2017/745 **RO** - Dispozitiv medical conform Regulamentului (UE) 2017/745 **SL** - Medicinski pripomoček skladen z Uredbo (EU) 2017/745 **SV** - Medicinsk utrustning i enlighet med förordning (EU) 2017/745 **SA** - جهاز طبي متوافق مع اللائحة (الاتحاد الأوروبي) 2017/745
MD	**IT** - Dispositivo medico **GB** - Medical device **FR** - Dispositif médical **DE** - Medizinisches Gerät **ES** - Dispositivo medico **PT** - Dispositivo médico **GR** - Ιατρική συσκευή **BG** - Медицинско изделие **CZ** - Lékařské zařízení **DA** - Medicinsk udstyr **FI** - Lääketieteellinen laite **HR** - Medicinski uređaj **HU** - Orvosi eszköz **LV** - Medicīnas ierīce **NL** - Medisch apparaat **PL** - Urządzenie medyczne **RO** - Dispozitiv medical **SL** - Medicinski pripomoček **SV** - Medicinsk apparat **SA** - جهاز طبي		**IT** - Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso **GB** - Warning: Read and follow the instructions (warnings) for use carefully **FR** - Avertissement : Lisez et suivez attentivement les instructions (avertissements) d'utilisation. **DE** - Warnung: Lesen und befolgen Sie die Gebrauchsanweisungen (Warnhinweise) sorgfältig **ES** - Advertencia: Lea y siga atentamente las instrucciones (advertencias) de uso. **PT** - Aviso: Leia e siga as instruções (avisos) de uso com atenção **GR** - Προειδοποίηση: Διαβάστε και ακολουθήστε προσεκτικά τις οδηγίες (προειδοποιήσεις) χρήσης **BG** - Предупреждение: Прочетете и следвайте внимателно инструкциите (предупрежденията) за употреба **CZ** - Upozornění: Pečlivě si přečtěte a dodržujte pokyny (varování) k použití **DA** - Advarsel: Læs og følg brugsanvisningen (advarslerne) omhyggeligt **FI** - Varoitus: Lue käyttöohjeet (varoitukset) huolellisesti ja noudata niitä **HR** - Upozorenje: Pažljivo pročitajte i slijedite upute (upozorenja) za uporabu **HU** - Figyelmeztetés: Olvassa el és kövesse figyelmesen a használati utasításokat (figyelmeztetéseket). **LV** - Brīdinājums: uzmanīgi izlasiet un ievērojiet lietošanas instrukcijas (brīdinājumus). **NL** - Waarschuwing: Lees en volg de instructies (waarschuwingen) voor gebruik zorgvuldig **PL** - Ostrzeżenie: Przeczytaj uważnie instrukcje (ostrzeżenia) dotyczące użytkowania i postępuj zgodnie z nimi **RO** - Avertisment: Citiți și urmați cu atenție instrucțiunile (avertismente) de utilizare **SL** - Opozorilo: natančno preberite in upoštevajte navodila (opozorila) za uporabo **SV** - Varning: Läs och följ instruktionerna (varningarna) för användning noggrant **SA** - تحذير: اقرأ واتبع التعليمات (التحذيرات) للاستخدام بعناية
EU REP	**IT** - Rappresentante autorizzato nella Comunità europea **GB** - Authorized representative in the European Community **FR** - Représentant autorisé dans la Communauté Européenne **DE** - Autorisierter Vertreter in der Europäischen Gemeinschaft **ES** - Representante autorizado en la Comunidad Europea **PT** - Representante autorizado na Comunidade Europeia **GR** - Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα **BG** - Упълномощен представител в Европейската общност **CZ** - Autorizovaný zástupce v Evropském společenství **DA** - Autoriseret repræsentant i Det Europæiske Fællesskab **FI** - Valtuutettu edustaja Euroopan yhteisössä **HR** - Ovlašteni predstavnik u Europskoj zajednici **HU** - Meghatalmazott képviselő az Európai Közösségben **LV** - Pilnvarotais pārstāvis Eiropas Kopienā **NL** - Gemachtigde vertegenwoordiger in de Europese Gemeenschap **PL** - Autoryzowany przedstawiciel na terenie Wspólnoty Europejskiej **RO** - Reprezentant autorizat în Comunitatea Europeană **SL** - Pooblaščeni zastopnik v Evropski skupnosti **SV** - Auktoriserad representant i Europeiska gemenskapen **SA** - الممثل المعتمد في المجموعة الأوروبية		**IT** - Non sterile **GB** - Not sterile **FR** - Non stérile **DE** - Nicht steril **ES** - No estéril **PT** - Não estéril **GR** - Όχι στείρα **BG** - Не е стерилен **CZ** - Ne sterilní **DA** - Ikke steril **FI** - Ei steriili **HR** - Nije sterilan **HU** - Nem steril **LV** - Nav sterils **NL** - Niet steriel **PL** - Nie sterylne **RO** - Nu este steril **SL** - Ni sterilno **SV** - Inte steril **SA** - غير معقمة
	IT - Importato da **GB** - Imported by **FR** - Importé par **DE** - Importiert von **ES** - Importado por **PT** - Importado por **GR** - Εισαγωγή από **BG** - Внесено от **CZ** - Dovezeno uživatelem **DA** - Importeret af **FI** - Tuoja **HR** - Uvezeno od strane **HU** - Importálta: **LV** - Importēja **NL** - Geïmporteerd door **PL** - Importowane przez **RO** - Importat de **SL** - Uvozil **SV** - Importerad av **SA** - مستورد بواسطة		**IT** - Data di fabbricazione **GB** - Manufacturing date **FR** - Date de fabrication **DE** - Herstellungsdatum **ES** - Fecha de fabricación **PT** - Data de fabricação **GR** - Ημερομηνία κατασκευής **BG** - Дата на производство **CZ** - Datum výroby **DA** - Fremstillingsdato **FI** - Valmistuspäivämäärä **HR** - Datum proizvodnje **HU** - Gyártási dátum **LV** - Ražošanas datums **NL** - Productiedatum **PL** - Data produkcji **RO** - Data fabricatiei **SL** - Datum izdelave **SV** - Tillverkningsdatum **SA** - تاريخ التصنيع
	IT - Leggere le istruzioni per l'uso **GB** - Read the instructions for use **FR** - Lire la notice d'utilisation **DE** - Lesen Sie die Gebrauchsanweisung **ES** - Lea las instrucciones de uso. **PT** - Leia as instruções de uso **GR** - Διαβάστε τις οδηγίες χρήσης **BG** - Прочетете инструкциите за употреба **CZ** - Přečtěte si návod k použití **DA** - Læs brugsanvisningen **FI** - Lue käyttöohjeet **HR** - Pročitajte upute za uporabu **HU** - Olvassa el a használati utasítást **LV** - Izlasiet lietošanas instrukcijas **NL** - Lees de gebruiksaanwijzing **PL** - Przeczytaj instrukcję użytkowania **RO** - Citiți instrucțiunile de utilizare **SL** - Preberite navodila za uporabo **SV** - Läs bruksanvisningen **SA** - اقرأ تعليمات الاستخدام		**IT** - Limite di umidità **GB** - Humidity limit **FR** - Limite d'humidité **DE** - Luftfeuchtigkeitsgrenze **ES** - Límite de humedad **PT** - Limite de umidade **GR** - Όριο υγρασίας **BG** - Ограничение на влажността **CZ** - Limit vlhkosti **DA** - Fugtgrænse **FI** - Kosteusraja **HR** - Granica vlažnosti **HU** - Páratartalom határérték **LV** - Mitruma robeža **NL** - Vochtigheidslimiet **PL** - Limit wilgotności **RO** - Limita de umiditate **SL** - Omejitev vlažnosti **SV** - Fuktighetsgräns **SA** - حد الرطوبة
	IT - Fabbricante **GB** - Manufacturer **FR** - Fabricant **DE** - Hersteller **ES** - Fabricante **PT** - Fabricante **GR** - Κατασκευαστής **BG** - производител **CZ** - Výrobce **DA** - Fabrikant **FI** - Valmistaja **HR** - Proizvođač **HU** - Gyártó **LV** - Ražotājs **NL** - Fabrikant **PL** - Producent **RO** - Producător **SL** - Proizvajalec **SV** - Tillverkare **SA** - الشركة المصنعة		**IT** - Limite di temperatura **GB** - Temperature limit **FR** - Limite de température **DE** - Temperaturgrenze **ES** - Límite de temperatura **PT** - Limite de temperatura **GR** - Όριο θερμοκρασίας **BG** - Температурна граница **CZ** - Teplotní limit **DA** - Temperaturgrænse **FI** - Lämpötilaraja **HR** - Ograničenje temperature **HU** - Hőmérséklet korlát **LV** - Temperatūras ierobežojums **NL** - Temperatuurlimiet **PL** - Limit temperatury **RO** - Limită de temperatură **SL** - Temperaturna omejitev **SV** - Temperaturgräns **SA** - حد درجة الحرارة

S. Jee ENTERPRISES
CANTT ROAD, Sialkot-51310-Pakistan
Made in Pakistan

Gima S.p.A
Via Marconi,1 - 20060 Gessate (Mi) - Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

EU REP
KHANG BAO INTERNATIONAL KFT
Auróra utca, 13. al. 1,
Budapest – 1084 Hungary

M26653-M-Rev.2.08.26