GIMA

PROFESSIONAL MEDICAL PRODUCTS

# SISTEMA REGGIARTO
# LIFTING POLE FOR BED
# POTENCE POUR LIT
# SISTEMA INCORPORADOR DE TRAPECIO PARA CAMA
# PÓLO DE LEVANTAMENTO DA CAMA
# UMREIFUNGSSYSTEM
# SISTEM DE FIXARE CU BENZI
# BANDNINGSSYSTEM
# SYSTEM MOCOWANIA
# ΣΥΣΤΗΜΑ ΠΡΟΣΔΕΣΗΣ

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**GIMA 27664**

*M27664-M-Rev.1-10.26*

**Anhui Longway Medical Technology Co., Ltd.**
No. 59 Lingji Rd., Industrial Park, Mingguang, Anhui, China
Made in China

REF **LW10207**

EU REP **Lotus NL B.V.**
Koningin Julianaplein 10, 1 e Verd, 2595AA,
The Hague, Netherlands

UK REP **Lotus global co ltd**
23 Maine Street, Reading, RG2 6AG,
England, United Kingdom

**Gima S.p.A.**
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
**www.gimaitaly.com**

MD

***Weight Capacity: 75kgs***

## INTENDED USE:

Longway Medical Accessories are Medical Devices intended to compensate for a disability or inability of the bedridden user. Support device with ground base, adaptable to any bed. Structure in tubular steel painted that can be positioned according to the patient's needs and triangle with strap to adjust its height.

**CAUTION (!):**

• In case of damage or malfunction of your product DOT NOT use the product and please contact the authorized and trained representative.
• Do not use the product for a purpose not indicated in this manual
• LONGWAY Medical declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
• LONGWAY Medical reserves the right to change the information contained in this document without previous notice
• There is a ± 3%. Deviation for the Dimensions of the Products.

## GENERAL SAFETY WARNINGS:

⚠ • If you don't read this user manual, it is preferable not to use this product or another available part. If you don't understand the usage or the precautions, please contact with the dealer or the proper technical person before you use the product because may be caused damage.
⚠ • Keep away the product from heat sources.
⚠ • Save this document for future reference.
⚠ • Be careful when children are nearby and do not allow the children to play with the product.
⚠ • Do not exceed the maximum payload.
⚠ • Special care is always required where there are moving parts that can cause limb entrapment and injury.
⚠ • Do not attempt to lift the product by any removable parts.
⚠ • Do not use the product seat for a different purpose from the intended one.
⚠ • The lifting poles are not intended to carry the full weight of an individual.

## DECLARATION OF CONFORMITY:

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 745/2017 and where appropriate, the standards and legislation referred to.

## BEFORE USE:

• Check if the frame is damaged in order to guarantee a safe use of the product. (There are no cracks or fractures in the frame).
• In case of damage do not use the product and contact your dealer for further instructions.
• Check, if the product is properly assembled and all screws are secure and well screwed on.
• Always check the condition of wear for the mechanical parts to ensure the absolute safe use of the product for people and objects

## ENVIRONMENTAL PROTECTION:

If one day you find that your product needs to be replaced or it is no longer working for you, consider protecting the environment:
1) Do not dispose your product along with the rest of the public waste (this is also the meaning of the shown recycling sign).
2) Contact your Public Authorities and they will instruct you of the Recycling centers to which your product must be disposed.
3) Correct disposal of your product helps the protection of the environment as well as the recycling of the product's components.


## MAINTENANCE, CLEANING & DISINFECTION INFORMATION:

To guarantee a safe use and adequate standard of hygiene, the user should perform these procedures before every use.
The user has to make certain of the structural integrity of the device and its components. For cleaning and sanitation procedures, follow the steps below:
1. Gloves should be worn.
2. Remove the evident dirt before to proceed with the most accurate disinfection.
3. Use water or mild detergent for the cleaning procedure.
4. Wipe with a dry and clean cloth.
5. DO NOT leave the product wet.
6. Do not use chemical cleaners to clean the frame. This could cause damage to the surface of the product
7. If you need to disinfect the product, use a common, mild disinfectant
8. Protect the item from scratches, cuts and punctures.
9. Lubrication: the product has designed for a minimal maintenance; however, check every six months the lubrication on the moving & mechanical parts, this would ensure an optimal duration.

## ASSEMBLING:

NOTE: The illustrations for the assembly are indicative to make the assembly easier. May will be differences from reality.
The support with LW10207 base is delivered disassembled and requires the following assembly operations:
1. Use 2*M10*115mm screws to fix the straight tube between the left and right leg tubes;

2. Insert the lifting pole in the space provided, and use plastic screw to fixed;
3. Insert the triangle with adjustable strap into the space provided.

## WARRANTY:

The product is guaranteed for a period of one (1) year from the purchasing date.
The warranty apply on factory defects & DOES NOT apply on damages resulting from misuse, poor maintenance, modification, overuse, or non-compliance with the product's operating instructions.
Also DOES NOT apply on the parts, that can wear out during use or over time. Damages or defects caused by natural disasters, unauthorized maintenance or repair, power supply problems (where applicable) & carriers are not covered by the warranty.
No service or spare part is applied unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the end user- customer or by dealer.

## *SYMBOL LEGEND*

| | |
|---|---|
|  | Manufacturer |
|  | Date of manufacture |
| REF | Product code |
| LOT | Lot number |
| EU REP | Authorized representative in the European community |
|  | Caution: read instructions (warnings) carefully |
|  | Follow instructions for use |
| CE | Medical Device compliant with Regulation (EU) 2017/745 |
| MD | Medical device |
|  | Keep away from sunlight |
|  | Keep in a cool, dry place |
|  | Imported by |
| UDI | Unique Device identifier |
| UK REP | Authorised representative in the United Kingdom |

**GIMA WARRANTY TERMS**
The Gima 12-month standard B2B warranty applies.