GIMA

PROFESSIONAL MEDICAL PRODUCTS

TELO PORTAFERITI 190X75 CM
CARRYING COT 190X75 CM
CIVIÈRE DE TRANSPORT 190X75 CM
PORTA-HERIDOS 190X75 CM
MACA DE LONA PARA TRANSPORTAR FERIDOS 190X75 CM
TRAGETUCH FÜR VERLETZTE 190X75 CM

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ENGLISH

INDICATIONS: the flexible carrying cot moves patient securely in restricted areas where standard stretchers have difficulties to be used. **It should not be used if there is a spine, backbone or pelvis suspected injury.**
SPECIFICATIONS: strong polyester fabric, PVC coated, with four secured handles per side to enable a comfortable patient carriage. The materials are waterproof and are easily washable and disinfectable, bacterial and fungus proof, non putrescent. Transparent to X rays.
INSTRUCTIONS: **should be used by skilled staff**. After laying the patient on the carrying cot,
four rescuers (as per fig. 1) should pull up contemporary on the lifting handles and bring the patient to destination. However if four are not available, two rescuers (as per fig. 2).
CLEANING INSTRUCTIONS: after using the carrying sheet, it should be accurately washed with water and soap. It can also be disinfected by a bacterial liquid solution.
STORAGE: no special instructions are requested. Better if kept in a dry place, out of dust and direct sunlight.
MAXIMUM LOAD: 150 kg.
Made according EN 1865-1:2010+A1:2015

GIMA WARRANTY TERMS
The Gima 12-month standard B2B warranty applies.

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