



Disposable Dermal curette

- Stainless steel blades with PP handle
- Ribbed handle for sure grip and control
- Size identification on each curette
- Used for curettage of tissue, mostly skin
- The curettes are sterile and single use
- Clinical areas: Dermatology, Surgery, Veterinary and Podiatry

Sizes:

MK402	2.0mm		MK405	5.0mm	
MK403	3.0mm		MK407	7.0mm	
MK404	4.0mm				

Packaging: Each curette in a blister pack and 20 pcs. in a dispenser box, 20 boxes per export carton.

Safe handling

- Store any new products in a dry, clean and safe place
- Avoid high temperature, high humidity, direct sunlight, soaking and preserve at room temperature
- Do not preserve at the chemical products' place and gas occurrence place
- Preserve at the status of non-break to package material and non-pinhole
- It is necessary to confirm the use period on each package and box (Shelf life 5 years)
- By confirming the use period on each package and box, dispose the products which expire the use period.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged
- Avoid contact of the delicate blades at tips with hard surfaces or other instruments
- Re-use prohibited
- Re-sterile prohibited
- Do not use after expiry date
- Do not use if it is considered that the products aseptic condition is spoiled by broken package



kai industries co.ltd, 1110 Oyana, Seki City, Gifu Pref., 501-3992, Japan



kai Europe GmbH, Kottendorfer Straße 5, 42697 Solingen, Germany

CE0197

Product complies with
Requirements of Directive
93/42EEC