




GIMA

PROFESSIONAL MEDICAL PRODUCTS

GONIOMETRO DA 7" ROBINSON 180° 180° ROBINSON GONIOMETER 7" GONIOMÈTRE 7" 180° ROBINSON 180° GONIÓMETRO ROBINSON 180° - 7"

Manuale d'uso - User manual - Manuel de l'utilisateur - Guía de uso



Tipo di acciaio - Steel Type
Type d'acier - Tipo de acero:
AISI 430

Le misure si intendono indicative
The measures are intended as indicative
Les mesures sont destinées à être indicatives
Las medidas tienen carácter indicativo

GIMA 24902

Fabbricante, Manufacturer, Fabricante, Fabricante:
JIMCO INDUSTRIES
4-KM Aimen Abad Road Near Akbar Abad Chowk Bhadal,
Sialkot-51310-Punjab - Pakistan - info@jimcoind.com
Made in Pakistan

REF RE-210-10

Importato da / Imported by / Importé par/ Importado por:
Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



GB - INSTRUCTIONS OF USE

Goniometers measure range of motion of finger joints (metacarpophalangeal and interphalangeal) and other small joints including toes. A goniometer measures range-of-motion (ROM) of various joints.

The 180° head is marked in 5° increments.

Arm measures in inches and centimeters.

Tension between the arms is controlled by the thumb knob.

INDICATION OF USE & PRECAUTIONS

Those devices are designed for professional use by qualified operators with necessary skills.

The use of the product involves direct contact with patient's skin; always wear protective gloves, do not use on dull skin, handle carefully to avoid accidental abrasion.





Clean and sterilize before each use; use neutral non-aggressive disinfectant solution (PH7), do not use solutions containing chloride ions or high alcohol content.

This product is intended to provide easy-to-read indicative measurements.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

Simboli / Symbols / Symboles / Símbolos

	<p>IT Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB Caution: read instructions (warnings) carefully FR Attention: lisez attentivement les instructions (avertissements) ES Precaución: lea las instrucciones (advertencias) cuidadosamente</p>		<p>IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil ES Conservar al amparo de la luz solar</p>	<div style="border: 1px solid black; padding: 2px; text-align: center;">LOT</div>	<p>IT Numero di lotto GB Lot number FR Numéro de lot ES Número de lote</p>
	<p>IT Leggere le istruzioni per l'uso GB Consult instructions for use FR Consulter les instructions d'utilisation ES Consultar las instrucciones de uso</p>		<p>IT Conservare in luogo fresco e asciutto GB Keep in a cool, dry place FR À conserver dans un endroit frais et sec ES Conservar en un lugar fresco y seco</p>	<div style="border: 1px solid black; padding: 2px; text-align: center;">REF</div>	<p>IT Codice prodotto GB Product code FR Code produit ES Código producto</p>