




GIMA

PROFESSIONAL MEDICAL PRODUCTS

GONIOMETRO DA 360 GRADI 14" **14" 360 DEGREE GONIOMETER** **GONIOMÈTRE 360 DEGRÉS 14"** **GONIÓMETRO 360 GRADOS - 14"**

Manuale d'uso - User manual - Manuel de l'utilisateur - Guía de uso



Tipo di acciaio - Steel Type
Type d'acier - Tipo de acero:
AISI 430

Incrementi di 1°
1° increments
Incréments de 1 degré
Intervalos de 1°

Le misure si intendono indicative
The measures are intended as indicative
Les mesures sont destinées à être indicatives
Las medidas tienen carácter indicativo

GIMA 24903

Fabbricante, Manufacturer, Fabricante, Fabricante:
JIMCO INDUSTRIES
4-KM Aimen Abad Road Near Akbar Abad Chowk Bhalal,
Sialkot-51310-Punjab - Pakistan - info@jimcoind.com
Made in Pakistan

REF RE-230-10

Importato da / Imported by / Importé par/ Importado por:
Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



GB - INSTRUCTIONS OF USE

A goniometer measures range-of-motion (ROM) of various joints. The 360° head is marked in 1° increments. Tension between the arms is controlled by the thumb knob. Goniometer has two opposing 180° scales and one 360° scale, all marked in 1° increments. Sturdy adjustment knob, controlling tension of arms.







INDICATION OF USE & PRECAUTIONS

Those devices are designed for professional use by qualified operators with necessary skills. The use of the product involves direct contact with patient's skin; always wear protective gloves, do not use on dull skin, handle carefully to avoid accidental abrasion. Clean and sterilize before each use; use neutral non-aggressive disinfectant solution (PH7), do not use solutions containing chloride ions or high alcohol content.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

Simboli / Symbols / Symboles / Simbolos

	<p>IT Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB Caution: read instructions (warnings) carefully FR Attention: lisez attentivement les instructions (avertissements) ES Precaución: lea las instrucciones (advertencias) cuidadosamente</p>	 <p>IT Conservare al riparo dalla luce solare GB Keep away from sunlight FR À conserver à l'abri de la lumière du soleil ES Conservar al amparo de la luz solar</p>		<p>IT Numero di lotto GB Lot number FR Numéro de lot ES Número de lote</p>
	<p>IT Leggere le istruzioni per l'uso GB Consult instructions for use FR Consulter les instructions d'utilisation ES Consultar las instrucciones de uso</p>	 <p>IT Conservare in luogo fresco e asciutto GB Keep in a cool, dry place FR À conserver dans un endroit frais et sec ES Conservar en un lugar fresco y seco</p>		<p>IT Codice prodotto GB Product code FR Code produit ES Código producto</p>