


Universal Holder DRACULA

Conveniently keeps boxes of disposable gloves, tissues, or wipes within reach when mounted near work area. Perfect for doctor's offices, laboratories, industrial facilities or wherever efficiency counts. The spring-loaded clamp is adjustable to fit many box sizes.

Installation:

Universal Holder comes ready to install with 2 adhesive strips on back, for mounting on smooth non-porous surfaces. Two screw holes have been provided, if stronger installation is required.



ger installation is required.

To insert box, push top of holder upwards with box as shown. The "teeth" will clamp down on box and hold it securely in place.