

INSTRUCTIONS FOR USE

INFUSION SET

PURPOSE OF USE

Medical device for the administration of medicinal substances.

INSTRUCTIONS FOR USE:

Before use, check the expiry date.

Check the integrity of the packaging before use.

Open the package and use immediately.

Do not use if the packaging is damaged.

In case of reuse, infections may develop: consult a doctor.

Prepare the area carefully and aseptically.

Remove the device from the packaging.

Proceed with administration.

At the end of the process, dispose of as hospital waste.

LATEX FREE:

The device is latex free to prevent any risk of allergic reaction.

STERILIZATION METHOD:

Sterilized by Ethylene Oxide (EtO), residual amount of gas in accordance with the regulations.

DISPOSAL:

Disposal must be made in accordance with the applicable government regulations.

Incineration or landfill disposal must be carried out under controlled conditions.

The product that did not come in to contact with infectious agents or other hazardous materials, cannot be considered dangerous waste.

CONSERVATION:

Store in a cool dry place and far away from the any source of heat at a temperature.

These Instructions for Use note must be preserved by the user for the duration of the DM.



Health&Safety

RAYS SPA

via Francesco Crispi - 60027 Osimo (AN)

Tel. 071 2868935 - Fax 071 2868911 - www.rays.it - info@rays.it