



Disposable Biopsiblades

- Stainless steel blades with PP handle
- Sure grip handle
- Used for curettage of tissue, mostly skin
- The blades are sterile and single use
- Clinical areas: Dermatology, Gynecology, Podiatry

BB-01:




Packaging: Each blade is individually sealed in a sterile pouch and 20 pcs. in a dispenser box, 40 boxes per export carton

Safe handling

- Store any new products in a dry, clean and safe place
- Avoid high temperature, high humidity, direct sunlight, soaking and preserve at room temperature
- Do not preserve at the chemical products' place and gas occurrence place
- Do not place in or near any products containing chlorides as rusting may occur
- Preserve at the status of non-break to package material and non-pinhole
- It is necessary to confirm the use period on each package and box (Shelf life 5 years)
- By confirming the use period on each package and box, dispose the products which expire the use period.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged
- Avoid contact of the delicate blades at tips with hard surfaces or other instruments
- Re-use prohibited
- Re-sterile prohibited
- Do not use after expiry date
- Do not use if it is considered that the products aseptic condition is spoiled by broken package



kai industries co.ltd, 1110 Oyana, Seki City, Gifu Pref., 501-3992, Japan



kai Europe GmbH, Kottendorfer Straße 5, 42697 Solingen, Germany

CE0197

Product complies with
Requirements of Directive
93/42EEC