



GIMA

PROFESSIONAL MEDICAL PRODUCTS

BORSE E ZAINI EMERGENZA
EMERGENCY BAGS AND BACKPACKS
SACS ET SACS À DOS D'URGENCE
BOLSOS Y MOCHILAS EMERGENCIA
ERSTE-HILFE-TASCHEN UND ERSTE-HILFERUCKSÄCKE
BOLSAS E MOCHILAS DE EMERGÊNCIA
ΤΣΑΝΤΕΣ ΚΑΙ ΣΑΚΙΔΙΑ ΕΠΕΙΓΟΥΣΑΣ ΑΝΑΓΚΗΣ
СПЕШНИ ЧАНТИ И РАНИЦИ
NOUZOVÉ TAŠKY A BATOHY
NØDTASKER OG RYGSAKKE
HÄDAABI KOTID JA SELJAKOTID
HÄTÄLAUKUT JA REPUT
TORBE I RUKSAK ZA HITNU POMÓČ
VÉSZTÁSKA ÉS HÁTIZSÁK
AVARINIAI KREPŠIAI IR KUPRINĖS
ÅRKÅRTAS SOMAS UN MUGURSOMAS
NØDSEKKER OG RYGGSEKER
NOODTASSEN EN RUGZAKKEN
TORBY I PLECAKI AWARYJNE
GENŢI ŞI RUCI DE URGENŢĂ
NÚDZOVÉ TAŠKY A BATOHY
TORBE IN NAHRBTNIKI ZA SILI
NÖDSÄCKAR OCH RYGGSÄCKAR
ئىراوطلار تالاحل رەھز جورخو بىئاقح

لصانع

Manufacturer/Fabbricante/Fabricant/Fabricante/Hersteller/Fabricante/Výrobce/Proizvodáč/
Κατασκευαστής/Gamintojas Fabrikant/Producent/Producent/Tilverkare/Gyártó/Valmistaja/
производитель/Fabrikant/Tootja/Ražotājs/Výrobca/Proizvajalec/Produsent/مصنع
QUANZHOU SKYSOAR BAGS CO.,LTD
No. 163 Beiqingdong Road, Beifeng Industrial Area, Quanzhou 362000, China
info@skysoars.cn

Importato da / Imported by / Importé de / Importado de / Importado por / Eingeführt von / Εισαγωγή
από / Importat de / Importowane przez / Importováno / Importerad av / Maahantuoja / Uvezao /
Importálja / Importované z / Importeret af / Uvozil / Importované z / Geïmporteerd door / Importert av /
Importavo / :
Importē / Внесен от / Imporditud asukohast / استه براد
Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China

FEATURES

GIMA bags and backpacks are designed to make first aid interventions more efficient. They assure safe transportation and the instruments, medications, products or medical and first aid devices are stored in practical compartments. Bags and backpacks are supplied empty

PRESCRIPTIONS

The product can be used only by qualified people, aware of the hygienic limitations due to the use of the product in places where there could be polluting or bacteriologically infected elements. GIMA cannot be held responsible for damages caused to the content or that the contents of the bag may cause to the bags or backpacks.

The product, or any of its parts, cannot be used for purposes other than those it was designed for.

Do not leave the product unattended, in the reach of children or unqualified personnel.

Do not use the product if it is damaged, contact dealer.

Avoid any temporary repair.

During use be careful not to damage the zippers.

When the zippers are being opened and closed verify that they do not get jammed into other materials and always open or close them completely.

If any liquids or powders should leak out of bottles, clean and disinfect the bag or backpack, because they could damage or contaminate the bag itself or the products inside



Do not fill up too much. The weight and shape of the content must fit the type of bag or backpack. Always verify that the weight of the content does not exceed the weight indications of the bag model or the possibility of the user. For any indications regarding maximum weights that can be lifted and moved manually, refer to the local regulations in force..

UNPACKING



The elements of the packaging (paper, cellophane, metal staples, adhesive tape, etc.) may cut and/or hurt if not handled with care. Remove them with appropriate tools and do not leave with unauthorized people. The same goes for the tools used to open the packaging (scissors, knives, etc.).

The first operation to do after opening the packaging is to control all the parts of the product; verify that all the components needed are in place and verify their conditions: movement of zippers, fastening of clips, Velcro or strips, resistance of handles.

MAINTENANCE

Emptying all the pockets and shaking out any loose debris. then, use a mild detergent and water to spot clean any stains or dirty areas. for a more thorough cleaning, you can hand wash the entire bag using a mild detergent and warm water. Gently scrub the fabric with a soft brush or cloth, and then rinse thoroughly. Finally, hang the bag to air dry. Avoid using harsh chemicals or bleach, as these can damage the polyester fabric

The bags and backpacks cannot be put in an autoclave

Comply with the laws in the country of use to properly dispose of the product.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.