



seca

seca 203





GB Guarantee

We guarantee defects attributable to poor materials or workmanship for 24 months from the date of purchase.

The guarantee is invalidated if the device is opened by unauthorised parties.

Important: guarantee valid only in conjunction with proof of purchase.



WHR	✓	⚠
♂	< 1.0 *	≥ 1.0 *
♀	< 0.85 *	≥ 0.85 *

*Obesity: preventing and managing the global epidemic. Report of a WHO Consultation. WHO technical report series 894 (1999: Geneva Switzerland).

Technische Daten • Technical Data • Fiche technique • Specifiche tecniche • Datos técnicos • Dane techniczne • Dados técnicos

MAX	205 cm
d	1 mm
	±5 mm

CE Declaration of Conformity

seca gmbh & co. kg hereby declares that the product meets the terms of the applicable European directives. The unabridged declaration of conformity can be found at: www.seca.com.

	$\varnothing = 70 \text{ mm}$ $h = 28 \text{ mm}$
	~ 75 g
	0.5 - 0.9: ±0.03 0.9 - 1.2: ±0.04 1.2 - 1.6: ±0.05 1.6 - 2.0: ±0.06
	10 °C - 40 °C

	-25 °C - 65 °C
--	----------------

Reinigung • Cleaning • Nettoyage • Pulizia Limpieza • Czyszczenie • Limpeza

Clean the measuring device as required using a domestic cleaning agent or commercially-available disinfectant. Follow the manufacturer's instructions.

Medical Measuring
Systems and Scales
since 1840

seca gmbh & co. kg
Hafenweg 33a/31a/31b/31c
22089 Hamburg - Germany
Telephone +49 40 20 00 00 00
Fax +49 40 20 00 00 50
info@seca.com

All contact data under www.seca.com

seca
group

seca deutschland seca russia seca suomi
seca france seca mexico seca america latina
seca united kingdom seca australia seca asia pacific
seca north america seca brazil
seca schweiz seca middle east
seca zhong guo seca india