

ENGLISH

VELCRO BEDSIDE PAD
WARNINGS

Before carrying out any operation, visually check the integrity of the device. If the product shows any abnormalities, immediately contact the point of sale. Erregi sas shall not be held responsible for the improper use of the device.

In the event of a serious accident involving the use of the device, please report the incident to the manufacturer and to the competent authority of the member state where the user and/or patient lives.

DESCRIPTION AND COMPOSITION

Art 27637: Side rail pad with velcro made of water-repellent material 100% polyester support class 1 flame-retardant acrylic/pu. PE ratio 75%, resin 25%. Polyurethane foam padding.

INDICATIONS






Device used to prevent accidental bumps against the rigid side rails of the bed and to ensure safety in people with motor and control difficulties, supplied in different types of fabrics and shapes.

CONTRAINDICATIONS

No contraindications were found except in cases of ascertained sensitivity to the component materials.

WASHING INSTRUCTIONS

ART 27637

			
WASH AT MAX 40°	DO NOT BLEACH	DO NOT IRON	DO NOT TUMBLE DRY
	USE ANY SOLVENT EXCEPT TETRACHLORETHYLENE		

PADDING:

WASH WITH A SOFT DAMP CLOTH

DISINFECTION











Use medical surgical products.

STORAGE

Store away from sources of light, heat and humidity. Avoid direct exposure to sunlight.

DISPOSAL

Dispose of the device in compliance with the current regulations on environmental protection and recycling.

	Product code		Medical Device
	Lot number		Manufacturer
	Medical Device compliant with Regulation (EU) 2017/745		Consult instructions for use
	Keep away from sunlight		Keep in a cool, dry place
	Date of manufacture		Not made with natural rubber latex

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

REF 18/12/I (Gima 27637)

 **Erregi S.a.s.**
Via Adua 40
21040 Vedano Olona (Va)
commerciale@ausilimedici.it
www.ausilimedici.it
MADE IN ITALY



M27637-Rev.1-02.24