

REF

KDB-504E (GIMA 27704)


 Zhongshan Kangdebao(KDB) Rehabilitation Equipment Co.,Ltd.
 Card 1, No.5, 24# Longcheng Road, Dongsheng Town, Zhongshan City, Guangdong Province, China
 Made in China

EC REP


SUNGO Europe B.V.
 Fascinatio Boulevard 522, Unit 1.7, 2909VA
 Capelle aan den IJssel, The Netherlands


Gima S.p.A.
 Via Marconi, 1 - 20060 Gessate (MI) Italy
 gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



NAME
Bedside rail
Model number: KDB-504E

INTENDED USE:
 The device is used as gripping object for preventing the patient, elderly and frail person from slipping over when getting in or out of bed.

Model	Specification	Load-bearing	Material	Picture	Foldable	Height adjustable
KDB-504E	Bedside rail	125KG	Iron		YES	YES

It is warned the users who have allergic reactions to certain latex proteins that feet contain natural rubber

CARE AND MAINTAIN

- To clean, wipe with damp cloth;
- Any additional and replacement parts should only be obtained from manufacturer.

PACKAGE

Seal the Bedside Rail in a plastic bag, then put it in carton.

STORAGE














The packed Bedside Rail should be stored in a ventilated and dry place, no corrosive gas and well ventilated and clean room. To avoid damp, strong light and keep away from heat sources during transportation.

DISPOSAL

Please dispose the product after use to comply with local regulation.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

LABELS, PACKAGING LOGO DESIGN:

Symbol	Meaning	Symbol	Meaning
	Medical device compliant with Regulation (EU) 2017/745		Medical Device
	Manufacturer		Authorized representative in the European Community
	Manufacturing date		Store away from sunlight
	Lot Number		Non-sterile
	Warning: Read and follow the instructions (warnings) for use carefully		Product code
	Read the instructions for use		Store in a cool, dry place
	Imported by		

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.