



# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**SOSTEGNO PER SACCHE URINA  
CON COPERTURA RIUTILIZZABILE**

**HANGER HOLDER WITH REUSABLE COVER FOR URINE  
BAG**

**CROCHET AVEC ENVELOPPE RÉUTILISABLE POUR  
POCHE À URINE**

**GANCHO CON COBERTURA REUTILIZABLE PARA BOLSA  
ORINA**

**SUPORTE DE CABIDE COM CAPA REUTILIZÁVEL PARA  
SACO DE URINA**

**BÜGELHALTER MIT WIEDERVERWENDBARER  
HÜLLE FÜR URINBEUTEL**

**WIESZAK DO WORKÓW NA MOCZ Z POKRYCIEM  
WIELOKROTNEGO UŻYTKU**

**ΚΡΕΜΑΣΤΡΑ ΜΕ ΕΠΑΝΑΧΡΗΣΙΜΟΠΟΙΟΥΜΕΝΟ  
ΚΑΛΥΜΜΑ ΓΙΑ ΟΥΡΟΣΥΛΛΕΚΤΗ**



#### ENGLISH

Hanger holder with reusable cover for urine bag

**Instruction for use:**

- Place the urine bag inside the cover.
- Make sure that the scale printed on the urine bag appears through the window.
- When using urine bag with bottom valve, insert the outlet valve through the bottom hole.
- Place the hanger in the relative holes to hang the urine bag.

**REF 28698**

Importato da / Imported by / Importé par / Importado por / Importado por /  
Eingeführt von / Importowane przez / Εισαγωγή από:

**Gima S.p.A.**

Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com

[www.gimaitaly.com](http://www.gimaitaly.com)

Made in Israel

Fabbricante / Manufacturer / Fabricant n/ Fabricante / Hersteller/ Fabricante /  
Προδουcent/ Κατασκευαστής:

**Mizra Medical Ltd.**

21 Hasadna Street, Industrial Park, Holon, 5881411, Israel  
dror@mizramedical.com