



---

**ENGLISH**

## Instruction Manual

Thermo Feel is a cylindrical device, one end of which is made of steel and the other of plastic. Thanks to the different thermal properties of the two materials, which are placed alternately in contact with the skin (see figure), the patient can be asked whether or not they perceive the temperature variation. Thermo Feel is not a Medical Device, but it is only an indicative support tool for an initial screening of diabetic patients potentially at risk of neuropathy. The results provided are not useful for diagnostic purposes.

### **GIMA WARRANTY TERMS**

The Gima 12-month standard B2B warranty applies

**REF 31278**



Distribuito da / Distributed by / Distribué par / Distribuido por  
/ Distribuído por / Vertrieben von / Διανέμεται από / Distribuit  
de / Dystrybuowane przez / Distributer / Forgalmazza /  
Išplatino / Разпространяван от / توزيع

**Gima S.p.A.**

Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com

[www.gimaitaly.com](http://www.gimaitaly.com)

Made in China