



---

**ENGLISH**

## Instruction Manual

Thermo feel is a cylindrical device, one end of which is made of steel and the other of plastic. The different thermal properties of the two materials, placed alternately in contact with the skin (see figures) allows the doctor to ask the patient if he/she feels the change in temperature.

### **GIMA WARRANTY TERMS**

The Gima 12-month standard B2B warranty applies



Fabbricante/Manufacturer/Fabricant/Fabricante/Fabricante  
Hersteller/Producent/Κατασκευαστής/Gyártó/Gamintojas  
Producător/Proizvajalec/производитель/الشركة المصنعة

**SAMSON INNOVATION CORPORATION LTD**

ROOM F, 8TH FLOOR, HOP HING INDUSTRIAL BUILDING, NO.704  
CASTLE PEAK ROAD, LAI CHI KOK, KOWLOON, HONG KONG  
manager@hksic.com

**REF TF001D**

Importato da/Imported from/Importé de/Importado de/Importado  
de Importiert von/Importowane z/Εισαγωγή από/Importált innen  
Importuota iš/Importat din/Uvoženo iz/Вносот/مستورد من

**Solution Group S.r.l.**

Via Carlo Sormani, 1 - 20095 Cusano Milanino (Mi) - Italy  
solutiongroup@legalmail.it - www.solutiongroup.it  
Made in China

Distribuito da / Distributed by / Distribué par / Distribuido por  
/ Distribuido por / Vertrieben von / Διανέμεται από / Distribuit de  
/ Dystrybuowane przez / Distributer / Forgalmazza / Išplatino /  
Разпространяван от / توزيع

**Gima S.p.A.**

Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

**Gima 31278**