



GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. – Via Marconi, 1 – 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 – fax. 199 400 403
Export: tel. +39 02 953854209/221/225 – fax. +39 02 95380056
gima@gimaitaly.com – export@gimaitaly.com
www.gimaitaly.com

ORL CURETTES FOR WAX REMOVAL



User manual

ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

REF 31470



GIMA S.p.A.
Via Marconi, 1
20060 Gessate (MI)
Italy





Double-Ended Disposable Ear Curettes

Double-Ended Disposable Ear Curettes feature the functionality of metal curettes with the comfort, safety and convenience of disposables. Double-ended disposable ear curettes feature the functionality of metal curettes

- Made of polypropylene, these ear curettes are as effective as metal curettes and virtually unbreakable.
- Rounded tips reduce trauma and lessen the risk of lacerations.
- Eliminate expensive and time-consuming sterilization and the risk of cross contamination with these disposable ear curettes. ·

Double endend curette;designed for normal curettage in a variety of patients / Infant spoon designed for smaller ears and may be used through an otoscope. /· Loop smooth flexible tip that gives the user exceptional control while offering the patient more comfort.