

IMMOBILISATION BELTS

ENGLISH

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

Size: 5x213 cm
 This product is used in stretchers during patient transfer.
 Loading Capacity: More than 200 Kgs.











OPERATION

Unfasten the safety strap, and fix the patient on a stretcher, then fasten the buckle, adjust the length of strap for securing the patient properly.

PRODUCT STORAGE AND MAINTENANCE

After each use, please use clear water or neutral detergent to clean the dirt on the surface. Store it for subsequence use when it is dry. If there are bloodstains or body fluids left, please disinfect the product with disinfectant, and clean up.

SYMBOLS

	Product code
	Caution: read instructions (warnings) carefully
	Medical Device
	Consult instructions for use
	Manufacturer
	Date of manufacture
	Medical Device compliant with Regulation (EU) 2017/745
	Keep in a cool, dry place
	GB Keep away from sunlight
	Lot number

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

REF 34012 - 34013 - 34015

 **Gima S.p.A.**
 Via Marconi, 1 - 20060 Gessate (MI) Italy
 gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
 Made in China

