



CE
EN1865

TAVOLA SPINALE SPINE BOARD
TABLA ESPINAL SPINE BOARDS
PLAN DUR

PVS SpA
Cassina de Pecchi – I

ENGLISH

Specification : Spinal board is made in one piece, high density polyethylene construction is impervious to all body fluids. Complete x-ray translucent. Lightweight and buoyant. Is thin enough to accommodate in the ambulance. Compatible with all head immobilization devices

Spinal boards are to be used by skilled staff. Avoid any contact between the surface and the open wound.

Cleaning : wash with water and soap or bactericidal liquid

DATI DIMENSIONALI - SIZES - DIMENSIONS ABMESSUNGEN - DIMENSIONES	
PORTATA - LOADING CAPACITY - CAPACITE DE CHARGE - BELASTBARKEIT - ALCANCE	150 KG.
LUNGHEZZA - LENGHT - LONGUEUR - ERSTRECKUNGSLAENGE - LONGITUD	184 cm.
LARGHEZZA - WIDTH - LARGEUR - BREITE - ANCHURA	41 cm
SPESSORE - TICKNESS - EPAISSEUR - DICKE - ESPESOR	5 cm.
PESO - WEIGHT - POIDS - GEWICHT - PESO	7,20 KG.