

**G** The **BOSTON** splint is resistant even though very light weighted . They are easily adaptable to an injured arm or leg. The drawings printed on the splint shows how to use them. They are X ray translucent, washable and can be disinfected, and coated with a non-toxic material. The **BOSTON** splint can be re-used, but has to be cleaned and kept as instructed. For more security it is indicated that the **BOSTON** splint should not be re used if in contact with body fluid or blood. Store in a clean and dry place, far from heat and sunrays. Has to be protected in a bag or case if to be re-used. Wash with water and soap or bacterial cleaning solution. Dry they well. Made in PRC.