

# ***INSTRUCTIONS FOR USE***

## ***NON-STERILE ADHESIVE PLASTERS***

**PRODUCT TRADE NAME:** PAP02NBSC

**MANUFACTURER:** RAYS S.p.A. – Via Francesco Crispi 26 – 60027 – Osimo (AN)

**STORAGE:**

Store in a cool, dry place away from heat sources.

**SINGLE USE INDICATION:** Single use medical device

**INTENDED USE:**

Medical device for dressings.

**INSTRUCTIONS FOR USE:**

Wash your hands thoroughly

Open the package.

Adapt the adhesive part to clean, dry skin.

Renew the dressing if necessary.

Do not apply directly to damaged skin.

In case of reuse, infections may develop: consult a doctor.

**LATEX FREE:**

The device is latex-free to avoid the risk of allergic reactions.

**DISPOSAL:**

Disposal must be carried out in compliance with current national regulations.

Incineration or landfill disposal must be carried out under controlled conditions. The

product that has not meet infectious agents or other risky materials cannot be

considered hazardous waste.

*These Instructions for Use note must be preserved by the user for the duration of the DM.*



**Health&Safety**

**RAYS SPA**

via Francesco Crispi - 60027 Osimo (AN)

Tel. 071 2868935 - Fax 071 2868911 - [www.rays.it](http://www.rays.it) - [info@rays.it](mailto:info@rays.it)