

INSTRUCTIONS FOR USE

NON-STERILE ADHESIVE PLASTERS

PRODUCT TRADE NAME: PAP03NBSC

MANUFACTURER: RAYS S.p.A. – Via Francesco Crispi 26 – 60027 – Osimo (AN)

STORAGE:

Store in a cool, dry place away from heat sources.

SINGLE USE INDICATION: Single use medical device

INTENDED USE:

Medical device for dressings.

INSTRUCTIONS FOR USE:

Wash your hands thoroughly

Open the package.

Adapt the adhesive part to clean, dry skin.

Renew the dressing if necessary.

Do not apply directly to damaged skin.

In case of reuse, infections may develop: consult a doctor.

LATEX FREE:

The device is latex-free to avoid the risk of allergic reactions.

DISPOSAL:

Disposal must be carried out in compliance with current national regulations.

Incineration or landfill disposal must be carried out under controlled conditions. The product that has not meet infectious agents or other risky materials cannot be considered hazardous waste.

These Instructions for Use note must be preserved by the user for the duration of the DM.



Health&Safety

RAYS SPA

via Francesco Crispi - 60027 Osimo (AN)

Tel. 071 2868935 - Fax 071 2868911 - www.rays.it - info@rays.it