

**USE WARNINGS**  
**CE IN-VITRO DIAGNOSTIC DEVICE - For professional healthcare use**

- In-vitro diagnostic devices: test tubes, containers for biological liquids collection, Petri dishes, microbiology spreaders, sample cups for autoanalyzer. Only for collection of specimens, other than blood, from humans. Do not use the product for any other use which would be considered improper. Only use the device with suitable accessories.
- Use protective gloves.
- Do not place the device near flame.
- Do not use the product if it is out of date, damaged or if the package is open.
- Keep out of the reach of children.
- Store in a dry place, away from sources of heat and direct sunlight.
- Dispose of the product in accordance with current legislation in force, in a biohazardous waste container.
- Before using the product with particular substances, check the table of material resistance/compatibility.
- **Sterile product: sterility is guaranteed until the packaging is opened. Sterilization is performed by Beta rays.**
- The device is disposable. The reuse of the device will entail the loss of all performances and features.
- Any serious incident occurring in relation to the device must be reported to the manufacturer and to the Competent Authority of the State where the patient resides.



F.L. Medical s.r.l. Unipersonale – Via Enrico Mattei, 20- 35038 Torreglia (PD) – Italy  
Phone +390495211475 – Fax +390495212566 – mail info@flmedical.com

AL536\_Use Warnings IVD\_FI Medical



Rev. 2 del 12.04.2022