



AZIENDA CON SISTEMA DI GESTIONE
PER LA QUALITÀ CERTIFICATO DA DNV
= UNI EN ISO 9001:2008 =

VAIMAS SRL

ART. 90120150 NORMAL CHEST BELTVCA 150
ART. 90120155 CHEST BELT WRAP PADDED VCA 155
ART. 90120160 CHEST BELT WITH SUSPENDERS PADDED VCA 160



MANUAL

PRODUCT COMPLIANT **CE**

DIRETTIVA 93/42 CEE
AUSILI MEDICI CLASSE I

DESCRIPTION

Class 1 medical devices manufactured by VAIMAS SRL Marano Vicentino in compliance with Directive 93/42 EEC.

INDICATION

To prevent the patient sitting on an wheelchair from sliding forward, with sturdy straps and adjustable shoulder straps.
Easy to clean.

MAINTENANCE

Do not machine wash or high temperatures (over 30 ° C) for not compromising materials, to dry away from heat sources. Store and transport in a dry place at room temperature.

WARNINGS

The device meditate should be prescribed by a doctor and applied by an orthopedic technician or by a qualified professional. The orthotist is the professional responsible for the application and reference information for its safe and proper use in accordance with individual needs.

Any further changes in the structure or regulation and use of the device must be prescribed by a physician and performed by an orthopedic technician or qualified personnel. In case of doubt on the application of the device, consult an orthopedic technician or qualified personnel.

The device is designed to be well tolerated by the patient.

SAFETY



- Use only for the purposes intended by the manufacturer.
- The manufacturer declines any liability for damage caused by incorrect or not provided in this manual.
- Any use other than that intended by the manufacturer may cause damage.
- Any modification or refurbishing must be approved by the manufacturer.

WARNING: The information contained in this document is subject to change without notice.