



TITANOX s.r.l. - via Canove 2/A - Loc. Canove dé  
Biazzì - 26038 Torre dé Pícenardi (CR) - Italy -  
www.titanox.it [info@titanox.it](mailto:info@titanox.it) MADE IN ITALY

## Instruction manual “Mayo table”

Stainless steel base, 5 spokes with castors diam. 50mm. Frame and removable tray made of stainless steel. Adjustable height from 90 to 115 cm.

Tray dimensions: 70x45x1,7h cm.

All models are produced eliminating sources of potential risk to the user, should unintended collisions occur, such as cutting edges or sharp corners.

These instruments are intended to be used by medical practitioners who are specially trained on how to use and care of them. Use for purposes other than what it is intended for is not allowed. The incorrect use, poor or inappropriate maintenance can rapidly lead to deterioration of the instruments.

- Check the right tightening torque of all mobile parts before and during use.
- For cleaning, do not use abrasive substances, acids, alcohol, chlorine-based cleaners, disinfectants, and acetone: the manufacturer company cannot be held responsible for the damage that is caused by using materials that could damage the product surface or corrosive chemicals during cleaning.



TITANOX s.r.l. - via Canove 2/A - Loc. Canove dé  
Biazi - 26038 Torre dé Picenardi (CR) - Italy -  
[www.titanox.it](http://www.titanox.it) [info@titanox.it](mailto:info@titanox.it) MADE IN ITALY










## TITANOX WARRANTY TERMS

Titanox 12-month warranty applies from final invoice date.

**All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.**


**REF M600480**



	<b>IT</b> Fabbricante <b>ENG</b> Manufacturer <b>FR</b> Fabricant
	<b>IT</b> Attenzione: Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso <b>ENG</b> Caution: read instructions (warnings) carefully <b>FR</b> Attention: lisez attentivement les instructions (avertissements)
	<b>IT</b> Leggere le istruzioni per l'uso <b>ENG</b> Consult instructions for use <b>FR</b> Consulter les instructions d'utilisation
	<b>IT</b> Conservare al riparo dalla luce solare <b>ENG</b> Keep away from sunlight <b>FR</b> Á conserver à l'abri de la lumière du soleil
	<b>IT</b> Conservare in luogo fresco ed asciutto <b>ENG</b> Keep in a cool, dry place <b>FR</b> Á conserver dans un endroit frais et sec
	<b>IT</b> Codice prodotto <b>ENG</b> Product code <b>FR</b> Code produit
	<b>IT</b> Numero di lotto <b>ENG</b> Lot number <b>FR</b> Numéro de lot
	<b>IT</b> Dispositivo medico <b>ENG</b> Medical Device <b>FR</b> Dispositif médical
	<b>IT</b> Il dispositivo è conforme ai requisiti del regolamento sui dispositivi medici (UE) 2017/745 ed è classificato come dispositivo medico di classe I <b>ENG</b> The product complies with the requirements of the regulation on medical devices (EU) 2017/745 and is classified as a Class I medical device <b>FR</b> Le produit est conforme aux exigences du règlement sur les dispositifs médicaux (UE) 2017/745 et est classé comme dispositif médical de classe I



Manufacturer: TITANOX s.r.l. - via Canove 2/A -  
Loc. Canove de' Biazzi - 26038 Torre de' Picenardi  
(CR) - Italy - [www.titanox.it](http://www.titanox.it) MADE IN ITALY



MANTENERE IN ASSE CON TRAVERSO SUPERIORE  
KEEP IN AXIS WITH THE UPPER CROSSPIECE  
MANTENIR DANS L'AXE DE LA TRAVERSE SUPERIEURE

M600480