



FIRST AID FOR BURNS


91cm x 76cm (36" x 30")
FIRE/TRAUMA BLANKET FOR FIRST AID TREATMENT OF BURN INJURIES

FIRST AID FOR BURNS

- STOPS THE BURNING PROCESS
• REMOVES HEAT FROM BURN SITE
• EXTINGUISHES SMALL FIRES

CONTENTS: One gel soaked Fire/Trauma Blanket in pouch.

INSTRUCTIONS: Remove blanket from inner bag. Cover and conform to burned area. Obtain medical attention as soon as possible.



WARNING: Electricity must be shut off before applying blanket. Chemical burns must be thoroughly flushed with water before using blanket. This product is designed for temporary use only. If packaging has been damaged, DO NOT USE.

Water-Jel Technologies, LLC
50 Broad Street, Carlstadt, NJ 07072 USA
Safeguard Technologies Ltd
Willow Grove, Delgany, Co Wicklow A63 XY89, Ireland
PART OF SAFEGUARD MEDICAL
www.burnfree.com MEDICAL DEVICE MADE IN USA

Storage limit temperature icon (15°C to 25°C), Latex Free icon, CE mark, Lot icon, Barcode, REF #SB-9176 04/2021

Updated Company Info

Updated temp icon and rev date. CE slightly larger. Added Latex Free & Do not Reuse

ARTWORK SPECIFICATIONS
CLIENT INFO
CLIENT NAME: WATER-JEL
PRODUCT: BURN FREE 3630 DRESSING
ITEM NUMBER: BFT 7/6.00.018
FILE INFO
FILE NAME: WJ BF 9176.00.018 3630
PROGRAM: ILLUSTRATOR CC
MATCH COLORS IN PROCESS:
PANTONE REFLEX/BLUE C 100%
PANTONE YELLOW 485 C
PANTONE BLACK C 2835 C
PANTONE 485 C
PANTONE BLACK C
CYAN DIELINE DOES NOT PRINT

UNWIND CHART/COPY POSITION
DESCRIPTION: 1183 Edgewater Ave • Ridgefield, NJ 07657
OPERATOR:
# OF COLORS:
DATE: 8/25/16
ITEM #:
TPO#:
REFERENCE#:
LABEL SIZE: 12.5 X 8.25
DIE SPEC:
RECTANGLE
DATE: 8/25/16
WARNING: These colors represented on this proof are not necessarily actual colors, variations in tone and intensity will occur depending on actual material used.
DO NOT PRINT
DIE VINYL
IMPORTANT: This art has been checked and proofed for accuracy. It is the responsibility of the client to make all final approvals before the release of this art. PLEASE DOUBLE CHECK FOR ACCURACY.
X
Date:
PLEASE SIGN AND RETURN FOR COPY CORRECTNESS, COPY POSITION AND COLOR BREAKS
PLEASE NOTE: PMS COLORS ON PROOF ARE NOT EXACT, PLEASE CONSULT A PMS BOOK FOR COLOR ACCURACY