



FIRST AID FOR BURNS


244cm x 152cm (96" x 60")
FIRE/TRAUMA BLANKET FOR FIRST AID TREATMENT OF BURN INJURIES

FIRST AID FOR BURNS

- STOPS THE BURNING PROCESS
- REMOVES HEAT FROM BURN SITE
- EXTINGUISHES SMALL FIRES

CONTENTS: One gel soaked Fire/Trauma Blanket in pouch.

INSTRUCTIONS: Remove blanket from inner bag. Cover and conform to burned area. Obtain medical attention as soon as possible.



WARNING: Electricity must be shut off before applying blanket. Chemical burns must be thoroughly flushed with water before using blanket. This product is designed for temporary use only. If packaging has been damaged, DO NOT USE.

Water-Jel Technologies, LLC
 50 Broad Street, Carlstadt, NJ 07072 USA

EC REP Safeguard Technologies Ltd
 Willow Grove, Delgany, Co Wicklow A63 XY89, Ireland

PART OF **SAFEGUARD MEDICAL**

www.burnfree.com MEDICAL DEVICE **MADE IN USA**

REF #SB-244152 04/2021

15°C 25°C

STORAGE LIMIT TEMPERATURE

5 025646 090746

Updated Company Info

Updated temp icon and rev date.
 CE slightly larger. Added Latex Free & Do not Reuse

UNWIND CHART / COPY POSITION

1183 Edgewater Ave • Ridgefield, NJ 07657
 Tel: (201) 941-6611 • Fax: (201) 941-6307
 e-mail: trend@trendlabel.org

DESCRIPTION: 1183
 OPERATOR: _____
 # OF COLORS: _____

DATE: 8/25/16
 ITEM #: _____
 TPO#: _____
 REFERENCE#: _____

WARNING: These colors represented on this proof are not necessarily actual colors. Variations in tone and intensity will occur depending on actual material used.

DO NOT PRINT
 DIE VINYL
 DO NOT PRINT

IMPORTANT: This art has been checked and proofed for accuracy. It is the responsibility of the client to make all final approvals before the release of this art. PLEASE DOUBLE CHECK FOR ACCURACY.

X (Please sign and return for copy correctness, copy position and color breaks)
 PLEASE NOTE: PMS COLORS ON PROOF ARE NOT EXACT, PLEASE CONSULT A PMS BOOK FOR COLOR ACCURACY

ARTWORK SPECIFICATIONS

c-squared design inc. 333 Fairview Avenue, Westwood, NJ 07675
 Phone: 201.664.2212 • Fax: 201.664.5566

CLIENT INFO
 CLIENT NAME: WATER-JEL
 PRODUCT: BURN FREE 9690 DRESSING
 NET WT: n/a
 UPO OR PRODUCT CODE: 5025646090111
 ITEM NUMBER: BF244152.00.016
 FIG NO.: _____

FILE INFO
 FILE NAME: NU BF244152.00.016.9690
 PROGRAM: ILLUSTRATOR CC
 BACK-UP: INCLUDED IN FOLDER
 FONTS: CONVERTED TO OUTLINE

MATCH COLORS IN PROCESS:

	PANTONE REFLEX/BLUE C 100%
	PANTONE YELLOW 485 C
	PANTONE 485 C
	PANTONE BLACK C
	PANTONE 2895 C
	CYAN DIELINE DOES NOT PRINT

This is a digital color mechanical file. No traps have been made. It is the responsibility of the printer to identify to the firm separator/color make the required press specifications.

2: Design has made every effort to ensure the accuracy of the information included in this file. Final responsibility for accuracy of all information, spelling, and legal approvals remains with the client.