







CUSCINO ROTANTE ROTATING SEAT CUSHION

Manuale d'uso - User's manual

Ausilio per agevolare l'entrata o uscita dall'auto, dal letto e sulla sedia.
To assist someone to turn easily in a car, from the bed or on the chair.



ATTENZIONE: usare molta cautela in caso di utilizzo su superfici non piane o scivolose. Non indicato per pazienti affetti da gravi patologie motorie. Si raccomanda di controllare periodicamente la robustezza del prodotto. Per evitare eventuali rischi tenere lontano dalla portata dei bambini e controllare sempre lo stato.

Materiale: parte superiore in velluto, interno in nylon rivestito in silicone, fondo in PVC antiscivolo

Imbottitura: schiuma EPE

ATTENTION: use caution on slippery surfaces. Not indicated for patients with serious motor pathology. Check the status periodically. Keep it away from children to avoid any possible risks. Please keep under control before use.

Material :Upper with velour material ,inner with silicon coated nylon ,bottom with non-slip PVC

Padded: EPE foam

E' necessario segnalare qualsiasi incidente grave verificatosi in relazione al dispositivo medico da noi fornito al fabbricante e all'autorità competente dello Stato membro in cui si ha sede. All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.

ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

REF FC150026 (GIMA 43040)



 **HuiZhouShi YiHuKang Trading Co. Ltd.**
No.105 ,JinFuYuan, NanEr Rd., XinNa , XinWuXiaoZu ,Bogang Village, XiaoJinKou Street,
HuiCheng District ,Huizhou City,Guangdong, China
Made in China



Gima S.p.A.
Via Marconi, 1 20060 Gessate (Mi) - Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



MedPath GmbH
Mies-van-der-Rohe-Strasse 8
80807 Munich, Germany

