

135mm

25mm

13mm

GiSafe Mask - Colours

- Surgical Disposable 3-PLY Face Mask with ear loops - Type IIR
- Mascherina chirurgica monouso a 3 veli con elastici - Tipo IIR
- Masque chirurgical jetable à 3 couches avec élastiques - Type IIR
- Mascarilla quirúrgica desechable a 3 capas con elásticos - Tipo IIR
- Máscara Facial Cirúrgica Descartável 3-camadas com alças para orelha - Tipo IIR
- 3-lagige Einweg-OP-Schutzmaske mit Gummibändern - Typ IIR
- Χειρουργική μάσκα προσώπου μίας χρήσης 3-στρωσεων - με θηλιές αυτιών - Τύπου IIR
- قناع طبي للوجه يستخدم لمرة واحدة 3 طبقات - مع حلقتين للأذن - نوع IIR



10 PIECES PEÇAS
PEZZI STÜCK
UNITÉS TEMAXIA
PIEZAS قطعة

BFE
>98%



EN14683:2019+AC(2019)

www.gimaitaly.com



GiSafe Mask - Colours

- Non-woven material
• Class I medical device
• Comfortable and pliable nose piece
• Fiber glass free
- Materiale non tessuto
• Dispositivo medico di classe I
• Stringinaso modellabile
• Priva di fibre di vetro
- Matériau non tissé
• Dispositif médical de classe I
• Pinze-nez modelable
• Sans fibres de verre
- Material no tejido
• Producto sanitario de clase I
• Faja para nariz modelable
• Sin fibra de vidrio
- Material não tecido
• Dispositivo médico Classe I
• Peça para nariz confortável e flexível
• Sem fibra de vidro
- Vlies
• Medizinprodukt der Klasse I
• Verformbares bequemes Nasenstück
• Glasfiberfrei
- Μη υφασμένου υλικού
• Ιατροτεχνολογικό προϊόν κατηγορίας I
• Άνετη και εύκαμπτη στερέωση μύτης
• Χωρίς ίνες γυαλιού
- مادة خام غير صوفية من
• جهاز طبي فئة I
• قطعة للأنف مريةة وقابلة للطي
• خال من الألياف الزجاجية



生产厂商: 湖北飞腾防护用品有限公司
生产地址: 湖北省仙桃市纺织工业园区管沟社区18号
执行标准: EN14683:2019+AC(2019)
使用年限: 3年
生产日期:



8 023279 206883



000002



REF ELFM3600 (20688)

LOT



Eldo Care Co., Limited
Rm 1902, Easey Comm Bld 253-261
Hennessy Rd. Wanchai, Hong Kong
Made in China

EC REP Lotus NL B.V.
Koningin Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands

Imported by / Importato da:
Gima S.p.A. Via Marconi, 1
20060 Gessate (MI) Italy
gima@gimaitaly.com
export@gimaitaly.com
www.gimaitaly.com



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



S09553 Rev. 1 - 10/21 Size 13,5x25 cm