



URINE SCREEN 10



URINE SCREEN 11

Parametri / Parameter

Descrizione Description	Codice prodotto Product code	Contenuto Content	Glucosio Glucose	Acido ascorbico Ascorbic Acid	Chetoni Ketones	Proteine Protein	Valori pH pH-Value	Sangue Blood	Nitriti Nitrite	Leucociti Leucocytes	Peso specifico Spec. Gravity	Bilirubina Bilirubin	Urobilinogeno Urobilinogen
URINE SCREEN 10	24073 / 93120	100	■		■	■	■	■	■	■	■	■	■
URINE SCREEN 11	24074 / 93100	100	■	■	■	■	■	■	■	■	■	■	■



References / Referenzen / Referenser / Referencer / Viitteet / Referanser / Referanslar / Referencias / Références / Riferimenti / Referências / Referentיות / Referençe / Reference / βιβλιογραφικές αναφορές / Справочный
 Referenzbereiche für Kinder und Erwachsene von Heil/Ehrhardt (Roche) [pH Referenz daraus entnommen]; oder alternativ aus „Textbook of Urinalysis and Body Fluids“ von Landy J. McBride: Kaplan L.A., Pesce A.J. Clinical chemistry. 3rd ed. St. Louis: The CV Mosby Company, 1996.

