

GYNAECOLOGICAL AND DELIVERY BEDS



Standard leg holders



Quality leg holders




27510

- 27500 DELUXE GYN BED - quality leg holders
 - 27501 GYNAECOLOGY BED - standard leg holders
- Strong structure in chrome-plated. Rectangular steel tube. Adjustable head and leg sections. Simulated leather, black colour, washable upholstery, self extinguishing. Load: 100 kg. Complete with s/s bowl Ø 31 cm and leg holders. Made in Italy.
- 26586 S/S BOWL Ø 31 cm - spare

- 27510 DELUXE DELIVERY BED - quality leg holders
- Same as 27500, but with three sections with independent controls and Trendelenburg position. Load: 120 kg.
- 27503 STANDARD LEG HOLDERS - couple - spare
 - 27504 QUALITY LEG HOLDERS - couple - spare
- Kit leg holders include 2 leg holders and arms.
- 27418 COUCH ROLL HOLDER

3- SECTION MULTIFUNCTIONAL AND GYNAECOLOGICAL BEDS



- 27502 3-SECTION MULTIFUNCTIONAL BED - blue
 - 27505 3-SECTION GYNAECOLOGICAL BED - blue
- Strong and stylish structure in oval tube 50x25 mm, painted with epossidic powder, metallic gray colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material Class 1 IM. Back-rest and leg-rest section are fully adjustable by self mechanical sytem. Gynaecological bed is supplied with 2 foam leg holders and stainless steel bowl Ø 32 cm. Maximum load: 130 kg - Weight: 36 kg. Made in Italy.



- 27508 STANDARD LEG HOLDER - couple - spare
 - 27509 QUALITY LEG HOLDER - couple
- Kit leg holder, include 2 leg holder and arms.
- 27511 COUCH ROLL HOLDER for 27502-5
 - 27919 DESK with 2 drawers - 130x71x80h cm
- Strong structure in oval tube 40x20 painted with epossidic powder, metal gray colour. White work top with rubber edge and round corner.
- 27979 FOOT STOOL - 2 steps



HEIGHT ADJUSTABLE GYNAECOLOGICAL BED

- 27506 HEIGHT ADJUSTABLE GYN BED - green
 - 27507 HEIGHT ADJUSTABLE GYN BED - blue
- Strong and solid structure in oval tube 50x25x2 mm painted with epossidic powder, white colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material class 1 IM.
- Full adjustable**
- Height, up to 95 cm, is adjustable by electric engine (6000 N power) with foot control, head section is adjusted by a gas cylinder (400 N power), central section by a gas cylinder (600 N power) up to 65°. This gyn bed is supplied with 2 quality foam leg holders, stainless steel tray and couch roll holder.
- Maximum load: 130 kg.
Weight: 70 kg Made in Italy.
- 27509 QUALITY LEG HOLDER - couple - spare
- Kit leg holder include 2 leg holder and arms. Suitable for 27505-27506-27507.

