


MICROALBUMIN URINE SCREENING TEST - KIDNEY DISEASE

• 24539 MICROALBUMIN URINE TEST - box of 1 test

An easy-to-use visual test to detect elevated levels of albumin in urine, an early warning sign of kidney disease. Early detection is key: two main causes of kidney disease are **diabetes & high blood pressure**.

- excellent correlation to lab test
- results in 10 minutes
- simple, easy to use urine test
- sensitive; detects 18 mg/L of albumin in urine (same as professional tests)
- contains everything needed for consumer to perform test
- CE approved as self-testing device
- American Diabetes Association and American Kidney Foundation recommend that those with diabetes test annually.



SELF TEST


COLON DISEASE TEST, FOR SELF TESTING AT HOME

SELF TEST

• 24531 COLON DISEASE SELF TEST

Simple home test for detecting the early warning signs of colorectal disease, even before symptoms appear. Contains 5 test pads, 1 positive control package, results card and 1 instruction manual.

- easy throw in the bowl test for hidden blood in the stool
- no laboratory processing required, results in 2 minutes
- the American Cancer Society recommends annual testing after age of 50
- with early detection and treatment colorectal cancer is curable



Step1: after a bowel movement drop a test pad into the toilet bowl.

Step2: observe the test pad. If blood is present a blue color will appear within two minutes

HELICOBACTER PYLORI TEST FOR ULCER SCREENING

PROFESSIONAL


• 24528 HELICOBACTER PYLORI TEST - box of 40 tests

Helicobacter Pylori is the leading cause of gastric ulcers, this easy to use test allows a quick diagnosis in about 10 minutes. Easy to perform and interpret, no additional equipment required. For professional use.


Diagnosis of H. pylori infection in adults over 18 years of age.

• 24529 HELICOBACTER PYLORI SELF TEST - box of 1 test

Same as 24528 but packed in single box including one automatic lancet, single buffer and user manual. This product is provided with self testing certification to be sold to end users in OTC market.



24528 CE



SELF TEST

PROSTATE SCREENING - PSA ONE-STEP TEST

• 24534 PSA TEST - box of 25 tests - for doctor/hospital use

Provides accurate results in 5-8 minutes.

- over 98% correlation to lab tests
- simple, easy to use serum or plasma test
- detects PSA as low as 4 ng/ml
- contains everything needed to perform test (except lancets see page 5)
- for professional use



CE 0197

24532

CHOLESTEROL TEST

SELF TEST

• 24521 CHOLESTEROL STRIPS - box of 2 pcs.

Inexpensive, CE marked strips for visual reading of cholesterol in whole blood.

No meter is needed to measure the results.

The special color chart allows an easier quantification of the result. Cholesterol is offered in a 2 tests tube package inclusive of fingertips, color chart and instructions for use. Box and instructions in English, German, Italian.



24521

ALLERGY TEST

• 24580 CAT ALLERGY TEST - box of 1 test - CE 0197

• 24581 DOG ALLERGY TEST - box of 1 test - CE 0459

Simple, reliable and clinically proven screening test for allergies:

- result in 15 minutes
- excellent correlation to lab test
- easy to use a fingerstick whole blood
- for self testing
- each test includes: test device, vial of test solution, fingerstick device, user manual

SELF TEST



24580