

DOCTOR CABINET



Metal tubular frame with fixed height 100 cm and flexible arm 38 cm.

• 27990 STUDIO LIGHT
Metal tubular frame with plastic base Ø 60 cm with 5 plastic wheels.


TECHNICAL SPECIFICATIONS
100 Watt lamp (not included)
Power supply: 230 V

GIMA code	Product	Description	Size	Weight
27920	DESK - with 2 drawers	Desk with 2 side drawers with locking top drawer and enamelled tubular steel frame	120x70x78	31
27403	EXAMINATION COUCH	Single-hinged examination table with chromed steel tubular frame. Comes complete with simulated leather mattress (grey), backrest adjustable through different positions. Delivered in kit form.	176x60x78	22
27984	WASTE BIN - enamelled steel	Made of enamelled steel with plastic inner bucket - capacity 7 l	Ø25xh 35 cm	2
27985	WASTE BIN - stainless steel	Made of AISI 304 stainless steel with plastic inner bucket - capacity 7 l	Ø25xh 35 cm	2
27989	COAT HANGER	Free-standing coat hanger with chromed steel frame with 8 knobs and umbrella holder		5
27990	LIGHT	Not medical. 100 W. Aluminium tubular frame with 5 spoke base on 50 mm castors. Adjustable height thanks to flexible arm. Lamp not included.		3

CABINET WITH 1, 2, 4 DOORS AND WALL CABINETS



Hanging cabinet with 2 lockable doors. Inside division in 2 compartments left side with 2 removable shelves, right side with 1 removable shelf.



Hanging cabinet to keep medicine with 2 glass tempered doors and 2 tempered glass shelves adjustable on master rack to be hung by two holes on the back

GIMA CABINETS Choice of cabinets with epoxy coated steel structure and tempered glass. Rack-mounted adjustable glass shelves. Rubber feet.

TECHNICAL SPECIFICATIONS	27901	27905	27903	27911	27912
Door:	1	2	4	2	2
Adjustable shelves:	3	4	3	2	2
Weight kg:	30	47	45	16	18
Size cm:	53x36x144 h	75x38x169 h	75x38x169 h	60x30x60 h	80x33x45 h

- 27901 CABINET - 1 door
- 27905 CABINET - 2 doors
- 27903 CABINET - 4 doors
- 27911 WALL CABINET - glass sliding door
- 27912 WALL CABINET - painted door