


DOCTOR LIFE AIR COMPRESSION THERAPY SOLUTIONS

All human body is supplied by peripheral circulation of blood and lymph

Doctor Life system helps patients to recover from their disease by improving blood and lymph circulation and metabolism

Doctor Life system consists of air chambers and an air pump offering air compression massage to legs, arms or waist.




- 28430 DL 2002 D BASIC with 2 legs (98x38 cm) - 4 chambers 1 single and 1 double hose and cable
- This system can operate with 3 compression modes.

TECHNICAL SPECIFICATIONS

- pressure sensor
- digital timer (10, 20, 30, 40 min)
- digital pressure gauge (10-180 mmHg)
- one touch connection hose plug
- body size: 190x250x210 mm
- body weight: 3 kg

- 28435 LEG 93x32 cm
- 28434 LEG EXTENSION ZIPPER - 10 cm extra extension
- 28440 DOUBLE HOSE - to use two leg cuffs at the same time - spare
- 28436 ARM CUFF unique size - 86x29 cm
- 28437 WAIST CUFF 130x38 cm



3 compression modes

Medical indications

- DVT (Deep Vein Thrombosis)
- Lymph edema
- Post-paralytic
- Limb convulsion
- Venosclerosis
- Chronic venous insufficiency (CVI)
- Post-mastectomy
- Post-trauma
- Spinal cord injury
- Varicose veins
- Sciatica
- Sprains
- Diabetes
- Cellulite

Home care indications

- Lymph edema
- Diseases caused by interrupted blood circulation
- Fat dissolution
- Pregnant woman's limb treatment
- Benumbed feet or hands
- Cold feet or hands
- Settlement of stress
- Exercise effects
- Insomnia
- Healthier body and elastic skin

Drainage, Circulation & Massage

SEMI-PROFESSIONAL MODEL (Home-Hospital)

Universal cable included is suitable for leg, arm and waist cuff



Waist cuff 130x38 cm (28437)




Possibility to use 2 leg cuffs in the same time with double hose (code 28440)

HIGH CLASS PROFESSIONAL MODEL


- 28450 DL850 PROFESSIONAL SYSTEM with 2 legs (102x35 cm), 1 arm - 6 chambers (87x25 cm), double hose and cable
- This system can operate with 8 compression modes. Every chamber can have a different pressure or can work only with 1 selected chamber. Graphic LCD display. Inflation and interval time control by two separate keys.

TECHNICAL SPECIFICATIONS


- Pressure sensor
- Use time: 0-90 minutes
- Hold time: 0-6 seconds
- Interval time: 0-19 seconds
- Body size: 370x390x190 mm
- Body weight: 13 kg



Standard arm (28454)



Standard leg (28453)



- 28453 LEG CUFF 92x31 cm
- 28458 DOUBLE HOSE - to use two leg cuffs at the same time - spare
- 28459 LEG EXTENSION ZIPPER - 10 cm extra extension
- 28454 ARM CUFF - 6 CHAMBER 87x25 cm - spare
- 28455 WAIST CUFF 128x51 cm