

## NEBULIZERS FOR PROFESSIONAL AND HOSPITAL USE



28089

Nasaljet is a quick effective nasal irrigation device for rinsing the nasal and rhino-pharyngeal cavities



28090

CONTINUOUS USE COMPRESSOR

## PROFESSIONAL NEBULIZERS

Piston compressor nebulizers for professional and home care use. They embody all the dictates of modern aerosol therapy techniques. Reliable and simple to use, also by inexperienced users. Manual: GB, FR, IT, DE, ES. Made in Italy.

GIMA code	PROFESSIONAL NEBULIZER	Voltage
28134	HOSPINEB PROFESSIONAL	230 V - 50 Hz
28136	HOSPINEB PROFESSIONAL	230 V - 60 Hz
28137	HOSPINEB PROFESSIONAL	110 V - 60 Hz
28090	COMP-A NEB PROFESSIONAL	230 V - 50 Hz
28089	NASALJET	-
28145	FILTER - spare	-

• 28105 CORSIA PROFESSIONAL NEBULIZER - piston

• 28092 CORSIA PROFESSIONAL NEBULIZER - piston - without manometer

Piston-type compressor nebulizers designed for continuous use. High performance with any type of drug, ideal for intensive hospital and clinic use. Recommended for atomizing antibiotics and bronchodilator drugs. Manufactured with high thermal and electric

CE 0476



28105

CONTINUOUS USE COMPRESSOR



28134

CE 0434


Metal box

## TECHNICAL SPECIFICATIONS

	HOSPINEB	COMP-A NEB
Operating voltage:	110 V - 60 Hz 230 V - 60 Hz	230 V - 50 Hz 110 V - 60 Hz on request 230 V - 60 Hz on request
Size:	37x10xh 24 cm	19.6x35.7xh 18.5 cm
Max pressure:	350 kPa (3.5 bar)	280 kPa
Noise level:	60 dBA	55 dBA
Operating pressure:	0-130 kPa (0-1.3 bar)	0-100 kPa
Weight:	2.6 kg	2.5 kg
Air Flow:	15 l/min	15 l/min
MDD class:	IIa	IIa

## STANDARD ACCESSORIES

Fasterjet nebulizer with valve system  
Nosepiece irrigation (polypropylene)  
Mouthpiece (polypropylene) with valve (silicone)  
Adult and child mask (PVC DEHP free)  
Air tube (PVC DEHP free) and spare filter.  
Nasaljet (only with 28090).



insulation plastic chassis in accordance with the latest European safety regulations. The oil-free piston compressor has long durability and is equipped with the highly efficient HI-FLO jet nebulizer to guarantee quick and accurate drug delivery. Easy to handle and transport.

• 28104 FILTER - spare - box of 2

28092



## TECHNICAL SPECIFICATIONS

Operating Voltage:	220-230 V - 50/60 Hz
Power consumption:	184 VA
Max pressure:	250 kPa (2.5 bar)
Air flow:	16 l/min
Operating pressure:	130 kPa (1.30 bar)
Operating air flow:	5.2 l/min at 130 kPa
Nebulization rate:	0.50 ml/min with 4 ml solution
MMAD:	3.25 µm
GSD:	3.45
Size:	255x190xh 105 mm
Weight:	2.2 kg
Noise level:	approx 57 db
Made in Italy	

## STANDARD ACCESSORIES

HI-FLO Jet bulb  
Mouthpiece  
Adult and child mask  
Air tube and air filter

## ATOMIZER 2 SPEEDS, PISTON NEBULIZER FOR FAMILY USE

• 28094 ATOMIZER NEBULIZER - 2 speeds

Piston-type compressor with 2 speeds for family use. With accessory compartment separated from cord compartment for improved hygiene. Made in Italy.

• 28093 FILTER - spare - box of 3

## TECHNICAL SPECIFICATIONS

Operating voltage:	230 V - 50 Hz
Operating cycle:	continuous
Operating pressure:	speed 1: 45 kPa - speed 2: 60 kPa
Air Flow:	speed 1: 8 l/min - speed 2: 11 l/min
Nebulization rate:	speed 1: 0.2 ml/min - speed 2: 0.35 ml/min
Particle size:	80% < 5 µm
Max pressure:	1.85 bar / 185 kPa
Weight:	1.4 kg
Noise level:	56 db
Class of risk:	IIa



Standard accessories

CE 0434

28094

## STANDARD ACCESSORIES

Nebjet bulb with connector  
Adult and child mask  
Mouth piece and nose piece  
Air tube and 3 air filter



Accessories compartment