



## "KARMAN" AND "BERKELEY" INTRAUTERINE SUCTION CURETTES




Length 24 cm with connector, 20,5 cm without

GIMA code	"KARMAN" STYLE CURETTES FLEXIBLE - SINGLE USE - STERILE	Minimum order
29860	Karman Curette N° 6 - 6 mm x 24 cm	box of 20
29861	Karman Curette N° 7 - 7 mm x 24 cm	box of 20
29862	Karman Curette N° 8 - 8 mm x 24 cm	box of 20
29863	Karman Curette N° 9 - 9 mm x 24 cm	box of 20
29864	Karman Curette N° 10 - 10 mm x 24 cm	box of 20
29865	Karman Curette N° 11 - 11 mm x 24 cm	box of 20
29866	Karman Curette N° 12 - 12 mm x 24 cm	box of 20



Length 20 cm with connector, 17 cm without

GIMA code	"BERKELEY" STYLE CURETTES CURVED - RIGID - SINGLE USE - STERILE	Minimum order
29870	Berkeley Curette N° 7 - Ø 7 mm x 20 cm	box of 20
29871	Berkeley Curette N° 8 - Ø 8 mm x 20 cm	box of 20
29872	Berkeley Curette N° 9 - Ø 9 mm x 20 cm	box of 20
29873	Berkeley Curette N° 10 - Ø 10 mm x 20 cm	box of 20
29874	Berkeley Curette N° 11 - Ø 11 mm x 20 cm	box of 20
29875	Berkeley Curette N° 12 - Ø 12 mm x 20 cm	box of 20



20 cm long CE 0120

• **29918 SUCTION SYRINGE - one-way self-blocking**  
Special syringe with adaptor suitable for use with Karman and Berkeley currettes of any brand.  
One-way self-blocking valve for safety and easier use.  
Single use and sterile, 60 ml capacity. **minimum order: 10 pcs.**

## AMNIOSCOPE: TUBES AND SET



To light amnioscope tubes you can use GIMA light source and Fibre Optic cables see page 165

• **29970 AMNIOSCOPE SET** includes 29971/2/3/5 and 30840  
Set of three amnioscopes tubes with shutter 200 mm length, 12, 16 and 20 mm diameter, and a clamp for fibre optic light. All tubes are graduated in centimetres, have a rounded end and are sterilizable with any method.

To light amnioscope set you can use GIMA 150 light source (code 30801), adapter (code 30830), F.O. cable (code 30818). Single tube need in addition adapter for clamp (code 30840) and clamp (29975).

- 29971 AMNIOSCOPE TUBE 12x24x200 mm
- 29972 AMNIOSCOPE TUBE 16x24x200 mm
- 29973 AMNIOSCOPE TUBE 20x24x200 mm
- 29975 ANGLED CLAMP FOR LIGHT CONNECTION

## OBSTETRIC STETHOSCOPES

- 32501 SHORT OBSTETRIC STETHOSCOPE - beech-wood
  - 32502 LONG OBSTETRIC STETHOSCOPE - beech-wood
- These stethoscopes are hand-craft manufactured, made in beech-wood accurately selected and sesasoned for at least 3 years. They have a perfect finish with a special transparent varnish. Packed in single box. Made in Germany.
- 32505 PINARD STETHOSCOPE - aluminium



32505 length 13 cm 32501 length 17 cm 32502 length 17 cm

## SCHULTZE SALPINGOGRAPH



**Riester**

CE 0124

- **29780 SALPINGOGRAPH**  
Instrument for hystero-salpingography and tubal insufflation, Prof. Gunter K. F. Schultze from Berlin type. Made in Germany.
- **29784 HYSTERO SYRINGE - 10 cc - LUER LOCK - spare**

**HYSTEOSALPINGOGRAPH: STANDARD ACCESSORIES**  
Watch for hystero Rack-bar      Hystero syringe 10 cc Luer Lock  
3 interchangeable cones (small-medium-large)

## SUCTION CUPS

- 29980 SILICONE SUCTION CUP - Ø 60x210 mm
  - 29981 SILICONE SUCTION CUP - Ø 70x210 mm
- Made in medical grade silicone

CE 0434

