

MAGNIFYING FLUORESCENT COLD LIGHTS - 3 DIOPTRES



CE
93/42

CEI 62-5

30740



30750
Magnifying glass
Ø 120 mm

Power supply unit integrated in the base with double safety fuse, lead retainer with spiral protection and bipolar light switch with green indicator light.


30751

SOLENORD
Suggested for dermatology. It is equipped with a biconvex optical magnifying glass diameter 120 mm. The cold fluorescent light allows prolonged use. Light source protection in plexiglass. Spring arm and 5 castors trolley allow easy adjustability. Fluorescent round lamp 230 V. Made in Italy.

GIMANORD
Magnifying glass self-illuminating fluorescent cold light 3 dioptrés (0.75X). Suggested for dermatology, it is equipped with a biconvex optical magnifying glass Ø 120 mm. The cold fluorescent light allows prolonged use. Polycarbonate screen for protection of the light source. Available also with wall support. Made in Italy.

Model	SOLENORD GIMANORD
Type of lighting	fluorescent
Consumption (Watt)	22
Life of the bulbs (hours)	6,000
Light intensity at 50 cm (lux)	550
Light intensity at 1 m (lux)	275
Colour temperature ± 5% at 1 m (K)	4 550
External diameter of the reflector (Ø mm)	230
Temperature increase at 1 m (°C)	0
Norms SOLENORD	IEC-601-1
Norms GIMANORD	IEC-598-2-25
Marks SOLENORD	CE 93/42
Marks GIMANORD	CE 2006/95
Weight (kg)	7
Operating voltage V-Hz	230 - 50
Adjustment height with castors base: min-max cm	60-180

GIMA code	MAGNIFYING LIGHTS
30740	SOLENORD on trolley
30741	SOLENORD wall
30750	GIMANORD desk model - with clamp
30751	GIMANORD on trolley
SPARES AND FIXING SYSTEMS	
30767	Spare fluorescent lamp - 22W
30771	Spare plexiglass shield
30769	Trolley - spare (A)
30774	Wall connection (B)
30775	Rail connection with clamp (C)
30776	Table connection (D)



LED MAGNIFYING LAMP - 5 DIOPTRES

• **30737 LOUPE LED MAGNIFYING LIGHT**
Smartly designed cold light LED lamp with 5-dioptres magnification. 5 wheels base and articulated arm allow easy movement. Manual in GB, FR, IT, ES, PT, DE Made in Europe.

MIMSAL

CE
93/42




TECHNICAL SPECIFICATIONS

Type of light:	LED
Dioptrés:	5
Lens:	Ø 12.7 cm
Operating voltage:	230 V - 50 Hz
Class II	
Electrical security:	UNI EN 60601-1 (1990)
Electromagnetic compatibility:	EN 60601-1-2 (1993)

Height
50+50 (arm) +
95 cm (trolley)

30737

• **30736 DERMALED MAGNIFYING LAMP**
Cold light LED lamp with 5-dioptres magnifying lens and brightness adjustment. It comes with a high stability 5 wheels base and an articulated arm for easy movement and fixing of the lamp. Suggested for dermatology and beauty treatments. Made in China.



TECHNICAL SPECIFICATIONS

Type of light:	LED
Voltage:	8 W
Dioptrés:	5
Operating voltage:	100-240 V - 50/60 Hz
Net/Gross weight:	10-12 kg

30736

