

SIGMA F.O. DIAGNOSTIC SETS

SIGMA otoscopes and ophthalmoscopes are high quality pocket diagnostic instruments designed for professional use in ENT and ophthalmology. Light weight and manufactured from impact resistant re-enforced plastic.

Supplied in deluxe rigid plastic cases or pouch (31524) with 14 disposable specula (7 with Ø 2.5 mm, 7 with Ø 4 mm).

- 31524 SIGMA F.O. OTOSCOPE 2.5 V - pouch
- 31525 SIGMA F.O. OTOSCOPE 2.5 V - case
- 31526 SIGMA F.O. LED OTOSCOPE - LED light 60,000 LUX - bulb estimated life 70,000 hours - case

Manufactured with high quality materials and top-notch workmanship ensures excellent reliability in every examination.

Battery life with Led bulb (31526) 6 time longer than with traditional bulb

- distal fiber optic illumination with super bright Xenon/Halogen light or LED light.
- swivel type window with 3X magnification lens for sharp images and minimal reflection
- insufflation port for pneumatic testing of tympanic mobility.

Connector port (code 31532) and insufflation bulb (code 31500) should be ordered separately.

- 31527 SIGMA F.O. OPHTHALMOSCOPE 2.5 V - case

Top quality ophthalmoscope with 5 different apertures for almost all ophthalmic examinations. High quality optics in mini format makes this instrument ideal for house calls and rounds.

- enhanced xenon/halogen super bright light
- 5 different apertures

Large spot: for normal fundoscopy

Small spot: for reduced reflection with small pupils

Semi circle: for reduced reflection with small pupils

Fixation star: for determination of central or eccentric fixation, suitable for examining children


Red-free filter: for contrast-enhancement when assessing fine vascular disorders

- range of lenses from -20 D to +20 D

- 31528 SIGMA F.O. OTO-OPHTHALMOSCOPE SET - case
- 34485 OTOSCOPE BULB - xenon - halogen 2.5V
- 31531 OPHTHALMOSCOPE BULB - xenon - halogen 2.5V
- 31532 PNEUMATIC CONNECTOR PORT
- 31500 INSUFFLATION BULB



31524



Range of lenses from -20D to +20D

+	1	2	3	4	6	8	10	15	20
-	1	2	3	4	6	8	10	15	20



31527



31525



31525



31527



31528

HEINE BETA® 200 RETINOSCOPE

- 31745 HEINE BETA 200 RETINOSCOPE - slit - 2.5V - C-034.10.118

With ParaStop® for finding the neutralization point (US pat. 5,859,687): precise, easy selection of a parallel beam

- dustproof: maintenance free
- single control for vergence and rotation: comfortable operation
- Integrated polarisation filter: eliminates stray light and reflexes for a brighter pupilling reflex
- streak or spot retinoscope: simply by changing the bulb
- XHL (xenon halogen technology): bright, white light. Bright fundus reflex, easy recognition of neutralization point
- ergonomic shape: protects the examiner's orbit from stray light



5 YEARS WARRANTY

