

# NON WOVEN SURGICAL SPECIALISTIC GOWNS - STERILE

## SURGICAL SPECIALISTIC GOWNS

Specialistic, water repellent, breathable gowns for a barrier effect to fluids and bacteria

- cotton cuffs 8 cm
- ultrasonic welding
- neck velcro closure
- 4 waist straps for internal and external closure

### Pre folded

For easy programmed dressing

### Welded reinforcement (only 21015, 21016)

Water-proof, breathable reinforcement on front and on sleeves ultrasonic welded.

All gowns include two dry-paper towels.

Conforms to UNI EN ISO 13795.

Made in Italy.



21020-21023

**HIGH QUALITY**  
50 g/m<sup>2</sup>

GIMA code	SURGICAL GOWNS STERILE	Minimum order
<b>SURGICAL GOWNS 35 g</b>		
21010	Size M - length 115 cm	box 30
21011	Size L - length 130 cm	box 30
21012	Size XL - length 135 cm	box 30
<b>SURGICAL GOWNS 50 g</b>		
21020	Size M - length 120 cm	box 30
21021	Size L - length 130 cm	box 30
21022	Size XL - length 140 cm	box 30
21023	Size XXL - length 150 cm	box 30
<b>REINFORCED GOWNS 35 g</b>		
21015	Size L - length 130 cm	box 30
21016	Size XL - length 135 cm	box 30



## DISPOSABLE GOWNS

### DISPOSABLE SURGICAL GOWNS - green

Disposable non-woven gowns, suitable when comfort and light protection are required.

Colour: green

Grams: 25 g/m<sup>2</sup>

One size: 110x136 cm

With cotton cuffs (26100 only).

26100



GIMA code	NON WOVEN GOWNS	Minimum order
26099	Disposable Gowns - sterile*	box of 50
26100	Disposable Gowns - non sterile	bag of 10

\*Sterilized by ETO

## POLYTHENE GOWNS

### DISPOSABLE POLYTHENE GOWNS - green

Surgical polythene gowns with neck and waist laces and cotton cuff.

Made in waterproof non-woven fabric.

One size 115x140 cm


Available sterile or non sterile.

Made in Italy.




GIMA code	POLYTHENE GOWNS	Minimum order
26092	Poly Gowns - sterile	box of 50
26093	Poly Gowns - non sterile	box of 100

## NON-WOVEN UNIFORM, JACKET AND PANTS



- 26076 UNIFORM COAT+PANT - S
  - 26077 UNIFORM COAT+PANT - M
  - 26078 UNIFORM COAT+PANT - L
  - 26079 UNIFORM COAT+PANT - XL
- minimum order: box of 50 - non sterile  
Short sleeve jacket "V" neck, with 3 pockets.  
Pants with rear pocket, closure by laces.



- 26023 JACKET WITH BUTTON - L
  - 26024 JACKET WITH BUTTON - XL
  - 26027 JACKET AND PANTS - L
  - 26028 JACKET AND PANTS - XL
- minimum order: box of 50 - non sterile  
Jacket long sleeve, neck and cuffs made of elastic polyester. Front opening with snap buttons. 2 front pockets.  
Pants with laces and rear pocket.

## COLONOSCOPY PANTS




- 21030 COLONOSCOPY PANTS
- minimum order: box of 50 - non sterile  
Colonoscopy pants, fit all sizes, length 110 cm.  
Elastic waist, rear opening with velcro fixing.  
Made in water-repellent non-woven SMS 50g/m<sup>2</sup>.

## DISPOSABLE KITS FOR PATIENT/VISITOR DRESSING - NON STERILE

- 26017 KIMONO - blue
- box of 100

Opaque non-woven SMS kimono 40 g/m<sup>2</sup> with short-sleeves. Front opening and closing at the waist with laces.  
One size: length 115 cm, width 140 cm.



- 26018 PATIENT DRESSING KIT - box of 25
- Kit includes kimono, round cap with elastic band and shoe-covers.



- 26019 VISITOR DRESSING KIT - box of 25
- Kit includes gown for visitor, round cap with elastic band, three layers mask with ear loops, polyethylene shoe-covers.