



3M TEGADERM™ AND TEGADERM™+PAD FILM DRESSINGS




Transparent film dressings for skin and wound protection. Frame style allows customization of shape and size to fit any site. Unique sacral shape is designed for challenging and difficult-to-dress areas. Sterile waterproof bacterial barrier, highly breathable.

Sterile, waterproof bacterial barrier composed by a non-adherent absorbent pad and an highly breathable transparent dressing.

GIMA code	3M code	3M TEGADERM™ STERILE - transparent film dressing	Minimum order
34900	1624W	TEGADERM™ - 6x7 cm - sterile	box of 100 pcs.
34901	1626W	TEGADERM™ - 10x12 cm - sterile	box of 50 pcs.
34902	1627	TEGADERM™ - 10x25 cm - sterile	box of 20 pcs.
34903	1628	TEGADERM™ - 15x20 cm - sterile	box of 10 pcs.

GIMA code	3M code	3M TEGADERM™ STERILE - transparent film dressing + PAD	Minimum order
34895	3582	TEGADERM™ + PAD 5x7.2 cm - sterile	box of 50 pcs.
34896	3586	TEGADERM™ + PAD 9x10 cm - sterile	box of 25 pcs.
34897	3589	TEGADERM™ + PAD 9x15 cm - sterile	box of 25 pcs.
34898	3591	TEGADERM™ + PAD 9x25 cm - sterile	box of 25 pcs.

3M TEGADERM™ I.V. AND ADVANCED WOUND DRESSINGS




Advanced securement dressing: with strong hold and gentle removal to help minimize the pain, anxiety and costs of catheter dislodgement. Waterproof, sterile.

Sterile hydrocolloid dressing with highly breathable cover impermeable to liquids, bacteria and viruses.

GIMA code	3M code	3M TEGADERM™ I.V. - STERILE advanced I.V. securement dressing	Minimum order
34945	1633	TEGADERM™ I.V. - 7x8.5 cm - fixing strips	box of 100
34946	1635	TEGADERM™ I.V. - 8.5x10.5 cm - fixing strips	box of 50

GIMA code	3M code	3M TEGADERM™ HYDROCOLLOID	Minimum order
34968	90002	HYDROCOLLOID™ 10x10 cm square	box of 5
34969	90007	HYDROCOLLOID™ 16x17 cm sacral	box of 6



Fast wicking, non-swelling polyurethane foam with a highly breathable film backing which prevents wound exudate strike-through and is a barrier to outside contamination while the dressing remains intact without leakage.

Antimicrobial, highly absorbent wound dressing composed of sterile calcium alginate, carboxymethylcellulose and an ionic silver complex. It releases antimicrobial ionic silver which very effectively reduces the microbial load caused by a broad spectrum of microorganism, including antibiotic resistant bacteria (MRSA and VRE). Sterile highly conformable to be used on exuding wounds.

GIMA code	3M code	3M TEGADERM™ FOAM - non adhesive or adhesive dressing	Minimum order
34954	90601	TEGADERM™ FOAM - 10x10 cm	box of 10
34955	90604	TEGADERM™ FOAM - 9x9 cm for tracheostomy	box of 10
34957	90619	TEGADERM™ FOAM ADHESIVE for heel and elbow - 14x14 cm	box of 5

GIMA code	3M code	3M TEGADERM™ ALGINATE and ALGINATE AG	Minimum order
34964	90112	TEGADERM™ ALGINATE - 10x10 cm	box of 10
34965	90312	TEGADERM™ ALGINATE AG -10x10 cm with silver	box of 10


• 34974 CAVILON™ NO STING BARRIER FILM 28 ml - 3346E - box of 12
Alcohol-free liquid barrier film that dries quickly to form a breathable, transparent coating on the skin. Designed to protect intact or damaged skin from urine, feces, body fluids, tape trauma, and friction. 72 hours protection.



• 34972 TEGADERM™ HYDROGEL 15g - 91110 - box of 10 single dose tubes
Wound filler, preservative-free, latex-free, easy to apply. Formulated to help provide a moist wound healing environment.



3M RESTON™ FOAM PAD



Self-adhering open cell foam designed to protect skin from splints, casts and prostheses. Can be used as a pressure reducing device. It can also be used as a padding for beds and chairs.

3M COBAN™ ELASTIC WRAP



Self-adherent cohesive elastic wrap

GIMA code	3M code	3M RESTON™ FOAM PAD	Minimum order
34920	1560	RESTON™ 30x20x1.1 cm	box of 10
34921	1561	RESTON™ 30x20x2.5 cm	box of 5
34922	1563	RESTON™ ROLL 10cmx5mx4mm	1 roll

GIMA code	3M code	3M COBAN™ SELF ADHERENT ELASTIC WRAP - LATEX	Minimum order
34923	1583	COBAN™ 7.5 cm x 2 to 4.57 m	box of 24
34924	1584	COBAN™ 10 cm x 2 to 4.57 m	box of 18