


AN INTEGRATED SYSTEM OF TREVIRA® CURTAINS



A COMPLETE SEPARATION SYSTEM

A complete system to create separations inside hospital rooms or outpatient wards.

System includes

4 STRUCTURES - include rings

- adjustable wall arms manual or electric
- window curtain rails
- cabin curtain rail
- trolleys with 1 or 2 curtains

3 SIZES TREVIRA® CURTAINS 6 COLOURS - do not include rings

Window curtains rails

Available in 2 size: 160 cm (45492) need curtains 175xh 145 cm and 210 cm (45493) need curtains 225x145 cm

Cabin curtain rails

Size: 80x80 need curtains 225xh 180 cm

Swing arm curtains - manual

Swing arm curtain is adjustable horizontally (180°) and can be bent down to 90° (45501 only) when not in use.

Telescopic tubular arm extension: open 210 cm, closed 88 cm (45500) or 100 cm (45501).

Rings allow curtains to slide smoothly.

Swing arm curtains - electric

Adjustable horizontally (180°) by electrical control.

Telescopic tubular arm extension: open 188 cm, closed 90 cm.

GIMA code	STRUCTURE: ARM OR TROLLEY WITHOUT CURTAINS
45492	Window curtain rail 160 cm
45493	Window curtain rail 210 cm
45495	Cabin rail 80x80 cm
45500	Adjustable arm
45501	Adjustable foldable arm
45503	Electrical adjustable arm
45505	Trolley for 1 curtain - foldable
45506	Trolley for 2 curtains - foldable

Trolley curtains - 1 or 2 curtains

Trolley made of stainless steel ISO 304 with 4 revolving castors (2 lockable). It allows the use of 1 curtain 225xh 145 cm (45505) or 2 curtains (45506) 1 long 225xh 145 cm and 1 short 175xh 145 cm. Arm can be bent down to 90° when not in use and curtain remains tied up.

Trevira® curtains

Trevira® curtains are flame-proof (Class I), antiallergic, avoid bacterial growth, waterproof and highly resistant to abrasion. Thick weft to avoid dust infiltration and improve privacy, washable up to 60°C. Available in 6 colours to match any need.

Custom-made curtains are available on request.

Made in Italy.



TREVIRA CURTAINS			
FOR WALL ARMS AND CABIN	FOR TROLLEYS AND WINDOWS		COLOUR
	Long trolley arm/rail: 225xh145 cm	Short trolley arm/rail: 175xh145 cm	
45510	45520	45530	Light blue
45511	45521	45531	Blue
45512	45522	45532	White
45513	-	-	Beige
45514	-	-	Green
45515	45525	45535	Peach

Light blue

Blue

White

Beige

Green

Peach

SET OF 18 RINGS - SPARE	
45508	Open rings for window rails and cabin
45509	Closed rings for wall arms and trolleys