



Product Code	28402
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	1456069
CND	Z120601
EAN/UPC	8019781501097
UMDNS	13762
PARAF	936980604

## Description

### T-ONE MEDISPORT ELECTROTHERAPY - 4 channels

Professional equipment specifically designed for sport specialist users. Muscle stimulation specialized programs improve sport performance. It increases strength and helps recovery of sport injuries with specific treatment protocols.

#### Main features

- constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- four independent channels unit (up to 8 connectable electrodes)
- compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect;
- monophasic wave generator (only for iontophoresis programs)
- 2+2 users programs
- backlight graphic colour display
- user feedback: the device recognizes patient connection
- soft-touch easy to use keyboard
- on-line 5 languages guide (GB, FR, IT, ES, DE) on display
- 4,8 V rechargeable battery pack
- personalized programs for man and woman
- maximum intensity 240 mApp, with 1 mA step regulation

Size: 150x90x35 mm

## Technical Specifications

- BEAUTY (18 PROGRAMMES)
  - Firming up
  - Toning up
  - Definition
  - Modelling
  - Microlifting
  - Lipolysis abdomen
  - Lipolysis thighs
  - Lipolysis glutei and hips
  - Lipolysis arms
  - Tissue elasticity
  - Capillarization
  - Heaviness in legs
  - Face Microlifting Handy Trode
  - Breast Microlifting Handy Trode
  - Décolleté Microlifting Handy Trode
  - Face wrinkles Handy Trode
  - Forehead Wrinkles Handy Trode
  - Collagen

- REHA (10 PROGRAMMES)

- Ionophoresis L
- Ionophoresis M
- Ionophoresis H
- MENS Microcurrent
- Haematomas
- Oedema
- Tens sequential
- Tens Burst
- Atrophy prevention
- Atrophy

- TENS (16 PROGRAMMES)

- Tens eapid
- Tens Endorphinic
- Tens Max Values
- Inflammatory
- Neck pain
- Backache/Sciatic pain
- Sprains/Bruises
- Vascularization
- Muscle Relaxant
- Hand and wrist pain
- Plantar stimulation
- Epicondylitis
- Epitroclea
- Periarthritis
- Tens Burst
- Microcurrent

- NEMS (14 PROGRAMMES)

- Warming up
- Resistance
- Resistance strength
- Basic strength
- Fast strength
- Explosive strength
- Deep capillarization
- Muscle recovery
- Agonist/Antagonist
- Sequential tonic contractions
- Sequential phasis contractions
- Muscle Relaxant
- Deep massage
- EMS
- 30' Resistance
- 30' Resistance strength
- 20' Base strength
- 20' Fast strength
- 20' Explosive strength
- 20' Hypertrophy

- NEM

- Tens mem 1
- Tens mem 2
- Tens mem 3
- Tens mem 4
- Tens mem 5
- NEMS mem 1
- NEMS mem 2
- NEMS mem 3
- NEMS mem 4
- NEMS mem 5
- NEMS alternated
- Battery test

- 2+2

- Tens rapid
- Tens endorphinic
- firming up
- Toning up
- Definition
- Modelling
- Microlifting
- Muscle Relaxant
- Massage
- Capillarization

Inflammatory  
Resistance, resistance strength  
Basic strength, fast strength  
Explosive strength, hypertrophy  
Face Microlifting Handy Trode  
Breast microlifting Handy Trode  
Décolleté microlifting Handy Trode  
Face wrinkles Handy Trode  
Forehead wrinkles Handy Trode

## Standard accessories

4 connection cables (8 connectable electrodes)  
Kit of pre-gelled self-adhesive electrodes (square, rectangular, round)  
lonophoresis kit: elastic belt, 2 silicone electrodes, 2 sponges, carrying bag  
Internal Ni-Mh rechargeable battery pack and battery charger  
Electrode position manual and user manual (GB, FR, IT, ES). On request DE, PT